

Instructions for use

PRODUCT NAME: INFUSION SET FOR SINGLE USE

BRAND NAME: CHIRAPLUS - P / CHIRAPLUS - G

INTENDED PURPOSE:

- Infusion set for single use is intended for administration liquids and drugs into the circulation system by using of intravenous catheters and cannulas.

INTENDED USERS:

- Health care professionals with relevant knowledge of use these types of medical devices. Infusion sets are simple and well-known products.

INTENDED PATIENT POPULATION:

- The type of Infusion set for single use containing phthalates should not be used for treatment of pregnant or nursing women, infants and children.



CONTRAINDICATIONS:


- Not to be used other than as indicated in intended purpose.
- Administration of highly viscous fluids.
- Blood Transfusion.

INSTRUCTIONS FOR USE:

- Remove Infusion set from unit package, closed the roller clamp with flow regulator.
- Prepare container with infusion fluid.
- Remove protector from piercing spike. Insert the spike to its full length into the top of the Infusion fluid container. Suspend/hang the container.
- Squeeze the drip chamber and set fluid level.
- Hang the connector end in the groove provided in the roller clamp.
- Open the airvent cap of piercing spike.
- Open the clamp. The fluid starts running into the tubing and priming takes automatically up to the connector end. Close the clamp.
- Remove protective cap from the connector end and connect with venous access device.
- Open roller clamp and adjust flow of infusion with flow regulator.

CAUTIONS:

- Read instructions before use. The product should be used according to the Instructions for use.
- Check the product and package integrity prior to use and also the expiry date.
- Do not use if protective cap is loose or missing.
- Do not use if package is opened or damaged.
- For pressure or gravity feed, as marked on product packaging (G = Gravity, P=Pressure).
- Use the product immediately after opening the individual blister packing.
- Close the airvent during periods of interrupted infusion therapy.
- Do not re-sterilize. Discard the set after single use.
- This product contains phthalates (**DEHP**); as marked with symbol .
- This product contains phthalates DEHP. The plasticizer DEHP is present in the fluid pathway of the Infusion set tubing. Leaching of DEHP is possible on contact with lipids by administering lipophilic drugs. Lipophilic drugs should not be administered with this type of infusion set during the treatment of pregnant or nursing women, infants and children.
- Infusion set with Flashball contains **LATEX** and is marked with symbol .
- Reuse and cleaning of the product may alter their structural and mechanical properties. It may lead to infections or other illness/injury.
- The device should be replaced and dispose as per facility approved protocol or guidelines.
- Store at 0-40°C temperature, avoid excessive heat, protect from direct sunlight and moisture.
- **CHIRANA T. Injecta** will not responsible for any direct incident or consequential damages resulting from reuse or improper use of the product.
- Report any health incident associated with medical device to healthcare professional or manufacturer.



CHIRANA T. Injecta, a.s.

Nám. Dr. Schweitzera 194
916 01 Stará Turá
Slovak Republic
Tel.: +421/32/7709 991

www.t-injecta.com



Fragile, handle with care

Keep away from sunlight

Keep dry



This product meets the requirements of the Regulation 2017/745



Caution



Consult instructions for use



Do not re-use



Do not re-sterilize



Do not use if package is damaged



Catalogue number



Batch code



Date of manufacture



Use-by date



Sterilized using ethylene oxide




Single sterile barrier system



Single sterile barrier system with protective packaging outside



Distributor



Temperature limit



Manufacturer



Non-pyrogenic



Contains or presence of phthalate



Contains hazardous substances



Contains or presence of natural rubber latex



LATEX Free



Liquid filter with pore size



Drops per millilitre



Unique device identifier



Medical device