

CHIRANA® Insulin and Tuberculin syringes



According to EN, ISO standards | Quality management certified according to EN ISO 13485:2016 | CE marked

- ✓ needle permanently fixed to syringe body, integrated
- ✓ precision lancet point grinding, silicone treated tip for smooth penetration
- ✓ no dead space volume

| Type | Size | Size G/'' | Size metric | Pieces per Box | Pieces per Case | Item Code |
|--------------|------|--------------|-------------|----------------|-----------------|------------|
| INSULIN U40 | 1 ml | 29G x 1/2'' | 12,7 mm | 100 | 2000 | CHINS4129 |
| TUBERCULIN | | 29G x 1/2'' | 12,7 mm | 100 | 2000 | CHTUB0129 |
| | | Luer | | 100 | 2000 | CHTUB01 |
| INSULIN U100 | | 27G x 1/2'' | 12,7 mm | 100 | 2000 | CHINS0127 |
| | | 29G x 5/16'' | 8 mm | 100 | 2000 | CHINS01298 |
| | | 29G x 1/2'' | 12,7 mm | 100 | 2000 | CHINS0129 |

CHIRANA® Insulin and Tuberculin syringes

| Type | Size | Size G/'' | Size metric | Pieces per Box | Pieces per Case | Item Code |
|--------------|--------|--------------|-------------|----------------|-----------------|-------------|
| INSULIN U100 | 1 ml | 30G x 1/2'' | 12,7 mm | 100 | 2000 | CHINS0130 |
| | | 30G x 5/16'' | 8 mm | 100 | 2000 | CHINS01308 |
| | | 31G x 4/16'' | 6 mm | 100 | 2000 | CHINS01316 |
| | | Luer | | 100 | 2000 | CHINS01 |
| | 0,5 ml | 29G x 1/2'' | 12,7 mm | 100 | 2000 | CHINS00529 |
| | | 30G x 1/2'' | 12,7 mm | 100 | 2000 | CHINS00530 |
| | | 30G x 5/16'' | 8 mm | 100 | 2000 | CHINS005308 |
| | | 31G x 5/16'' | 8 mm | 100 | 2000 | CHINS005318 |
| | | 31G x 4/16'' | 6 mm | 100 | 2000 | CHINS005316 |
| | 0,3 ml | 29G x 1/2'' | 12,7 mm | 100 | 2000 | CHINS00329 |
| | | 30G x 1/2'' | 12,7 mm | 100 | 2000 | CHINS00330 |
| | | 31G x 5/16'' | 8 mm | 100 | 2000 | CHINS003318 |
| | | 31G x 4/16'' | 6 mm | 100 | 2000 | CHINS003316 |

- ✓ possibility of 10 pcs POLY BAG packaging
- ✓ individual OEM manufacturing in reasonable quantity available