

MEDOJECT® fine Insulin Pen needles



According to EN, ISO standards | Quality management certified according to EN ISO 13485:2016 | CE marked

- ✓ compatible with most of type A Insulin pens
- ✓ precision lancet point grinding, silicone treated tip for smooth penetration
- ✓ shorter needles reduce risk of a painful injection into the muscle
- ✓ comfortable injection experience
- ✓ sterile, individually packed

Size	Color code	Pieces per Box	Pieces per Case	Item Code
29G x10 mm	yellow	100	2000	CHPN2910
29G x 12,7 mm	pink	100	2000	CHPN29127
30G x 6 mm	dark blue	100	2000	CHPN3006
30G x 8 mm	light blue	100	2000	CHPN3008
31G x 4 mm	green	100	2000	CHPN3104
31G x 5 mm	dark violet	100	2000	CHPN3105

MEDOJECT® fine Insulin Pen needles

Size	Color code	Pieces per Box	Pieces per Case	Item Code
31G x 6 mm	dark blue	100	2000	CHPN3106
31G x 8 mm	light blue	100	2000	CHPN3108
32G x 4 mm	green	100	2000	CHPN3204
32G x 5 mm	dark violet	100	2000	CHPN3205
32G x 6 mm	dark blue	100	2000	CHPN3206
33G x 4 mm	green	100	2000	CHPN3304

- ✓ intended for Insulin application via most of type A Insulin pens
- ✓ made of polypropylene hub and stainless steel tube
- ✓ thin wall tube material