

## CHIRANA® 2-part disposable syringes



According to EN, ISO standards | Quality management certified according to EN ISO 13485:2016 | CE marked

- ✓ ergonomic design, smooth flow, highly transparent barrel, tight
- ✓ black extended graduation, easy to read, double stop lines, guaranteed 0 ml position
- ✓ designed tip for best connection, double stop lines, scratch-free for glove manipulation
- ✓ available in Luer tip
- ✓ made of polypropylene/polyethylene, rubber-free
- ✓ pyrogen-free, latex-free, DEHP-free, PVC-free, non-toxic
- ✓ sterile, individually blister packed

| Size  | Pieces per Box | Pieces per Case | Item Code LUER |
|-------|----------------|-----------------|----------------|
| 2 ml  | 100            | 2800            | CH002L         |
| 5 ml  | 100            | 1900            | CH005L         |
| 10 ml | 100            | 1200            | CH010L         |
| 20 ml | 80             | 960             | CH020L         |