

## CHIRAFLEX® IV Cannulas



According to EN, ISO standards | Quality management certified according to EN ISO 13485:2016 | CE marked

- ✓ intended for peripheral vascular application
- ✓ compatible with all standard medical devices for wide application
- ✓ catheter made from Teflon material, available in SAFETY version
- ✓ 50 pcs per box, 1000 pcs per case

Size G	Diameter (mm)	Length (mm)	Color code	Item Code WITH PORT WITH WINGS	Item Code WITHOUT PORT WITH WINGS	Item Code WITHOUT PORT WITHOUT WINGS
26	0,6	19	purple	CHFX0102619	CHFX0202619	CHFX0302619
24	0,7	19	yellow	CHFX0102419	CHFX0202419	CHFX0302419
22	0,9	25	blue	CHFX0102225	CHFX0202225	CHFX0302225
20	1,1	32	pink	CHFX0102032	CHFX0202032	CHFX0302032
18	1,3	45	green	CHFX0101845	CHFX0201845	CHFX0301845
17	1,5	45	white	CHFX0101745	CHFX0201745	CHFX0301745
16	1,7	45	grey	CHFX0101645	CHFX0201645	CHFX0301645
14	2,1	45	orange	CHFX0101445	CHFX0201445	CHFX0301445