

## CHIRAFLEX® Scalp vein sets



According to EN, ISO standards | Quality management certified according to EN ISO 13485:2016 | CE marked

✓ intended for shorter time infusion

✓ stainless steel needle with length of 20 mm (4/5"), 30 cm plastic tubing

Size G	Diameter (mm)	Color code	Item Code
27	0,4	grey	CH260001127
25	0,5	orange	CH260001125
23	0,6	blue	CH260001123
22	0,7	black	CH260001122
21	0,8	green	CH260001121
20	0,9	yellow	CH260001120
19	1,1	ivory	CH260001119
18	1,2	pink	CH260001118