

## CHIRAHem™ Transfusion set



According to EN, ISO standards | Quality management certified according to EN ISO 13485:2016 | CE marked

✓ intended for blood transfusion administration      ✓ all sets have 200 µm blood filter

Item	Chamber Volume (ml)	Packaging type	Pieces per Case	Item Code
<b>CHIRAHem</b> Tubing 1500mm Luer Lock, without flashball	17,0	blister	100	CHT314
<b>CHIRAHem DEHP free</b> Tubing 1500mm Luer Lock, without flashball	17,0	blister	100	CHT514
<b>CHIRAHem DEHP free w/o air ventil</b> Tubing 1500mm Luer Lock, without air ventil	17,0	blister	100	CHT515
<b>CHIRAHem DEHP free</b> Tubing 1500mm Luer Lock, with flashball and without air ventil	17,0	blister	100	CHT516
<b>CHIRAHem DEHP free</b> Tubing 1800mm Luer Lock, without flashball	17,0	blister	100	CHT517
<b>CHIRAHem DEHP free</b> Tubing 1500mm Luer Lock, with flashball and with air ventil	17,0	blister	100	CHT518
<b>CHIRAHem</b> Tubing 1500mm Luer Lock, with flashball	16,0	poly bag	300	CHT320
<b>CHIRAHem with needle</b> Tubing 1500mm Luer Lock, with flashball + needle 1,2x40mm (18Gx1 ½")	17,0	poly bag	300	CHT321
<b>CHIRAHem with metal spike</b> Tubing 1500mm Luer, with flashball + needle 1,2x40mm (18Gx1 ½")	17,0	poly bag	300	CHT322
<b>CHIRAHem with plastic spike</b> Tubing 1500mm Luer, with flashball + needle 1,2x40mm (18Gx1 ½")	17,0	poly bag	300	CHT323