

CHIRAPLUS™ Infusion set



According to EN, ISO standards | Quality management certified according to EN ISO 13485:2016 | CE marked

- ✓ all sets include 15 µm, ventilation port and luer lock adapter
- ✓ CHIRAPLUS G intended for gravity infusion includes flash ball
- ✓ CHIRAPLUS P intended for pressure infusion, does not include flash ball

| Item | Chamber Volume (ml) | Packaging type | Pieces per Case | Item Code |
|--------------------------------------|---------------------|----------------|-----------------|-----------|
| CHIRAPLUS P | 10,5 | blister | 100 | CHP336 |
| CHIRAPLUS P DEHP free | 10,5 | blister | 100 | CHP536 |
| CHIRAPLUS P 12 | 12,0 | blister | 100 | CHP306 |
| CHIRAPLUS P 12 DEHP free | 12,0 | blister | 100 | CHP506 |
| CHIRAPLUS P 190 | 10,5 | blister | 100 | CHP339 |
| CHIRAPLUS P 190 SAFETY | 10,5 | blister | 100 | CHP338 |
| CHIRAPLUS P DEHP Free Autofusion 180 | 15,0 | blister | 300 | CHP576 |

CHIRAPLUS™ Infusion set

| Item | Chamber Volume (ml) | Packaging type | Pieces per Case | Item Code |
|-------------------------------------|---------------------|----------------|-----------------|-----------|
| CHIRAPLUS G | 10,5 | blister | 100 | CHG331 |
| CHIRAPLUS G | 10,5 | poly bag | 500 | CHG317 |
| CHIRAPLUS G DEHP free | 10,5 | blister | 100 | CHG531 |
| CHIRAPLUS G 12 | 12,0 | blister | 100 | CHG301 |
| CHIRAPLUS G 12 DEHP free | 12,0 | blister | 100 | CHG501 |
| CHIRAPLUS G 12 DEHP free OPAQUE 190 | 12,0 | blister | 100 | CHG901 |
| CHIRAPLUS G | 8,5 | poly bag | 500 | CHG318 |
| CHIRAPLUS G DEHP free | 8,5 | poly bag | 500 | CHG518 |
| CHIRAPLUS G with needle | 8,5 | poly bag | 500 | CHG319 |
| CHIRAPLUS G with metal spike | 8,5 | poly bag | 500 | CHG316 |