

## CHIRAVAC™ Blood Collection needles



According to EN, ISO standards | Quality management certified according to EN ISO 13485:2012 | CE marked

- ✓ screw-on needle hub compatible with all available holders
- ✓ precision lancet point grinding, silicone treated tip for smooth penetration
- ✓ easy-to-open packaging and manipulation
- ✓ comfortable collection experience
- ✓ sterile, individually packed

Size metric	Size G/''	Color code	Pieces per Box	Pieces per Case	Item Code
0,60 x 25 mm	23G x 1''	blue	100	1000	CHBCN23100
0,60 x 38 mm	23G x 1 1/2''	blue	100	1000	CHBCN23112
0,70 x 25 mm	22G x 1''	black	100	1000	CHBCN22100
0,70 x 38 mm	22G x 1 1/2''	black	100	1000	CHBCN22112
0,80 x 25 mm	21G x 1''	green	100	1000	CHBCN21100
0,80 x 38 mm	21G x 1 1/2''	green	100	1000	CHBCN21112

## CHIRAVAC™ Blood Collection needles

Size metric	Size G/”	Color code	Pieces per Box	Pieces per Case	Item Code
0,90 x 25 mm	20G x 1”	yellow	100	1000	CHBCN20100
0,90 x 38 mm	20G x 1 1/2”	yellow	100	1000	CHBCN20112
1,20 x 25 mm	18G x 1”	pink	100	1000	CHBCN18100
1,20 x 38 mm	18G x 1 1/2”	pink	100	1000	CHBCN18112

- ✓ intended for blood collection with vacuum tubes
- ✓ made of polypropylene hub and stainless steel tube
- ✓ thin wall tube material
- ✓ color coded needle hub