





Surgical Sutures & Meshes

State of the art





When experience matters

Our surgical sutures are manufactured in line with state of the art technology. Each individual production process is subject to strict scientific quality control.

Continuous inspection of both mechanical and microbiological characteristics is given top priority.

We have more than 85 years of experience in manufacturing surgical sutures.

Our team understands the specific needs of medical professionals by providing them with access to the most complete portfolio of surgical sutures.

The result is the current modern product range of absorbable and non-absorbable sutures, which consists of more than 1500 individual surgical suture products.

Production takes place in accordance with GMP guidelines, and a comprehensive quality management system in line with ISO 13485 has been implemented.

All products have been CE marked since 1999.



Contents

Overview of Surgical Sutures

Absorbable	Page 6
Non-absorbable	Page 10
Labelling and Packaging	Page 12
Needle geometry	Page 16
Needle range	Page 19

Surgical Sutures

Absorbable Sutures	Page 27
Chirasorb rapid	Page 29
Chirlac rapid	Page 39
Monolac	Page 51
Chirasorb	Page 63
Chirasorb Plus	Page 83
Chirlac	Page 99
Polydox	Page 119

Non-absorbable Sutures	Page 131
Chiraflon	Page 133
Chiralen	Page 139
Silon monofilament	Page 153
Tervalon	Page 165
Silon braided/Orsilon braided	Page 177
Silk	Page 187
Steel Wire	Page 195

Special products

Fix-clip	Page 201
Minimal invasive surgery	Page 201
Ophthalmic - Microsurgery	Page 202
Eye-needles	Page 206
Eye-needles - Double spring	Page 207
Eye-needles - Regular	Page 213

Surgical Meshes

Partially absorbable meshes	Page 216
Capromesh meshes	Page 216
Non-absorbable meshes	Page 218
Chiralen meshes	Page 218
Chiralen meshes light	Page 218

Veterinary

Absorbable bottles	Page 222
Chirasorb	Page 222
Chirlac	Page 222
Polydox	Page 222

Non-absorbable bottles	
Silon monofilament	Page 223
Tervalon	Page 223
Silon braided	Page 223
Silk	Page 223

The CTI group	Page 225
---------------	----------

Surgical Sutures - Colour coding

Chirlac rapid

Chirasorb

Chirlac

Polydox

UK03
4



Surgical Sutures - Overview



Silk

Monolac

Tervalon

Chiralen

Chirasorb rapid

Silon braided/
Orsilon braided

Chiraflon

Silon monofilament

Chirasorb Plus

Steel Wire

Surgical Sutures - Absorbable threads

UK03
6

CHIRASORB rapid

Intended for use in Skin closure, Paediatric surgery, Episiotomies, Circumcision, Closure of oral mucosa and Ophthalmic surgery for conjunctival sutures.



Composition	Colour	Range	Characteristics
Braided and coated suture made from 90% glycolide and 10% L-lactide	Violet and Natural	6/0 to 2 USP	Absorbable, multifilament, coated, synthetic, high tensile strength, reduced capillary effect, short-term healing

CHIRLAC rapid

Intended for use in Gynaecology, Skin Closure, Ophthalmic, Odontology and Paediatric surgery.



Composition	Colour	Range	Characteristics
Braided and coated suture made of Polyglycolic acid (PGA)	Violet and Natural	6/0 to 4 USP	Absorbable, multifilament, coated, synthetic, high tensile strength, reduced capillary effect, short-term healing

MONOLAC

Intended for use in general soft tissue in all specialities, Skin Closure, Plastic, Gynaecology, Urology and Gastrointestinal surgery.



Composition	Colour	Range	Characteristics
Monofilament suture made of Glycolide and ε-caprolactone	Violet and Natural	6/0 to 2 USP	Absorbable, monofilament, synthetic, high tensile strength, excellent handling properties, mid-term healing

CHIRASORB

Intended for use in General surgery, Gynaecology, Ophthalmic, Skin Closure, Muscular layer and Ligation.



Composition	Colour	Range	Characteristics
Braided and coated suture made from 90% glycolide and 10% L-lactide	Violet and Natural	10/0 to 6 USP	Absorbable, multifilament, coated, synthetic, high tensile strength, reduced capillary effect, mid-term healing



CHIRASORB Plus

Antibacterial

Intended for use in General surgery, Gynaecology, Skin Closure, Muscular layer and Ligation.



Composition	Colour	Range	Characteristics
Antibacterial, braided and coated suture made from 90% glycolide, 10% L-lactide and chlorhexidine diacetate	Violet and Natural	6/0 to 4 USP	Absorbable, multifilament, anti-bacterial coating, synthetic, high tensile strength, reduced capillary effect, mid-term healing

CHIRLAC

Intended for use in General surgery, Gynaecology, Ophthalmic, Skin Closure, Muscular layer and Ligation.



Composition	Colour	Range	Characteristics
Braided and coated suture made of Polyglycolic acid (PGA)	Violet and Natural	9/0 to 6 USP	Absorbable, multifilament, coated, synthetic, high tensile strength, reduced capillary effect, mid-term healing

POLYDOX



Intended for use in general soft tissue approximation, Orthopaedics, Fascia Closure, Urology, Gynaecology and Plastic surgery. Recommended for long-term healing.






Composition	Colour	Range	Characteristics
Monofilament suture made of Poly-p-dioxanone	Violet and Natural	7/0 to 4 USP	Absorbable, monofilament, synthetic, high tensile strength, excellent knotting properties, long-term healing

Surgical Sutures - Absorbable threads

SHORT-TERM HEALING



MID-TERM HEALING




UK03
8




Surgical Sutures - Overview

MID-TERM HEALING



LONG-TERM HEALING



UK03

9



VITREX MEDICAL SUTURES	Strength Retention profile	Tensile Strength	Absorption time
CHIRASORB RAPID	SHORT-TERM HEALING 1 day 100% 5 days 50% 10 - 14 days 0%	12 days	42 days
CHIRLAC RAPID	SHORT-TERM HEALING 1 day 100% 5 days 60% 10 - 14 days 0%	12 days	42 days
MONOLAC	MID-TERM HEALING 1 day 100% 7 days 60% (violet) 50% (natural) 14 days 30% (violet) 20% (natural) 28 days 0%	30 days	90 - 120 days
CHIRASORB	MID-TERM HEALING 1 day 100% 14 days 75% 21 days 50% 30 days 25% 56 - 70 days 0%	35 days	56 - 70 days
CHIRASORB PLUS	MID-TERM HEALING 1 day 100% 14 days 75% 21 days 50% 56 - 70 days 0%	35 days	56 - 70 days
CHIRLAC	MID-TERM HEALING 1 day 100% 14 days 65% - 75% 21 days 45% - 50% 80 - 90 days 0%	28 - 35 days	80 - 90 days
POLYDOX	LONG-TERM HEALING 1 day 100% 45 days 50% 180 - 240 days 0%	90 - 98 days	180 - 240 days



Surgical Sutures - Non-absorbable threads

CHIRAFLON

Intended for use in all surgical specialities, especially Cardiovascular, Plastic, Reconstructive, Ophthalmic, Orthopaedics and Traumatology (optimal, temporary, easily removable suture).



Composition	Colour	Range	Characteristics
Monofilament suture made of Polyvinylidene Fluoride (PVDF) - similar to Teflon	Blue	8/0 to 3 USP	Non-absorbable, monofilament, synthetic, no capillary effect, very smooth, excellent knotting properties, easy to handle, physiologically inert

CHIRALEN

Intended for use in all surgical specialities, especially Cardiovascular, Plastic, Reconstructive, Ophthalmic, Orthopaedics and Traumatology (tendons).



Composition	Colour	Range	Characteristics
Monofilament suture made of Polypropylene	Blue and Black	10/0 to 2 USP	Non-absorbable, monofilament, synthetic, no capillary effect, smooth, non-porous surface, suitable for stitches in infected wounds

SILON monofilament

Intended for use in Skin Closure, Plastic, Reconstructive, General and Ophthalmic surgery.



Composition	Colour	Range	Characteristics
Monofilament suture made of Polyamide 6	Blue and Black	10/0 to 2 USP	Non-absorbable, monofilament, synthetic, very smooth, excellent knotting properties, easy to handle, loss of tensile strength after long-term implantation

TERVALON

Intended for use in Cardiovascular, Plastic, Reconstructive, General surgery and whenever non-absorbability and low-friction passage through tissues are needed.




Composition	Colour	Range	Characteristics
Braided and coated suture made of Polyester	Green and White	6/0 to 6 USP	Non-absorbable, multifilament, synthetic, high tensile strength, no capillary effect, well tolerated by the tissue, low friction

Surgical Sutures - Non-absorbable threads

SILON braided ORSILON braided

Intended for use in Skin Closure, Subcutaneous tissue, Transversely cut muscles, Fascia, Peritoneum, Gastrointestinal Tract and Blood vessels ligation.



UK03
11

Composition	Colour	Range	Characteristics
Braided suture made of Polyamide	White and Black	5/0 to 5 USP	Non-absorbable, multifilament, synthetic, high tensile strength, excellent knotting properties, loss of tensile strength after long-term implantation

SILK

Intended for use in Skin Closure, Subcutaneous tissue, Fascia, Peritoneum and Gastrointestinal Tract.



Composition	Colour	Range	Characteristics
Braided and coated suture made of Silk	Blue and Black	9/0 to 5 USP	Non-absorbable, multifilament, synthetic, no capillary effect, loss of tensile strength after long-term implantation

STEEL Wire

Intended for use in Abdominal wound closure, Sternal closure, Orthopaedic procedures and Tendons.



Composition	Colour	Range	Characteristics
Monofilament suture made of 316LVM Stainless steel	Steel	1 to 7 USP	Non-absorbable, monofilament, superior steel, flexible wire, excellent tissue penetration, high tensile strength


European Pharmacopoeia

The publication of the European Pharmacopoeia (Ph. Eur.) resulted in a standardization of the classification of threads and gauges of surgical suture material.

The Ph. Eur. specifies a decimal classification and metric coding for the gauges. It quotes a thread diameter of 1/10 mm, so that the thread gauge of 3.5 has a diameter of 0.35 mm.

The following table of the gauges gives a comparison of the metric and the USP gauge coding.

EP/ Metric	USP	Ø (mm)	Ø (mm) Absorbable monofilaments
0.2	10-0	0,020-0,029	
0.3	9-0	0,030-0,039	
0.4	8-0	0,040-0,049	
0.5	7-0	0,050-0,069	0,050-0,094
0.7	6-0	0,070-0,099	0,095-0,149
1	5-0	0,100-0,149	0,150-0,199
1.5	4-0	0,150-0,199	0,200-0,249
2	3-0	0,200-0,249	0,250-0,339
3	2-0	0,300-0,349	0,340-0,399
3.5	0	0,350-0,399	0,400-0,499
4	1	0,400-0,499	0,500-0,570
5	2	0,500-0,599	0,571-0,610
6	3+4	0,600-0,699	
7	5	0,700-0,799	
8	6	0,800-0,899	




Labelling and Packaging

UK03
12



Surgical Sutures - Overview



LOT number YPPTTRR
 Y - Shelf life in years
 PP - Production number
 TT - Week of production
 RR - Year of production

Identification label

For quick registration

The label contains:

Name of material

Batch / LOT number

Expiration date

Catalogue number

Manufacturer

CHIRASORB
LOT 5080814
2022-08 LV 0224
CHIRANA T. Injecta





Packaging

The packaging of CTI's surgical suture material is stable, user-friendly and functional. It provides optimal mechanical and bacterial safety. The sutures are ready for use in a sterile foil or aluminium pack.

Each primary sales unit contains 24 individual packs, whereas some contain 12 individual packs (Loop sutures or Microsurgery sutures) or 20 individual packs (Chirapack sutures).

CTI short threads

Pre-cut length surgical sutures.

CTI atraumatic sutures

Pre-cut length surgical sutures, armed with a drilled end needle, ready for use. Removable needles are marked with an A behind the needle description.


CTI Chirapack

2.5 or 3 meters of suture material on a thread spool, ready for use. Each sales unit contains 20 individual packs.



How to unpack the suture - Alufoil

UK03
14



1. Open the non-sterile outer package



2. Transfer to sterile area



3. Tear off the opening system, indicated by the dotted line



4. Take immediate, direct hold of the needle. Pull out the suture towards the right




How to unpack the suture - Blisterfoil


UK03
15




Surgical Sutures - Overview




1. Open the non-sterile outer package



2. Transfer to sterile area



3. Tear off the opening system, indicated by the dotted line



4. Take immediate, direct hold of the needle. Pull out the suture towards the right




Surgical Sutures - Needle geometry

UK03
16



Surgical Sutures - Overview



Taper point
(round body)



Trapezoid spatulated




Conventional cutting
(inside cutting edge)



Tapercut



Reverse trapezoid spatulated




Reverse cutting
(outside cutting edge)




Blunt point





Needle curvature	Needle Body	Needle Point	Speciality
V - 1/4	R - Round body	T - Tapercut	S - Slim, rectangular body
H - 1/2	S - Reverse cutting	N - Blunt	A - Take-off needle
D - 3/8	SC - Conventional cutting	M - Micropoint	X - Extra strong needle
F - 5/8	LZ - Spatula		Black - Black needle
J - Hook	RF - Round, flattened body		
A - Asymmetric			
K - Ski			
G - Straight			




Surgical Sutures - Needle geometry

Needle curvature

UK03
18



Surgical Sutures - Overview




F	5/8 circle	Anal region, Nasal cavity, Urogenital tract, Pelvic region or other confined spaces.
H	1/2 circle	Gastrointestinal tract, Respiratory tract, Urogenital tract, Cardiovascular, Fascia, Muscle, Nasal cavity, Oral cavity, Eye, Pelvis, Skin, Tendons, Subcutaneous fat.
D	3/8 circle	Cardiovascular, Aponeurosis, Nerves, Dura-mater, Fascia, Periostium, Pleura, Urogenital tract, Myocardium, Skin, Eye, Muscles.
V	1/4 circle	Ophthalmology
J	Hook	Urology, laparoscopy
G	Straight	Nasal cavity, Gastro-intestinal tract, Oral cavity, Skin, Tendons, Larynx, Vessels.
K	Ski	Laparoscopy, Skin. Allows easy removal of fragments and easy needle holder manipulation.
A	Asymmetric	Plastic, Aesthetics, Robotics, Endoscopy.

Needle/thread combinations

SINGLE ARMED

Single stitch
Continuous suture



DOUBLE ARMED




Mainly for anastomosis of blood vessels
or other anatomical ducts

LOOP

Abdominal wall and tendons

LIGATURES

Surgical Sutures - Needle range

HR 	HR 	HR 
½ circle, round body needle	½ circle, round body needle	½ circle, round body needle
HR 8 		
HR 10 	HR 35 	HR 50 
HR 13 	HR 37 	HR 60 
HR 15 	HR 40 	HR 65 
HR 16 	HR 43 	HR 75 
HR 17 	HR 45 	
HR 18 	HR 48 	
HR 22 		
HR 25 		
HR 26 		
HR 26 Black 		
HR 27 		
HR 30 		

Surgical Sutures - Needle range

UK03
20



Surgical Sutures - Overview

HR s

½ circle, round body needle
- slim

HR s 18



HR s 22 Black



HR s 27 Black



HRX

½ circle, round body needle
- extra strong

HRX 22



HRX 26



HRX 27



HRX 30



HRX 35



HRX 36



HRX 37



HRX

½ circle, round body needle
- extra strong

HRX 40



HRX 45




HRX 48



HRX 50



HRX 70



Surgical Sutures - Needle range


HRT

½ circle, round body
needle, tapercut




HRT

½ circle, round body
needle, tapercut




HRTX

½ circle, round body
needle - extra strong, tapercut



HRN

½ circle, round body
needle - blunt point



Surgical Sutures - Needle range

UK03
22



Surgical Sutures - Overview

HRF 	DR 	DRT 
<p>1/2 circle, flattened body needle, taper point</p>	<p>3/8 circle, round body needle</p>	<p>3/8 circle, round body needle, tapercut</p>
<p>HRF 13 </p>	<p>DR 6 </p>	<p>DRT 10 </p>
<p>HRF 17 </p>	<p>DR 8 </p>	<p>DRT 10 Black </p>
<p>HRF 20 </p>	<p>DR 10 </p>	<p>DRT 11 </p>
<p>HRF 22 </p>	<p>DR 12 </p>	<p>DRT 13 </p>
<p>HRF 26 </p>	<p>DR 13 </p>	<p>DRT 17 </p>
<p>HRF 31 </p>	<p>DR 15 </p>	<p>DRN </p>
	<p>DR 17 </p>	<p>3/8 circle, round body needle, blunt point</p> <p>DRN 60 </p>
	<p>DR 18 </p>	<p>FR </p>
	<p>DR 21 </p>	<p>5/8 circle, round body needle</p> <p>FR 27 </p>
	<p>DR 30 </p>	<p>FR 38 </p>
	<p>DR 31 </p>	

Surgical Sutures - Needle range

GR ○	JRT ☺	HS ▼
<p>Straight, round body needle</p> <p>GR 20</p>  <p>GR 45</p>  <p>GR 60</p> 	<p>Hook shaped, round body needle, tapercut</p> <p>JRT 30</p> 	<p>½ circle, reverse cutting needle</p> <p>HS 12</p>  <p>HS 15</p>  <p>HS 18</p>  <p>HS 21</p>  <p>HS 22</p>  <p>HS 25</p>  <p>HS 26</p>  <p>HS 27</p>  <p>HS 30</p>  <p>HS 36</p> 
<p>ART ☺</p> <p>Asymmetric, round body needle, tapercut</p> <p>ART 30</p>  <p>ART 30 Black</p> 	<p>J-RX ○</p> <p>Hook shaped, round body needle</p> <p>J-RX 22</p>  <p>J-RX 30</p> 	
<p>KR ○</p> <p>Ski needle, round body</p> <p>KR 27</p> 		

Surgical Sutures - Needle range


UK03
24



Surgical Sutures - Overview


HS ▼

½ circle, reverse cutting needle





DS ▼

3/8 circle, reverse cutting needle



DS ▼

3/8 circle, reverse cutting needle



Surgical Sutures - Needle range

UK03
25

DSC ▲	DSM ▼	GS ▼
<p>3/8 circle, conventional cutting needle</p> <p>DSC 15 </p> <p>DSC 16 </p> <p>DSC 22 </p> <p>DSC 25 </p> <p>DSC 26 </p> <p>DSC 40 </p> <p>DSC 60 </p> <p>DSC 90 </p>	<p>3/8 circle, reverse cutting needle, micropoint</p> <p>DSM 6.5 </p> <p>DSM 12 </p> <p>DSM 15 </p> <p>DSM 19 </p> <p>DSM 25 </p>	<p>Straight, reverse cutting needle</p> <p>GS 13 </p> <p>GS 50 </p> <p>GS 60 </p> <p>GS 65 </p>
<p>DSMC ▲</p>		
<p>3/8 circle, conventional cutting needle, micropoint</p> <p>DSMC 16 </p> <p>DSMC 19 </p> <p>DSMC 20 </p>		<p>GS 75 </p>



Surgical Sutures - Needle range

Microsurgery

UK03
26



Surgical Sutures - Overview

DR

3/8 circle, round body needle

DR 3.75



DR 5



DR 5 Black



DR 5.1



DR 6



DR 6.5



DSM

3/8 circle, reverse cutting needle, micropoint

DSM 6.5



DSM 8



VSM

1/4 circle, reverse cutting needle, micropoint

VSM 8



HLZ

1/2 circle, spatulated

HLZ 8



DLZ

3/8 circle, spatulated

DLZ 6.2



DLZ 6.4



VLZ

1/4 circle, spatulated


VLZ 8



GLZ

Straight, spatulated

GLZ 16



Absorbable

Chirasorb rapid

Chirlac rapid

Monolac

Chirasorb

Chirasorb Plus

Absorbable

Chirlac

Polydox

Absorbable

Absorbable






ole Ab-
e Absorb
sorbab
able Ab
ole Abs
Absorb
orbable
able Ab
e Absor
sorbab
orbable
ble Abs
Absorb
sorbab
able Ab
e Absor
Absorb
orbable
ole Abs
e Absor
sorbab
able Ab
ole Abs
Absorb
orbable
able Ab
e Absor
Absorb
orbable
ble Abs



**Chirasorb
rapid**

Chirasorb

Chirasorb
rapid



CHIRASORB rapid

CHIRASORB rapid is a short-term absorbable, coated, synthetic, multifilament suture made of 90% glycolide and 10% L-lactide. The substances contained in the coating and suture are nonantigenic and noncollagenous. Sterile Chirasorb rapid is available in violet and natural.

Absorption begins as a loss of tensile strength without appreciable loss of mass. The product elicits a minimal acute inflammatory reaction in tissues.

Progressive loss of tensile strength and absorption of Chirasorb rapid occurs by hydrolysis, when the polymer is degraded to glycolic and lactic acids, which are subsequently absorbed and metabolised by the body.



UK03
31

Absorbable - Chirasorb rapid

Indication: Skin closure (e.g. intracutaneous, subcutaneous), Paediatric surgery, Episiotomies, Circumcision, Closure of oral mucosa and Ophthalmic surgery for conjunctival sutures.

Contraindication: Chirasorb rapid should not be used in tissues subject to expansion, stretching or distension, or which might require long-term mechanical support.

Sizes: 6/0 to 2 USP

Structure: Multifilament

Absorption method: Hydrolysis

Strength Retention Profile: SHORT-TERM HEALING (12 days)

1 Day	100%
5 Days	50%
10-14 Days	0%

Mass Absorption profile: 42 days

CHIRASORB rapid violet ligatures

UK03
32



Absorbable - Chirasorb rapid

Pre-cut threads

EP	USP	Length	Cat. No
1	5/0	3 x 45 cm	LR 0012
1.5	4/0	3 x 45 cm	LR 0013
2	3/0	3 x 45 cm	LR 0014
3	2/0	3 x 45 cm	LR 0015
3.5	0	3 x 45 cm	LR 0016
4	1	3 x 45 cm	LR 0017
5	2	3 x 45 cm	LR 0018
1	5/0	6 x 45 cm	LR 0042
1.5	4/0	6 x 45 cm	LR 0043
2	3/0	6 x 45 cm	LR 0044
3	2/0	6 x 45 cm	LR 0045
3.5	0	6 x 45 cm	LR 0046
4	1	6 x 45 cm	LR 0047
1	5/0	2 x 75 cm	LR 0022
1.5	4/0	2 x 75 cm	LR 0023
2	3/0	2 x 75 cm	LR 0024
3	2/0	2 x 75 cm	LR 0025

Pre-cut threads


EP	USP	Length	Cat. No
3.5	0	2 x 75 cm	LR 0026
4	1	2 x 75 cm	LR 0027
5	2	2 x 75 cm	LR 0028
2	3/0	1 x 150 cm	LR 0084
3	2/0	1 x 150 cm	LR 0085
3.5	0	1 x 150 cm	LR 0086

Chirapack - 20 pcs./box

EP	USP	Length	Cat. No
1.5	4/0	1 x 300 cm	LR 0003
2	3/0	1 x 300 cm	LR 0004
3	2/0	1 x 300 cm	LR 0005
3.5	0	1 x 300 cm	LR 0006
4	1	1 x 300 cm	LR 0007
5	2	1 x 300 cm	LR 0008

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
	1	5/0	1 x 75 cm	Violet	LR 0254	
	1.5	4/0	1 x 75 cm	Violet	LR 0255	
	1.5	4/0	1 x 75 cm	Violet	LR 0256	
	2	3/0	1 x 75 cm	Violet	LR 0257	
	1.5	4/0	1 x 75 cm	Violet	LR 0261	
	2	3/0	1 x 75 cm	Violet	LR 0262	
	3	2/0	1 x 75 cm	Violet	LR 0263	
	2	3/0	1 x 75 cm	Violet	LR 0266	
	3	2/0	1 x 75 cm	Violet	LR 0267	
	3.5	0	1 x 75 cm	Violet	LR 0268	
	3.5	0	1 x 75 cm	Violet	LR 0269	
	2	3/0	1 x 75 cm	Violet	LR 0270	
	3	2/0	1 x 75 cm	Violet	LR 0271	
	3.5	0	1 x 90 cm	Violet	LR 0372	

CHIRASORB rapid

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
HR 37 	3	2/0	1 x 45 cm	Violet	LR 0273	
HR 40 	3	2/0	1 x 75 cm	Violet	LR 0275	
	3.5	0	1 x 75 cm	Violet	LR 0276	
	4	1	1 x 75 cm	Violet	LR 0277	
HR 48 	3	2/0	1 x 90 cm	Violet	LR 0375	
	3.5	0	1 x 90 cm	Violet	LR 0376	
HR 60 	4	1	1 x 90 cm	Violet	LR 0295	
HR 60 	3.5	0	1 x 120 cm	Violet	LR 1281	
HRX: ½ circle, round body needle - extra strong						
HRX 27 	4	1	1 x 75 cm	Violet	LR 0432	
HRX 35 	3	2/0	1 x 90 cm	Violet	LR 0380	
	3.5	0	1 x 90 cm	Violet	LR 0381	
	4	1	1 x 90 cm	Violet	LR 0382	
HRX 40 	3.5	0	1 x 90 cm	Violet	LR 0385	
	4	1	1 x 90 cm	Violet	LR 0386	
	5	2	1 x 90 cm	Violet	LR 0387	
HRX 45 	4	1	1 x 90 cm	Violet	LR 0390	
HRX 50 	4	1	1 x 90 cm	Violet	LR 0395	
	5	2	1 x 90 cm	Violet	LR 0396	

CHIRASORB rapid

UK03
34



Absorbable - Chirascorb rapid

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRT: ½ circle, round body needle - tapercut						
	3.5	0	1 x 75 cm	Violet	LR 0593	
	3	2/0	1 x 90 cm	Violet	LR 0597	
	3.5	0	1 x 90 cm	Violet	LR 0598	
	4	1	1 x 90 cm	Violet	LR 0599	
	3	2/0	1 x 120 cm	Violet	LR 0697	
3.5	0	1 x 120 cm	Violet	LR 0698		
HRTX: ½ circle, round body needle - tapercut, extra strong						
	4	1	1 x 120 cm	Violet	LR 0729	
HS: ½ circle, reverse cutting needle						
	1.5	4/0	1 x 75 cm	Violet	LR 0423	
	2	3/0	1 x 75 cm	Violet	LR 0424	
	1.5	4/0	1 x 75 cm	Violet	LR 0300	
	2	3/0	1 x 75 cm	Violet	LR 0301	
	4	1	1 x 90 cm	Violet	LR 0316	
	4	1	1 x 75 cm	Violet	LR 0457	
	5	2	1 x 75 cm	Violet	LR 0458	
	3.5	0	1 x 75 cm	Violet	LR 0466	
	4	1	1 x 75 cm	Violet	LR 0476	
	5	2	1 x 75 cm	Violet	LR 0477	

CHIRASORB rapid

Needle	EP	USP	Length	Colour	Cat. No	Other information
DS: 3/8 circle, reverse cutting needle						
DS 12 	1 1.5	5/0 4/0	1 x 45 cm 1 x 45 cm	Violet Violet	LR 0202 LR 0203	
DS 15 	1.5 2	4/0 3/0	1 x 45 cm 1 x 45 cm	Violet Violet	LR 0206 LR 0207	
DS 19 	1 1.5 2 3	5/0 4/0 3/0 2/0	1 x 75 cm 1 x 75 cm 1 x 75 cm 1 x 75 cm	Violet Violet Violet Violet	LR 0210 LR 0211 LR 0212 LR 0213	
DS 25 	1.5 2 3	4/0 3/0 2/0	1 x 75 cm 1 x 75 cm 1 x 75 cm	Violet Violet Violet	LR 0223 LR 0224 LR 0225	
DS 30 	2 3	3/0 2/0	1 x 75 cm 1 x 75 cm	Violet Violet	LR 0234 LR 0235	

UK03
35



Absorbable - Chirasorb rapid



CHIRASORB rapid natural ligatures

UK03
36



Absorbable - Chirasorb rapid

Pre-cut threads

EP	USP	Length	Cat. No
1	5/0	3 x 45 cm	LQ 0012
1.5	4/0	3 x 45 cm	LQ 0013
2	3/0	3 x 45 cm	LQ 0014
3	2/0	3 x 45 cm	LQ 0015
3.5	0	3 x 45 cm	LQ 0016
1.5	4/0	6 x 45 cm	LQ 0043
2	3/0	6 x 45 cm	LQ 0044
3	2/0	6 x 45 cm	LQ 0045
3.5	0	6 x 45 cm	LQ 0046

Chirapack - 20 pcs./box

EP	USP	Length	Cat. No
1.5	4/0	1 x 300 cm	LQ 0003
2	3/0	1 x 300 cm	LQ 0004
3	2/0	1 x 300 cm	LQ 0005
3.5	0	1 x 300 cm	LQ 0006

CHIRASORB rapid

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
 HR 13	1	5/0	1 x 45 cm	Natural	LQ 0252	
	1.5	4/0	1 x 45 cm	Natural	LQ 0253	
 HR 18	1	5/0	1 x 75 cm	Natural	LQ 0255	
	1.5	4/0	1 x 75 cm	Natural	LQ 0256	
	2	3/0	1 x 75 cm	Natural	LQ 0257	
 HR 22	1.5	4/0	1 x 75 cm	Natural	LQ 0261	
	2	3/0	1 x 75 cm	Natural	LQ 0262	
	3	2/0	1 x 75 cm	Natural	LQ 0263	
	3.5	0	1 x 75 cm	Natural	LQ 0264	
 HR 27	2	3/0	1 x 75 cm	Natural	LQ 0266	
 HR 40	3.5	0	1 x 75 cm	Natural	LQ 0276	
	3	2/0	1 x 90 cm	Natural	LQ 0375	
	3.5	0	1 x 90 cm	Natural	LQ 0376	
 HR 50	3.5	0	1 x 90 cm	Natural	LQ 1392	

CHIRASORB rapid

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRX: ½ circle, round body needle - extra strong						
 HRX 36	3	2/0	1 x 90 cm	Natural	LQ 0383	
 HRX 40	3.5	0	1 x 90 cm	Natural	LQ 0385	
HRT: ½ circle, round body needle - tapercut						
 HRT 37	4	1	1 x 90 cm	Natural	LQ 0599	
 HRT 37	3	2/0	1 x 120 cm	Natural	LQ 0697	
HRF: ½ circle, flattened body needle, taper point						
 HRF 17	1.5	4/0	1 x 75 cm	Natural	LQ 9256	
 HRF 22	1.5	4/0	1 x 75 cm	Natural	LQ 9261	
 HRF 31	3	2/0	1 x 75 cm	Natural	LQ 9269	
HS: ½ circle, reverse cutting needle						
 HS 21	2	3/0	1 x 75 cm	Natural	LQ 0301	

UK03
37



Absorbable - Chirasorb rapid

CHIRASORB rapid

UK03
38



Absorbable - Chirisorb rapid

Needle	EP	USP	Length	Colour	Cat. No	Other information
DS: 3/8 circle, reverse cutting needle						
DS 15 	1.5	4/0	1 x 75 cm	Natural	LQ 1216	
	2	3/0	1 x 75 cm	Natural	LQ 1217	
DS 19 	1.5	4/0	1 x 45 cm	Natural	LQ 0216	
	1	5/0	1 x 75 cm	Natural	LQ 0210	
	1.5	4/0	1 x 75 cm	Natural	LQ 1226	
	2	3/0	1 x 75 cm	Natural	LQ 1227	
DS 25 	2	3/0	1 x 75 cm	Natural	LQ 0224	
	3	2/0	1 x 75 cm	Natural	LQ 0225	
DS 30 	2	3/0	1 x 75 cm	Natural	LQ 0234	
	3	2/0	1 x 75 cm	Natural	LQ 0235	
	3	2/0	1 x 90 cm	Natural	LQ 0239	
DSC: 3/8 circle, conventional cutting needle						
DSC 22 	2	3/0	1 x 45 cm	Natural	LQ 8227	
DSM: 3/8 circle, reverse cutting needle - micropoint						
DSM 12 	1	5/0	1 x 45 cm	Natural	LQ 0202	Advanced cutting needle
	1.5	4/0	1 x 45 cm	Natural	LQ 1203	
DSM 15 	1	5/0	1 x 45 cm	Natural	LQ 0206	
	1.5	4/0	1 x 45 cm	Natural	LQ 0207	
DSM 19 	1.5	4/0	1 x 75 cm	Natural	LQ 0212	
	2	3/0	1 x 75 cm	Natural	LQ 0213	
DSM 25 	2	3/0	1 x 75 cm	Natural	LQ 0221	
GS: Straight, reverse cutting needle						
GS 60 	2	3/0	1 x 75 cm	Natural	LQ 0236	
	3	2/0	1 x 75 cm	Natural	LQ 0237	




**Chirlac
rapid**

Chirlac rapid

Chirlac

Chirlac rapid



CHIRLAC rapid

CHIRLAC rapid is a short-term absorbable, coated, synthetic, multifilament suture made of low molecular weight polyglycolic acid. The substances contained in the coating and suture are nonantigenic and noncollagenous. Chirlac rapid is available in violet and natural.

Progressive loss of tensile strength and absorption of Chirlac rapid occurs by hydrolysis, where the polymer degrades to glycolic acid, which is subsequently absorbed and metabolised by the body.

Absorption begins as a loss of tensile strength without appreciable loss of mass. The product elicits a minimal acute inflammatory reaction in tissues.



UK03
41



Absorbable - Chirlac rapid

Indication: Gynaecology, Ophthalmic (e.g. conjunctiva), Odontology (e.g. oral mucosa), Paediatric and Skin Closure (e.g. intracutaneous, subcutaneous).

Contraindication: Chirlac rapid should not be used in tissues subject to expansion, stretching or distension, or which might require long-term mechanical support longer than seven days. The suture is not for use in ligation, cardiovascular or neurological tissues.

Sizes: 6/0 to 4 USP

Structure: Multifilament

Absorption method: Hydrolysis

Strength Retention Profile: SHORT-TERM HEALING (12 days)

1 Day	100%
5 Days	60%
10-14 Days	0%

Mass Absorption profile: 42 days

CHIRLAC rapid violet ligatures

UK03
42



Pre-cut threads

EP	USP	Length	Cat. No
1	5/0	3 x 45 cm	PR 0012
1.5	4/0	3 x 45 cm	PR 0013
2	3/0	3 x 45 cm	PR 0014
3	2/0	3 x 45 cm	PR 0015
3.5	0	3 x 45 cm	PR 0016
4	1	3 x 45 cm	PR 0017
5	2	3 x 45 cm	PR 0018
1	5/0	6 x 45 cm	PR 0042
1.5	4/0	6 x 45 cm	PR 0043
2	3/0	6 x 45 cm	PR 0044
3	2/0	6 x 45 cm	PR 0045
3.5	0	6 x 45 cm	PR 0046
4	1	6 x 45 cm	PR 0047
5	2	6 x 45 cm	PR 0048
1	5/0	2 x 75 cm	PR 0022
1.5	4/0	2 x 75 cm	PR 0023
2	3/0	2 x 75 cm	PR 0024
3	2/0	2 x 75 cm	PR 0025
3.5	0	2 x 75 cm	PR 0026
4	1	2 x 75 cm	PR 0027
5	2	2 x 75 cm	PR 0028
3.5	0	2 x 90 cm	PR 0097
2	3/0	1 x 150 cm	PR 0084
3	2/0	1 x 150 cm	PR 0085
3.5	0	1 x 150 cm	PR 0086

Chirapack - 20 pcs./box

EP	USP	Length	Cat. No
1.5	4/0	1 x 300 cm	PR 0003
2	3/0	1 x 300 cm	PR 0004
3	2/0	1 x 300 cm	PR 0005
3.5	0	1 x 300 cm	PR 0006
4	1	1 x 300 cm	PR 0007
5	2	1 x 300 cm	PR 0008



Absorbable - Chirlac rapid



CHIRLAC rapid

UK03
43

Absorbable - Chirlac rapid

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
HR 13 	1	5/0	1 x 75 cm	Violet	PR 0254	
	1.5	4/0	1 x 75 cm	Violet	PR 0255	
HR 17 	1.5	4/0	1 x 75 cm	Violet	PR 0258	
HR 18 	1	5/0	1 x 75 cm	Violet	PR 0259	
	1.5	4/0	1 x 75 cm	Violet	PR 0256	
	2	3/0	1 x 75 cm	Violet	PR 0257	
HR 22 	1.5	4/0	1 x 75 cm	Violet	PR 0261	
	2	3/0	1 x 75 cm	Violet	PR 0262	
	3	2/0	1 x 75 cm	Violet	PR 0263	
HR 27 	1.5	4/0	1 x 75 cm	Violet	PR 0265	
	2	3/0	1 x 75 cm	Violet	PR 0266	
	3	2/0	1 x 75 cm	Violet	PR 0267	
	3.5	0	1 x 75 cm	Violet	PR 0268	
HR 30 	3.5	0	1 x 90 cm	Violet	PR 0168	
HR 35 	2	3/0	1 x 75 cm	Violet	PR 0270	
	3	2/0	1 x 75 cm	Violet	PR 0271	
	3.5	0	1 x 90 cm	Violet	PR 0372	
HR 40 	3	2/0	1 x 75 cm	Violet	PR 0275	
	3.5	0	1 x 75 cm	Violet	PR 0276	
	4	1	1 x 75 cm	Violet	PR 0277	
	3	2/0	1 x 90 cm	Violet	PR 0375	
HR 48 	3.5	0	1 x 90 cm	Violet	PR 0376	
	4	1	1 x 90 cm	Violet	PR 0377	
	3	2/0	1 x 120 cm	Violet	PR 1375	
HR 48 	3.5	0	1 x 90 cm	Violet	PR 0294	
HR 50 	3.5	0	1 x 120 cm	Violet	PR 1293	
HR 60 	3.5	0	1 x 120 cm	Violet	PR 1281	

CHIRLAC rapid

UK03
44



Absorbable - Chirlac rapid

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRX: ½ circle, round body needle - extra strong						
HRX 27 	4	1	1 x 75 cm	Violet	PR 0432	
HRX 35 	3 3.5 4	2/0 0 1	1 x 90 cm 1 x 90 cm 1 x 90 cm	Violet Violet Violet	PR 0380 PR 0381 PR 0382	
HRX 40 	3.5 4 5	0 1 2	1 x 90 cm 1 x 90 cm 1 x 90 cm	Violet Violet Violet	PR 0385 PR 0386 PR 0387	
	4	1	1 x 120 cm	Violet	PR 0384	
HRX 45 	4	1	1 x 90 cm	Violet	PR 0390	
HRX 50 	4 5	1 2	1 x 90 cm 1 x 90 cm	Violet Violet	PR 0395 PR 0396	
HRT: ½ circle, round body needle - tapercut						
HRT 25 	3	2/0	1 x 75 cm	Violet	PR 0584	
	3.5	0	1 x 75 cm	Violet	PR 0593	
HRT 37 	3 3.5 4	2/0 0 1	1 x 90 cm 1 x 90 cm 1 x 90 cm	Violet Violet Violet	PR 0597 PR 0598 PR 0599	
	3 3.5	2/0 0	1 x 120 cm 1 x 120 cm	Violet Violet	PR 0697 PR 0698	

CHIRLAC rapid

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRTX: ½ circle, round body needle - tapercut, extra strong						

UK03
45

HRTX 50



4

1

1 x 120 cm

Violet

PR 0729

HS: ½ circle, reverse cutting needle

HS 18



1.5
2

4/0
3/0

1 x 75 cm
1 x 75 cm

Violet
Violet

PR 0423
PR 0424

HS 21



1.5
2

4/0
3/0

1 x 75 cm
1 x 75 cm

Violet
Violet

PR 0300
PR 0301

HS 37



4

1

1 x 90 cm

Violet

PR 0316

HS 40



4
5

1
2

1 x 75 cm
1 x 75 cm

Violet
Violet

PR 0457
PR 0458

HS 48



3.5

0

1 x 75 cm

Violet

PR 0466

HS 55



4
5

1
2

1 x 75 cm
1 x 75 cm

Violet
Violet

PR 0476
PR 0477



Absorbable - Chirlac rapid

CHIRLAC rapid

UK03
46

Absorbable - Chirlac rapid

Needle	EP	USP	Length	Colour	Cat. No	Other information
DS: 3/8 circle, reverse cutting needle						
DS 12 	1 1.5	5/0 4/0	1 x 45 cm 1 x 45 cm	Violet Violet	PR 0202 PR 0203	
DS 15 	1.5 2	4/0 3/0	1 x 45 cm 1 x 45 cm	Violet Violet	PR 0206 PR 0207	
	1	5/0	1 x 75 cm	Violet	PR 0228	
DS 19 	1 1.5 2 3	5/0 4/0 3/0 2/0	1 x 75 cm 1 x 75 cm 1 x 75 cm 1 x 75 cm	Violet Violet Violet Violet	PR 0210 PR 0211 PR 0212 PR 0213	
DS 25 	1.5 2 3	4/0 3/0 2/0	1 x 75 cm 1 x 75 cm 1 x 75 cm	Violet Violet Violet	PR 0223 PR 0224 PR 0225	
DS 30 	2 3	3/0 2/0	1 x 75 cm 1 x 75 cm	Violet Violet	PR 0234 PR 0235	
DSM: 3/8 circle, reverse cutting needle - micropoint						
DSM 19 	1	5/0	1 x 75 cm	Violet	PR 0185	Advanced cutting needle



CHIRLAC rapid natural ligatures

Pre-cut threads

EP	USP	Length	Cat. No
1	5/0	3 x 45 cm	PQ 0012
1.5	4/0	3 x 45 cm	PQ 0013
2	3/0	3 x 45 cm	PQ 0014
3	2/0	3 x 45 cm	PQ 0015
3.5	0	3 x 45 cm	PQ 0016
1.5	4/0	6 x 45 cm	PQ 0043
2	3/0	6 x 45 cm	PQ 0044
3	2/0	6 x 45 cm	PQ 0045
3.5	0	6 x 45 cm	PQ 0046



UK03
47

Chirapack - 20 pcs./box

EP	USP	Length	Cat. No
1.5	4/0	1 x 300 cm	PQ 0003
2	3/0	1 x 300 cm	PQ 0004
3	2/0	1 x 300 cm	PQ 0005
3.5	0	1 x 300 cm	PQ 0006
4	1	1 x 300 cm	PQ 0007



Absorbable - Chirlac rapid

CHIRLAC rapid

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
 HR 13	1	5/0	1 x 45 cm	Natural	PQ 0252	
	1.5	4/0	1 x 45 cm	Natural	PQ 0253	
 HR 18	1	5/0	1 x 75 cm	Natural	PQ 0255	
	1.5	4/0	1 x 75 cm	Natural	PQ 0256	
	2	3/0	1 x 75 cm	Natural	PQ 0257	
 HR 22	1.5	4/0	1 x 75 cm	Natural	PQ 0261	
	2	3/0	1 x 75 cm	Natural	PQ 0262	
	3	2/0	1 x 75 cm	Natural	PQ 0263	
	3.5	0	1 x 75 cm	Natural	PQ 0264	
 HR 27	2	3/0	1 x 75 cm	Natural	PQ 0266	
 HR 40	3.5	0	1 x 75 cm	Natural	PQ 0276	
 HR 50	3	2/0	1 x 90 cm	Natural	PQ 0375	
	3.5	0	1 x 90 cm	Natural	PQ 0376	
 HR 50	3.5	0	1 x 90 cm	Natural	PQ 1392	

CHIRLAC rapid

UK03
48



Absorbable - Chirlac rapid

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRX: ½ circle, round body needle - extra strong						
HRX 36 	3	2/0	1 x 90 cm	Natural	PQ 0383	
HRX 40 	3.5	0	1 x 90 cm	Natural	PQ 0385	
HRT: ½ circle, round body needle - tapercut						
HRT 37 	3	2/0	1 x 90 cm	Natural	PQ 0597	
	3 3.5	2/0 0	1 x 120 cm 1 x 120 cm	Natural Natural	PQ 0697 PQ 0698	
HRF: ½ circle, flattened body needle, taper point						
HRF 17 	1.5	4/0	1 x 75 cm	Natural	PQ 9256	
HRF 31 	3	2/0	1 x 75 cm	Natural	PQ 9269	
HS: ½ circle, reverse cutting needle						
HS 21 	1.5 2	4/0 3/0	1 x 75 cm 1 x 75 cm	Natural Natural	PQ 0300 PQ 0301	




CHIRLAC rapid

Needle	EP	USP	Length	Colour	Cat. No	Other information
DS: 3/8 circle, reverse cutting needle						
DS 12 	1	5/0	1 x 45 cm	Natural	PQ 0203	
DS 15 	1.5 2	4/0 3/0	1 x 75 cm 1 x 75 cm	Natural Natural	PQ 1216 PQ 1217	
DS 19 	1 1.5 2	5/0 4/0 3/0	1 x 45 cm 1 x 45 cm 1 x 45 cm	Natural Natural Natural	PQ 0215 PQ 0216 PQ 0217	
	1.5 2	4/0 3/0	1 x 75 cm 1 x 75 cm	Natural Natural	PQ 1226 PQ 1227	
DS 25 	2 3	3/0 2/0	1 x 75 cm 1 x 75 cm	Natural Natural	PQ 0224 PQ 0225	
DS 30 	2 3	3/0 2/0	1 x 75 cm 1 x 75 cm	Natural Natural	PQ 0234 PQ 0235	
	3	2/0	1 x 90 cm	Natural	PQ 0239	
DSC: 3/8 circle, conventional cutting needle						
DSC 22 	2	3/0	1 x 45 cm	Natural	PQ 8227	

UK03
49



Absorbable - Chirlac rapid



CHIRLAC rapid

UK03
50



Absorbable - Chirlac rapid

Needle	EP	USP	Length	Colour	Cat. No	Other information
DSM: 3/8 circle, reverse cutting needle - micropoint						
DSM 12	1 1.5	5/0 4/0	1 x 45 cm 1 x 45 cm	Natural Natural	PQ 0202 PQ 1203	Advanced cutting needle
	1	5/0	1 x 75 cm	Natural	PQ 0198	
DSM 15	1 1.5	5/0 4/0	1 x 45 cm 1 x 45 cm	Natural Natural	PQ 0206 PQ 0207	Advanced cutting needle
						
DSM 19	1.5 2	4/0 3/0	1 x 75 cm 1 x 75 cm	Natural Natural	PQ 0212 PQ 0213	Advanced cutting needle
						
DSM 25	2	3/0	1 x 75 cm	Natural	PQ 0221	Advanced cutting needle
						
GS: Straight, reverse cutting needle						
GS 60	3	2/0	1 x 75 cm	Natural	PQ 0237	
						

Wono



Monolac

Monolac

Monolac

Monolac

Wono



MONOLAC

MONOLAC is a mid-term absorbable, synthetic, monofilament suture made of glycolide and ε-Caprolactone, available in natural colour for the closure of superficial wounds and in violet for all other indications, when a mid-term monofilament is needed.

During absorption the material elicited only a very slight inflammation.

The material is degraded into adipic acid, which is completely metabolised in the body. The process begins with a loss of tensile strength, followed by a loss of mass at a later stage.



UK03
53



Absorbable - Monolac

Indication: Soft tissue in all specialities, Skin Closure (e.g. intracutaneous, subcutaneous), Plastic, Gynaecology, Urology and Gastrointestinal surgery.

Contraindication: Monolac should not be used in tissues subject to expansion, stretching or distension, or which might require long-term mechanical support.

Sizes: 6/0 to 2 USP

Structure: Monofilament

Absorption method: Hydrolysis

Strength Retention

Profile:	MID-TERM HEALING (30 days)	
	1 Day	100%
	7 Days	60% violet/50% natural
	14 Days	30% violet/20% natural

Mass Absorption profile: 90 - 120 days

MONOLAC violet ligatures

UK03
54

Pre-cut threads

EP	USP	Length	Cat. No
3	2/0	1 x 150 cm	MV 0085
3.5	0	1 x 150 cm	MV 0086



MONOLAC

Absorbable - Monolac

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
 HR 13	0.7	6/0	1 x 75 cm	Violet	MV 0251	
	1	5/0	1 x 75 cm	Violet	MV 0252	
	1.5	4/0	1 x 75 cm	Violet	MV 0253	
 HR 17	1.5	4/0	1 x 75 cm	Violet	MV 0254	
 HR 18	1	5/0	1 x 75 cm	Violet	MV 0256	
	1.5	4/0	1 x 75 cm	Violet	MV 0257	
	2	3/0	1 x 75 cm	Violet	MV 0258	
 HR 22	1.5	4/0	1 x 75 cm	Violet	MV 0261	
	2	3/0	1 x 75 cm	Violet	MV 0262	
	3	2/0	1 x 75 cm	Violet	MV 0263	
 2 x HR 22	3	2/0	1 x 45 cm	Violet	MV 0963	
	1.5	4/0	1 x 75 cm	Violet	MV 0965	
	2	3/0	1 x 90 cm	Violet	MV 0964	
 HR 26	2	3/0	1 x 75 cm	Violet	MV 0166	
	3	2/0	1 x 75 cm	Violet	MV 0167	
 HR 27	1	5/0	1 x 75 cm	Violet	MV 0264	
	1.5	4/0	1 x 75 cm	Violet	MV 0265	
	2	3/0	1 x 75 cm	Violet	MV 0266	
	3	2/0	1 x 75 cm	Violet	MV 0267	
	2	3/0	1 x 90 cm	Violet	MV 0366	
 2 x HR 27	2	3/0	1 x 75 cm	Violet	MV 2266	
	1.5	4/0	1 x 90 cm	Violet	MV 2275	
 HR 30	2	3/0	1 x 75 cm	Violet	MV 0269	
	3	2/0	1 x 75 cm	Violet	MV 0273	
	3.5	0	1 x 75 cm	Violet	MV 0274	

MONOLAC

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
2 x HR 30 	3	2/0	1 x 90 cm	Violet	MV 1373	
HR 35 	3	2/0	1 x 75 cm	Violet	MV 0271	
	3.5	0	1 x 75 cm	Violet	MV 0272	
	2	3/0	1 x 90 cm	Violet	MV 0373	
	4	1	1 x 90 cm	Violet	MV 0376	
HR 40 	3.5	0	1 x 75 cm	Violet	MV 0276	
	4	1	1 x 75 cm	Violet	MV 0277	
	2	3/0	1 x 90 cm	Violet	MV 0377	
	4	1	1 x 90 cm	Violet	MV 0397	
	5	2	1 x 90 cm	Violet	MV 0398	
HR s: ½ circle, round body needle - slim						
HR s 18 	1.5	4/0	1 x 75 cm	Violet	MV 1256	Slim needle
HR s 22 	1.5	4/0	1 x 75 cm	Violet	MV 1261	Slim needle
	2	3/0	1 x 75 cm	Violet	MV 1262	
HR s 22 Black 	2	3/0	1 x 75 cm	Violet	MV 2262	Slim, black needle
HR s 27 	3	2/0	1 x 75 cm	Violet	MV 1267	Slim needle



MONOLAC

UK03
56



Absorbable - Monolac

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRX: ½ circle, round body needle - extra strong						
HRX 35	3	2/0	1 x 75 cm	Violet	MV 0280	
	3	2/0	1 x 90 cm	Violet	MV 0380	
	3.5	0	1 x 90 cm	Violet	MV 0381	
HRX 40	4	1	1 x 75 cm	Violet	MV 0286	
	4	1	1 x 90 cm	Violet	MV 0386	
	4	1	1 x 150 cm	Violet	MV 4530	Loop
HRX 45	4	1	1 x 90 cm	Violet	MV 0390	
	3.5	0	1 x 90 cm	Violet	MV 0394	
	4	1	1 x 90 cm	Violet	MV 0395	
	4	1	1 x 150 cm	Violet	MV 4540	Loop
HRT: ½ circle, round body needle - tapercut						
HRT 17	1.5	4/0	1 x 75 cm	Violet	MV 0561	
	2	3/0	1 x 75 cm	Violet	MV 0581	
	3	2/0	1 x 75 cm	Violet	MV 0582	
HRT 37	4	1	1 x 90 cm	Violet	MV 1594	
HRN: ½ circle, round body needle - blunt point						
HRN 50	4	1	1 x 90 cm	Violet	MV 0686	

MONOLAC

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRF: ½ circle, flattened body needle, taper point						
HRF 17 	1.5	4/0	1 x 75 cm	Violet	MV 9257	
HRF 22 	2	3/0	1 x 75 cm	Violet	MV 9262	
DR: 3/8 circle, round body needle						
DR 8 	0.7	6/0	1 x 45 cm	Violet	MV 0703	
DR 20 	2	3/0	1 x 75 cm	Violet	MV 0726	
DRT: 3/8 circle, round body needle - tapercut						
DRT 13 	0.7	6/0	1 x 45 cm	Violet	MV 0544	
FR: 5/8 circle, round body needle						
FR 27 	2 3	3/0 2/0	1 x 75 cm 1 x 75 cm	Violet Violet	MV 0609 MV 0610	
GR: Straight, round body needle						
GR 20 	1	5/0	1 x 45 cm	Violet	MV 0405	
	1.5	4/0	1 x 75 cm	Violet	MV 0406	
	2	3/0	1 x 75 cm	Violet	MV 0407	
HS: ½ circle, reverse cutting needle						
HS 12 	0.7 1	6/0 5/0	1 x 45 cm 1 x 45 cm	Violet Violet	MV 0301 MV 0302	
HS 18 	1 1.5	5/0 4/0	1 x 75 cm 1 x 75 cm	Violet Violet	MV 0312 MV 0313	



MONOLAC

UK03
58



Absorbable - Monolac

Needle	EP	USP	Length	Colour	Cat. No	Other information
DS: 3/8 circle, reverse cutting needle						
DS 12 	0.7	6/0	1 x 45 cm	Violet	MV 0201	
	1	5/0	1 x 45 cm	Violet	MV 0202	
DS 13 	1	5/0	1 x 75 cm	Violet	MV 0203	
	1	5/0	1 x 45 cm	Violet	MV 0103	
DS 15 	1	5/0	1 x 45 cm	Violet	MV 0205	
	1.5	4/0	1 x 45 cm	Violet	MV 0206	
DS 16 	1.5	4/0	1 x 45 cm	Violet	MV 1206	
	1	5/0	1 x 45 cm	Violet	MV 0214	
DS 19 	1.5	4/0	1 x 45 cm	Violet	MV 0215	
	2	3/0	1 x 45 cm	Violet	MV 0216	
	1	5/0	1 x 75 cm	Violet	MV 0210	
DS 25 	1.5	4/0	1 x 75 cm	Violet	MV 0211	
	2	3/0	1 x 75 cm	Violet	MV 0212	
	2	3/0	1 x 45 cm	Violet	MV 0122	
DS 30 	2	3/0	1 x 75 cm	Violet	MV 0223	
	3	2/0	1 x 75 cm	Violet	MV 0224	
DS 35 	2	3/0	1 x 90 cm	Violet	MV 1223	
	3.5	0	1 x 75 cm	Violet	MV 0227	
	3	2/0	1 x 75 cm	Violet	MV 0228	
DS 35 	4	1	1 x 75 cm	Violet	MV 1240	
	3.5	0	1 x 75 cm	Violet	MV 0229	

MONOLAC

UK03
59

Absorbable - Monolac

Needle	EP	USP	Length	Colour	Cat. No	Other information
DSM: 3/8 circle, reverse cutting needle - micropoint						
 DSM 12	0.7 1	6/0 5/0	1 x 45 cm 1 x 45 cm	Violet Violet	MV 0230 MV 0231	Advanced cutting needle
 DSM 15	1 1.5	5/0 4/0	1 x 45 cm 1 x 45 cm	Violet Violet	MV 0234 MV 0235	Advanced cutting needle
 DSM 19	1.5 2	4/0 3/0	1 x 45 cm 1 x 45 cm	Violet Violet	MV 0237 MV 0238	Advanced cutting needle
 DSM 25	2 1.5	3/0 4/0	1 x 75 cm 1 x 90 cm	Violet Violet	MV 0242 MV 0240	Advanced cutting needle
GS: Straight, reverse cutting needle						
 GS 50	2	3/0	1 x 75 cm	Violet	MV 0456	
 GS 60	2 3	3/0 2/0	1 x 75 cm 1 x 75 cm	Violet Violet	MV 0466 MV 0467	



MONOLAC

UK03
60



Absorbable - Monolac

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
HR 13 	1	5/0	1 x 75 cm	Natural	MN 0252	
HR 18 	1 1.5	5/0 4/0	1 x 75 cm 1 x 75 cm	Natural Natural	MN 0256 MN 0257	
HR 22 	2	3/0	1 x 75 cm	Natural	MN 0262	
HR 35 	3.5	0	1 x 75 cm	Natural	MN 0275	
FR: 5/8 circle, round body needle						
FR 27 	3	2/0	1 x 75 cm	Natural	MN 0610	
DS: 3/8 circle, reverse cutting needle						
DS 12 	0.7 1	6/0 5/0	1 x 45 cm 1 x 45 cm	Natural Natural	MN 0201 MN 0202	
DS 15 	0.7 1 1.5	6/0 5/0 4/0	1 x 45 cm 1 x 45 cm 1 x 45 cm	Natural Natural Natural	MN 0204 MN 0205 MN 0206	
DS 19 	1 1.5 2	5/0 4/0 3/0	1 x 75 cm 1 x 75 cm	Natural Natural	MN 0208 MN 0209	
DS 25 	0.7 1 1.5 2 3	6/0 5/0 4/0 3/0	1 x 45 cm 1 x 45 cm 1 x 45 cm 1 x 45 cm	Natural Natural Natural Natural	MN 0213 MN 0214 MN 0215 MN 0216	
DS 30 	1.5 2 3	4/0 3/0 2/0	1 x 75 cm 1 x 75 cm 1 x 75 cm	Natural Natural Natural	MN 0222 MN 0223 MN 0224	
DS 30 	2	3/0	1 x 75 cm	Natural	MN 0227	

MONOLAC



Needle	EP	USP	Length	Colour	Cat. No	Other information
DSM: 3/8 circle, reverse cutting needle - micropoint						
DSM 12 	0.7	6/0	1 x 45 cm	Natural	MN 0230	Advanced cutting needle
	1	5/0	1 x 45 cm	Natural	MN 0231	
	1.5	4/0	1 x 45 cm	Natural	MN 0232	
DSM 15 	1	5/0	1 x 45 cm	Natural	MN 0234	Advanced cutting needle
	1.5	4/0	1 x 45 cm	Natural	MN 0235	
DSM 19 	1.5	4/0	1 x 45 cm	Natural	MN 0237	Advanced cutting needle
	2	3/0	1 x 45 cm	Natural	MN 0238	
	1.5	4/0	1 x 75 cm	Natural	MN 0240	
DSM 25 	2	3/0	1 x 75 cm	Natural	MN 0248	Advanced cutting needle
	3	2/0	1 x 75 cm	Natural	MN 0249	
GS: Straight, reverse cutting needle						
GS 50  	2	3/0	1 x 75 cm	Natural	MN 0456	
	1.5	4/0	1 x 75 cm	Natural	MN 0465	
GS 60  	2	3/0	1 x 75 cm	Natural	MN 0466	



chirasorb

Chirasorb

chirasorb

chirasorb

Chirasorb



CHIRASORB

CHIRASORB is a mid-term absorbable, coated, synthetic, multifilament suture made of 90% glycolide and 10% L-lactide. The substances contained in the coating and suture are nonantigenic and noncollagenous.

Sterile Chirasorb is available in violet or natural. The product elicits a minimal acute inflammatory reaction in tissues.

Progressive loss of tensile strength and absorption of Chirasorb occurs by hydrolysis, when the polymer is degraded to glycolic and lactic acids, which are subsequently absorbed and metabolised by the body.



Indication: General surgery, Gynaecology, Ophthalmic, Skin Closure (e.g. intracutaneous, subcutaneous), Muscular layer and Ligation.

Contraindication: Chirasorb should not be used in tissues subject to expansion, stretching or distension, or which might require long-term mechanical support.

Sizes: 10/0 to 6 USP

Structure: Multifilament

Absorption method: Hydrolysis

Strength Retention Profile: MID-TERM HEALING (35 days)

1 Day	100%
14 Days	75%
21 Days	50%

Mass Absorption profile: 56 - 70 days

CHIRASORB violet ligatures

UK03
66



Absorbable - Chirasorb

Pre-cut threads

EP	USP	Length	Cat. No
1	5/0	3 x 45 cm	LV 0012
1.5	4/0	3 x 45 cm	LV 0013
2	3/0	3 x 45 cm	LV 0014
3	2/0	3 x 45 cm	LV 0015
3.5	0	3 x 45 cm	LV 0016
4	1	3 x 45 cm	LV 0017
5	2	3 x 45 cm	LV 0018
1	5/0	6 x 45 cm	LV 0042
1.5	4/0	6 x 45 cm	LV 0043
2	3/0	6 x 45 cm	LV 0044
3	2/0	6 x 45 cm	LV 0045
3.5	0	6 x 45 cm	LV 0046
4	1	6 x 45 cm	LV 0047
5	2	6 x 45 cm	LV 0048
2	3/0	10 x 45 cm	LV 0064
3	2/0	10 x 45 cm	LV 0065
3.5	0	10 x 45 cm	LV 0066
4	1	10 x 45 cm	LV 0067
5	2	10 x 45 cm	LV 0068
3	2/0	12 x 45 cm	LV 4565
5	2	12 x 45 cm	LV 4568
1	5/0	2 x 75 cm	LV 0022
1.5	4/0	2 x 75 cm	LV 0023
2	3/0	2 x 75 cm	LV 0024
3	2/0	2 x 75 cm	LV 0025
3.5	0	2 x 75 cm	LV 0026
4	1	2 x 75 cm	LV 0027
5	2	2 x 75 cm	LV 0028
1.5	4/0	5 x 75 cm	LV 0053
2	3/0	5 x 75 cm	LV 0054
3	2/0	5 x 75 cm	LV 0055



Pre-cut threads

EP	USP	Length	Cat. No
3.5	0	5 x 75 cm	LV 0056
4	1	5 x 75 cm	LV 0057
5	2	5 x 75 cm	LV 0058
2	3/0	6 x 75 cm	LV 0059
3	2/0	6 x 75 cm	LV 0051
3.5	0	6 x 75 cm	LV 0060
4	1	6 x 75 cm	LV 4557
1.5	4/0	1 x 150 cm	LV 0083
2	3/0	1 x 150 cm	LV 0084
3	2/0	1 x 150 cm	LV 0085
3.5	0	1 x 150 cm	LV 0089
4	1	1 x 150 cm	LV 0087
5	2	1 x 150 cm	LV 0090

Chirapack - 20 pcs./box

EP	USP	Length	Cat. No
2	3/0	1 x 250 cm	LV 1004
3	2/0	1 x 250 cm	LV 1005
3.5	0	1 x 250 cm	LV 1006
4	1	1 x 250 cm	LV 1007
5	2	1 x 250 cm	LV 1008
1.5	4/0	1 x 300 cm	LV 0003
2	3/0	1 x 300 cm	LV 0004
3	2/0	1 x 300 cm	LV 0005
3.5	0	1 x 300 cm	LV 0006
4	1	1 x 300 cm	LV 0007
5	2	1 x 300 cm	LV 0008
6	3 + 4	1 x 300 cm	LV 0009

CHIRASORB

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
HR 10 	0.7	6/0	1 x 45 cm	Violet	LV 0246	
	1	5/0	1 x 75 cm	Violet	LV 0248	
HR 13 	0.7	6/0	1 x 45 cm	Violet	LV 0251	
	1	5/0	1 x 45 cm	Violet	LV 0252	
	1.5	4/0	1 x 45 cm	Violet	LV 0253	
	2	3/0	1 x 45 cm	Violet	LV 0250	
	0.7	6/0	1 x 75 cm	Violet	LV 1251	
	1	5/0	1 x 75 cm	Violet	LV 0254	
	1.5	4/0	1 x 75 cm	Violet	LV 0255	
HR 16 	1	5/0	1 x 75 cm	Violet	LV 3265	
HR 18 	1.5	4/0	1 x 45 cm	Violet	LV 0258	
	2	3/0	1 x 45 cm	Violet	LV 0259	
HR 18 	0.7	6/0	1 x 75 cm	Violet	LV 0241	
	1	5/0	1 x 75 cm	Violet	LV 0242	
	1.5	4/0	1 x 75 cm	Violet	LV 0256	
	2	3/0	1 x 75 cm	Violet	LV 0257	
	3	2/0	1 x 75 cm	Violet	LV 0260	
HR 20 	1.5	4/0	1 x 75 cm	Violet	LV 1261	
	2	3/0	1 x 75 cm	Violet	LV 1259	
HR 22 	1.5	4/0	1 x 75 cm	Violet	LV 0261	
	2	3/0	1 x 75 cm	Violet	LV 0262	
	3	2/0	1 x 75 cm	Violet	LV 0263	
HR 22 	2	3/0	4 x 75 cm	Violet	LV 0322	
	3	2/0	4 x 75 cm	Violet	LV 0323	
	2	3/0	1 x 120 cm	Violet	LV 0363	



CHIRASORB

UK03
68



Absorbable - Chirasorb

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
HR 22 A 	2	3/0	6 x 45 cm	Violet	LV 0327	Take-off needle
HR 25 	5	2	1 x 75 cm	Violet	LV 0297	
HR 26 	3	2/0	1 x 45 cm	Violet	LV 0136	
	2	3/0	1 x 75 cm	Violet	LV 0535	
	3	2/0	1 x 75 cm	Violet	LV 0536	
	3.5	0	1 x 75 cm	Violet	LV 0540	
HR 26 A 	2	3/0	6 x 45 cm	Violet	LV 0347	Take-off needle
	1.5	4/0	1 x 45 cm	Violet	LV 0245	
	2	3/0	6 x 45 cm	Violet	LV 0537	
	3	2/0	6 x 45 cm	Violet	LV 0538	
HR 27 	1.5	4/0	1 x 75 cm	Violet	LV 0265	
	2	3/0	1 x 75 cm	Violet	LV 0266	
	3	2/0	1 x 75 cm	Violet	LV 0267	
	3.5	0	1 x 75 cm	Violet	LV 0268	
	4	1	1 x 75 cm	Violet	LV 1269	
		2	3/0	1 x 90 cm	Violet	LV 0366
HR 27 A 	3	2/0	6 x 45 cm	Violet	LV 0338	Take-off needle
	2	3/0	4 x 75 cm	Violet	LV 0639	Take-off needle


CHIRASORB

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
2 x HR 27 	2	3/0	1 x 90 cm	Violet	LV 2466	
	3	2/0	1 x 90 cm	Violet	LV 2467	
	2	3/0	1 x 45 cm	Violet	LV 0169	
HR 30 	2	3/0	1 x 75 cm	Violet	LV 0269	
	3	2/0	1 x 75 cm	Violet	LV 0273	
	3.5	0	1 x 75 cm	Violet	LV 0274	
	4	1	1 x 75 cm	Violet	LV 1275	
	2	3/0	1 x 90 cm	Violet	LV 0369	
	3	2/0	1 x 90 cm	Violet	LV 0364	
	3.5	0	1 x 90 cm	Violet	LV 0371	
	4	1	1 x 90 cm	Violet	LV 0372	
HR 35 	2	3/0	1 x 75 cm	Violet	LV 0270	
	3	2/0	1 x 75 cm	Violet	LV 0271	
	3.5	0	1 x 75 cm	Violet	LV 0272	
	4	1	1 x 75 cm	Violet	LV 1273	
	5	2	1 x 75 cm	Violet	LV 1272	
	2	3/0	1 x 90 cm	Violet	LV 0373	
	3	2/0	1 x 90 cm	Violet	LV 0374	
	3.5	0	1 x 90 cm	Violet	LV 0375	
	4	1	1 x 90 cm	Violet	LV 0376	
HR 37 	3.5	0	1 x 75 cm	Violet	LV 0579	
	4	1	1 x 75 cm	Violet	LV 0379	
	3.5	0	1 x 90 cm	Violet	LV 0378	



CHIRASORB

UK03
70

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
	4	1	1 x 45 cm	Violet	LV 0177	
	5	2	5 x 45 cm	Violet	LV 2250	
 <p>HR 40</p>	3	2/0	1 x 75 cm	Violet	LV 0275	
	3.5	0	1 x 75 cm	Violet	LV 0276	
	4	1	1 x 75 cm	Violet	LV 0277	
	3	2/0	1 x 90 cm	Violet	LV 1475	
	3.5	0	1 x 90 cm	Violet	LV 1476	
	4	1	1 x 90 cm	Violet	LV 1477	
	5	2	1 x 90 cm	Violet	LV 1478	
	3	2/0	1 x 120 cm	Violet	LV 1479	
 <p>HR 43</p>	5	2	1 x 75 cm	Violet	LV 0479	
 <p>HR 48</p>	4	1	1 x 75 cm	Violet	LV 1396	
	5	2	1 x 75 cm	Violet	LV 1400	
	3.5	0	1 x 75 cm	Violet	LV 1402	
	3.5	0	1 x 90 cm	Violet	LV 1403	
	4	1	1 x 90 cm	Violet	LV 1401	
	5	2	1 x 90 cm	Violet	LV 1399	
 <p>HR 48 A</p>	4	1	4 x 45 cm	Violet	LV 0793	Take-off needle





Absorbable - Chirasorb

CHIRASORB

UK03
71

Absorbable - Chirasorb

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						

HR 50 	5	2	5 x 45 cm	Violet	LV 1201	
	3.5	0	1 x 90 cm	Violet	LV 1392	
HR 60 	3 3.5	2/0 0	1 x 75 cm 1 x 75 cm	Violet Violet	LV 0278 LV 0279	
HR 65 	3 4 5	2/0 1 2	1 x 75 cm 1 x 75 cm 1 x 75 cm	Violet Violet Violet	LV 1283 LV 1285 LV 1286	
	5	2	1 x 150 cm	Violet	LV 1287	
HR 75 	5	2	1 x 75 cm	Violet	LV 0496	

HR s: ½ circle, round body needle - slim						
---	--	--	--	--	--	--

HR s 18 	1.5	4/0	1 x 75 cm	Violet	LV 1256	Slim needle
HR s 27 Black 	3	2/0	1 x 75 cm	Violet	LV 2267	Slim, black needle

CHIRASORB

UK03
72



Absorbable - Chirasorb

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRX: ½ circle, round body needle - extra strong						
HRX 22 	2	3/0	1 x 75 cm	Violet	LV 0420	
HRX 26 	2 3	3/0 2/0	1 x 75 cm 1 x 75 cm	Violet Violet	LV 0419 LV 0417	
HRX 27 	2 3 4	3/0 2/0 1	1 x 75 cm 1 x 75 cm 1 x 75 cm	Violet Violet Violet	LV 0429 LV 0430 LV 0432	
HRX 30 	2	3/0	1 x 75 cm	Violet	LV 0426	
HRX 35 	3 4	2/0 1	1 x 75 cm 1 x 75 cm	Violet Violet	LV 0280 LV 0282	
	3.5	0	1 x 90 cm	Violet	LV 0381	



CHIRASORB

UK03
73



Absorbable - Chirisorb

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRX: ½ circle, round body needle - extra strong						
 HRX 36	3.5	0	1 x 75 cm	Violet	LV 1281	
	4	1	1 x 75 cm	Violet	LV 1282	
	3	2/0	1 x 90 cm	Violet	LV 1380	
 HRX 40	3.5	0	1 x 75 cm	Violet	LV 0285	
	4	1	1 x 75 cm	Violet	LV 0286	
	3.5	0	1 x 90 cm	Violet	LV 0385	
	4	1	1 x 90 cm	Violet	LV 0386	
	5	2	1 x 90 cm	Violet	LV 0387	
	3.5	0	1 x 150 cm	Violet	LV 4035	Loop
	4	1	1 x 150 cm	Violet	LV 4040	Loop
 HRX 45	4	1	1 x 90 cm	Violet	LV 0390	
	5	2	1 x 90 cm	Violet	LV 0391	
 HRX 50	3.5	0	1 x 75 cm	Violet	LV 0294	
	4	1	1 x 75 cm	Violet	LV 0295	
	5	2	1 x 75 cm	Violet	LV 0296	
	4	1	1 x 90 cm	Violet	LV 0395	
	5	2	1 x 90 cm	Violet	LV 0396	
	5	2	1 x 110 cm	Violet	LV 0393	

CHIRASORB

UK03
74



Absorbable - Chirasorb

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRT: ½ circle, round body needle - tapercut						
HRT 17 	1.5	4/0	1 x 45 cm	Violet	LV 0566	
HRT 20 	1.5	4/0	1 x 75 cm	Violet	LV 0561	
	2	3/0	1 x 75 cm	Violet	LV 0562	
HRT 22 	2	3/0	1 x 75 cm	Violet	LV 0572	
HRT 25 	2	3/0	1 x 75 cm	Violet	LV 0581	
HRT 31 	3	2/0	1 x 75 cm	Violet	LV 0606	
HRT 37 	3	2/0	1 x 75 cm	Violet	LV 0592	
	3.5	0	1 x 75 cm	Violet	LV 0593	
	3	2/0	1 x 90 cm	Violet	LV 0597	
HRT 40 	3.5	0	1 x 90 cm	Violet	LV 0598	
	4	1	1 x 90 cm	Violet	LV 0599	
	3.5	0	1 x 75 cm	Violet	LV 0653	
HRN: ½ circle, round body needle - blunt point						
HRN 50 	4	1	1 x 90 cm	Violet	LV 0686	

CHIRASORB

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRF: ½ circle, flattened body needle, taper point						
HRF 17 	1 1.5	5/0 4/0	1 x 75 cm 1 x 75 cm	Violet Violet	LV 9242 LV 9256	
HRF 20 	1.5	4/0	1 x 75 cm	Violet	LV 9261	
HRF 26 	2 3	3/0 2/0	1 x 75 cm 1 x 75 cm	Violet Violet	LV 9266 LV 9267	
HRF 31 	3	2/0	1 x 75 cm	Violet	LV 9273	
DR: 3/8 circle, round body needle						
DR 15 	1.5	4/0	1 x 75 cm	Violet	LV 0724	
DR 21 	3	2/0	1 x 75 cm	Violet	LV 0744	
DR 30 	2	3/0	1 x 75 cm	Violet	LV 0746	
DRN: 3/8 circle, round body needle - blunt point						
DRN 60 	4	1	1 x 100 cm	Violet	LV 0699	

UK03
75



Absorbable - Chirasorb

CHIRASORB

UK03
76



Needle	EP	USP	Length	Colour	Cat. No	Other information
FR: 5/8 circle, round body needle						
FR 27	2	3/0	1 x 75 cm	Violet	LV 0609	
	3	2/0	1 x 75 cm	Violet	LV 0610	
	3.5	0	1 x 75 cm	Violet	LV 0611	
GR: Straight, round body needle						
GR 20	1.5	4/0	1 x 75 cm	Violet	LV 0402	
	3	2/0	1 x 90 cm	Violet	LV 0399	
KR: Ski, round body needle						
KR 25	3	2/0	1 x 20 cm	Violet	LV 2736	
	2	3/0	1 x 40 cm	Violet	LV 4630	Loop
KR 27	4	1	1 x 120 cm	Violet	LV 2390	
						
JRT: Hook shaped, round body needle - tapercut						
JRT 30	4	1	1 x 75 cm	Violet	LV 0694	
						
J-RX: Hook shaped, round body needle						
J-RX 30	3	2/0	1 x 75 cm	Violet	LV 1387	
						

CHIRASORB

Needle	EP	USP	Length	Colour	Cat. No	Other information
HS: ½ circle, reverse cutting needle						
 HS 18	2	3/0	1 x 75 cm	Violet	LV 0424	
	3	2/0	1 x 75 cm	Violet	LV 0425	
 HS 21	1.5	4/0	1 x 75 cm	Violet	LV 0300	
	2	3/0	1 x 75 cm	Violet	LV 0301	
 HS 25	2	3/0	1 x 75 cm	Violet	LV 0441	
	3	2/0	1 x 75 cm	Violet	LV 0442	
	3.5	0	1 x 75 cm	Violet	LV 0443	
 HS 26	3.5	0	1 x 75 cm	Violet	LV 1455	
 HS 30	3	2/0	1 x 75 cm	Violet	LV 0445	
	3.5	0	1 x 75 cm	Violet	LV 0446	
	4	1	1 x 75 cm	Violet	LV 0447	
 HS 36	3	2/0	1 x 75 cm	Violet	LV 0350	
	3.5	0	1 x 75 cm	Violet	LV 0351	
 HS 37	3.5	0	1 x 75 cm	Violet	LV 1451	
	4	1	1 x 75 cm	Violet	LV 1452	
 HS 40	5	2	1 x 90 cm	Violet	LV 0453	
	3	2/0	1 x 75 cm	Violet	LV 0455	
	3.5	0	1 x 75 cm	Violet	LV 0456	
 HS 40	4	1	1 x 75 cm	Violet	LV 0457	
	5	2	1 x 75 cm	Violet	LV 0458	
	3.5	0	1 x 90 cm	Violet	LV 0460	
	4	1	1 x 90 cm	Violet	LV 0461	
 HS 48	3.5	0	1 x 75 cm	Violet	LV 0466	
	4	1	1 x 90 cm	Violet	LV 0469	
	5	2	1 x 90 cm	Violet	LV 0470	



CHIRASORB

UK03
78



Absorbable - Chirasorb

Needle	EP	USP	Length	Colour	Cat. No	Other information
HSX: 1/2 circle, reverse cutting needle - extra strong						
HSX 26 	4	1	1 x 75 cm	Violet	LV 1456	
	5	2	1 x 75 cm	Violet	LV 1454	
HSX 48 	3.5	0	1 x 90 cm	Violet	LV 0464	
DS: 3/8 circle, reverse cutting needle						
DS 11 	0.5	7/0	1 x 45 cm	Violet	LV 0193	
DS 12 	0.7	6/0	1 x 45 cm	Violet	LV 0201	
	1	5/0	1 x 45 cm	Violet	LV 0202	
	1.5	4/0	1 x 45 cm	Violet	LV 0203	
DS 15 	0.7	6/0	1 x 45 cm	Violet	LV 0200	
	1	5/0	1 x 45 cm	Violet	LV 0205	
	1.5	4/0	1 x 45 cm	Violet	LV 0206	
	2	3/0	1 x 45 cm	Violet	LV 0207	
DS 19 	1.5	4/0	1 x 75 cm	Violet	LV 0209	
	1	5/0	1 x 45 cm	Violet	LV 0214	
	2	3/0	1 x 45 cm	Violet	LV 0216	
DS 24 	1.5	4/0	1 x 75 cm	Violet	LV 0211	
	2	3/0	1 x 75 cm	Violet	LV 0212	
	3	2/0	1 x 75 cm	Violet	LV 0213	
DS 25 	2	3/0	1 x 90 cm	Violet	LV 2223	
DS 26 	3	2/0	1 x 90 cm	Violet	LV 2224	
	2	3/0	1 x 45 cm	Violet	LV 1221	

CHIRASORB

UK03
79

Absorbable - Chirasorb

Needle	EP	USP	Length	Colour	Cat. No	Other information
--------	----	-----	--------	--------	---------	-------------------

DS: 3/8 circle, reverse cutting needle

 <p>DS 30</p>	3	2/0	1 x 75 cm	Violet	LV 0234	
	3.5	0	1 x 75 cm	Violet	LV 0235	
	3	2/0	1 x 90 cm	Violet	LV 1239	

DSC: 3/8 circle, conventional cutting needle

 <p>DSC 16</p>	1.5	4/0	1 x 75 cm	Violet	LV 8208	
---	-----	-----	-----------	--------	---------	--

GS: Straight, reverse cutting needle

 <p>GS 13</p>	1.5	4/0	1 x 45 cm	Violet	LV 1198	
--	-----	-----	-----------	--------	---------	--



CHIRASORB natural ligatures

UK03
80



Absorbable - Chirasorb

Pre-cut threads

EP	USP	Length	Cat. No
1	5/0	3 x 45 cm	LN 0012
1.5	4/0	3 x 45 cm	LN 0013
2	3/0	3 x 45 cm	LN 0014
3	2/0	3 x 45 cm	LN 0015
3.5	0	3 x 45 cm	LN 0016
2	3/0	10 x 45 cm	LN 0064
3	2/0	10 x 45 cm	LN 0065
1.5	4/0	2 x 75 cm	LN 0023
2	3/0	2 x 75 cm	LN 0024
3	2/0	2 x 75 cm	LN 0025

Chirapack - 20 pcs./box

EP	USP	Length	Cat. No
1.5	4/0	1 x 300 cm	LN 0003
2	3/0	1 x 300 cm	LN 0004
3	2/0	1 x 300 cm	LN 0005
3.5	0	1 x 300 cm	LN 0006

CHIRASORB

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
HR 13 	1	5/0	1 x 45 cm	Natural	LN 0252	
	1.5	4/0	1 x 45 cm	Natural	LN 0253	
	1.5	4/0	1 x 75 cm	Natural	LN 0153	
HR 18 	1.5	4/0	1 x 75 cm	Natural	LN 0256	
	2	3/0	1 x 75 cm	Natural	LN 0257	
HR 22 	2	3/0	1 x 75 cm	Natural	LN 0262	
	3	2/0	1 x 75 cm	Natural	LN 0263	
HR 27 	2	3/0	1 x 75 cm	Natural	LN 0266	
HR 30 	3.5	0	1 x 75 cm	Natural	LN 0274	

CHIRASORB

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRT: ½ circle, round body needle - tapercut						
HRT 17 	1.5 2	4/0 3/0	1 x 75 cm 1 x 75 cm	Natural Natural	LN 0561 LN 0562	
HRT 22 	2	3/0	1 x 75 cm	Natural	LN 0572	
DS: 3/8 circle, reverse cutting needle						
DS 12 	0.7 1	6/0 5/0	1 x 45 cm 1 x 45 cm	Natural Natural	LN 1201 LN 1202	
DS 15 	1.5 2	4/0 3/0	1 x 45 cm 1 x 45 cm	Natural Natural	LN 1206 LN 1207	
DS 16 	0.7	6/0	1 x 45 cm	Natural	LN 1208	
DS 19 	1 1.5 2	5/0 4/0 3/0	1 x 45 cm 1 x 45 cm 1 x 45 cm	Natural Natural Natural	LN 1215 LN 1216 LN 1217	
	1.5 2 3	4/0 3/0 2/0	1 x 75 cm 1 x 75 cm 1 x 75 cm	Natural Natural Natural	LN 1211 LN 1212 LN 1213	
DS 25 	2	3/0	1 x 45 cm	Natural	LN 1332	
	1.5 2 3	4/0 3/0 2/0	1 x 75 cm 1 x 75 cm 1 x 75 cm	Natural Natural Natural	LN 1231 LN 1232 LN 1233	
DS 30 	3 3.5	2/0 0	1 x 75 cm 1 x 75 cm	Natural Natural	LN 1243 LN 1244	



CHIRASORB

UK03
82



Absorbable - Chirasorb

Needle	EP	USP	Length	Colour	Cat. No	Other information
DSM: 3/8 circle, reverse cutting needle - micropoint						
DSM 12 	0.7	6/0	1 x 45 cm	Natural	LN 0201	Advanced cutting needle
	1	5/0	1 x 45 cm	Natural	LN 0202	
	1.5	4/0	1 x 45 cm	Natural	LN 0203	
DSM 15 	1	5/0	1 x 45 cm	Natural	LN 0206	Advanced cutting needle
	1.5	4/0	1 x 45 cm	Natural	LN 0207	
	2	3/0	1 x 45 cm	Natural	LN 0208	
DSM 19 	1.5	4/0	1 x 45 cm	Natural	LN 0216	Advanced cutting needle
	1.5	4/0	1 x 75 cm	Natural	LN 0212	
	2	3/0	1 x 75 cm	Natural	LN 0213	
DSM 25 	2	3/0	1 x 75 cm	Natural	LN 0223	Advanced cutting needle
	3	2/0	1 x 75 cm	Natural	LN 0224	

Chirasorb


Chirasorb Plus



Chirasorb

Chirasorb

Chirasorb



CHIRASORB Plus

Antibacterial

CHIRASORB Plus is a mid-term absorbable, antibacterial coated, synthetic, multifilament suture made of 90% glycolide and 10% L-lactide. The braided suture is covered with antibacterial coating, which consists of copolymer glycolide and lactide, calcium stearate and chlorhexidine diacetate.

Chlorhexidine is a cationic antiseptic with a wide range of antimicrobial effects. It is effective on both Gram-positive and Gram-negative bacteria; the mechanism of action is cytoplasmic membrane disruption.

Sterile Chirasorb Plus is available in violet or natural. The product elicits a minimal acute inflammatory reaction in tissues.

Progressive loss of tensile strength and absorption of Chirasorb Plus occurs by hydrolysis, when the polymer is degraded to glycolic and lactic acids, which are subsequently absorbed and metabolised by the body.



UK03

85



Absorbable - Chirasorb Plus

Indication: General surgery, Gynaecology, Skin Closure (e.g. intracutaneous, subcutaneous), Muscular layer and Ligation.

Contraindication: Chirasorb Plus should not be used in tissues subject to expansion, stretching or distension, or which might require long-term mechanical support. Chirasorb Plus should not be used in patients with known allergic reactions to Chlorhexidine.

Sizes: 6/0 to 4 USP

Structure: Multifilament

Absorption method: Hydrolysis

Strength Retention Profile: MID-TERM HEALING (35 days)

1 Day	100%
14 Days	75%
21 Days	50%

Mass Absorption profile: 56 - 70 days

CHIRASORB Plus violet ligatures



UK03
86



Absorbable - Chirasorb Plus

Pre-cut threads

EP	USP	Length	Cat. No
1	5/0	3 x 45 cm	CV 0012
1.5	4/0	3 x 45 cm	CV 0013
2	3/0	3 x 45 cm	CV 0014
3	2/0	3 x 45 cm	CV 0015
3.5	0	3 x 45 cm	CV 0016
4	1	3 x 45 cm	CV 0017
5	2	3 x 45 cm	CV 0018
1	5/0	6 x 45 cm	CV 0042
1.5	4/0	6 x 45 cm	CV 0043
2	3/0	6 x 45 cm	CV 0044
3	2/0	6 x 45 cm	CV 0045
3.5	0	6 x 45 cm	CV 0046
4	1	6 x 45 cm	CV 0047
5	2	6 x 45 cm	CV 0048
2	3/0	10 x 45 cm	CV 0064
3	2/0	10 x 45 cm	CV 0065
3.5	0	10 x 45 cm	CV 0066
4	1	10 x 45 cm	CV 0067
5	2	10 x 45 cm	CV 0068
1	5/0	2 x 75 cm	CV 0022
1.5	4/0	2 x 75 cm	CV 0023
2	3/0	2 x 75 cm	CV 0024
3	2/0	2 x 75 cm	CV 0025
3.5	0	2 x 75 cm	CV 0026
4	1	2 x 75 cm	CV 0027
5	2	2 x 75 cm	CV 0028
1.5	4/0	5 x 75 cm	CV 0053
2	3/0	5 x 75 cm	CV 0054
3	2/0	5 x 75 cm	CV 0055

Pre-cut threads

EP	USP	Length	Cat. No
3.5	0	5 x 75 cm	CV 0056
4	1	5 x 75 cm	CV 0057
5	2	5 x 75 cm	CV 0058
1.5	4/0	1 x 150 cm	CV 0083
2	3/0	1 x 150 cm	CV 0084
3	2/0	1 x 150 cm	CV 0085
3.5	0	1 x 150 cm	CV 0089
4	1	1 x 150 cm	CV 0087
5	2	1 x 150 cm	CV 0090

Chirapack - 20 pcs./box

EP	USP	Length	Cat. No
2	3/0	1 x 250 cm	CV 1004
3	2/0	1 x 250 cm	CV 1005
3.5	0	1 x 250 cm	CV 1006
4	1	1 x 250 cm	CV 1007
5	2	1 x 250 cm	CV 1008
1.5	4/0	1 x 300 cm	CV 0003
2	3/0	1 x 300 cm	CV 0004
3	2/0	1 x 300 cm	CV 0005
3.5	0	1 x 300 cm	CV 0006
4	1	1 x 300 cm	CV 0007
5	2	1 x 300 cm	CV 0008
6	3 + 4	1 x 300 cm	CV 0009

CHIRASORB Plus

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
 HR 10	1	5/0	1 x 75 cm	Violet	CV 0248	
	1.5	4/0	1 x 45 cm	Violet	CV 0252	
 HR 13	2	3/0	1 x 45 cm	Violet	CV 0250	
	0.7	6/0	1 x 75 cm	Violet	CV 1251	
	1	5/0	1 x 75 cm	Violet	CV 0254	
 HR 18	1.5	4/0	1 x 45 cm	Violet	CV 0258	
	2	3/0	1 x 45 cm	Violet	CV 0259	
	1	5/0	1 x 75 cm	Violet	CV 0242	
	1.5	4/0	1 x 75 cm	Violet	CV 0256	
 HR 22	2	3/0	1 x 75 cm	Violet	CV 0257	
	3	2/0	1 x 75 cm	Violet	CV 0260	
	1.5	4/0	6 x 45 cm	Violet	CV 0526	
	2	3/0	6 x 45 cm	Violet	CV 0527	
	1.5	4/0	1 x 75 cm	Violet	CV 0261	
 HR 27	2	3/0	1 x 75 cm	Violet	CV 0262	
	3	2/0	1 x 75 cm	Violet	CV 0263	
	2	3/0	4 x 75 cm	Violet	CV 0322	
	3	2/0	4 x 75 cm	Violet	CV 0323	
	1.5	4/0	1 x 45 cm	Violet	CV 0245	
	2	3/0	6 x 45 cm	Violet	CV 0537	
 2 x HR 27	3	2/0	6 x 45 cm	Violet	CV 0538	
	1.5	4/0	1 x 75 cm	Violet	CV 0265	
	2	3/0	1 x 75 cm	Violet	CV 0266	
	3	2/0	1 x 75 cm	Violet	CV 0267	
 2 x HR 27	3.5	0	1 x 75 cm	Violet	CV 0268	
	2	3/0	1 x 90 cm	Violet	CV 0366	
 2 x HR 27	3	2/0	1 x 90 cm	Violet	CV 0367	
	2	3/0	1 x 90 cm	Violet	CV 2466	
 2 x HR 27	3	2/0	1 x 90 cm	Violet	CV 2467	




CHIRASORB Plus

UK03
88

Absorbable - Chirasorb Plus

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
 <p>HR 30</p>	2	3/0	1 x 45 cm	Violet	CV 0169	
	2	3/0	1 x 75 cm	Violet	CV 0269	
	3	2/0	1 x 75 cm	Violet	CV 0273	
	3.5	0	1 x 75 cm	Violet	CV 0274	
	4	1	1 x 75 cm	Violet	CV 1275	
	3.5	0	1 x 90 cm	Violet	CV 0371	
	<hr/>					
 <p>HR 35</p>	2	3/0	1 x 75 cm	Violet	CV 0270	
	3	2/0	1 x 75 cm	Violet	CV 0271	
	3.5	0	1 x 75 cm	Violet	CV 0272	
	4	1	1 x 75 cm	Violet	CV 1273	
	5	2	1 x 75 cm	Violet	CV 1272	
	3	2/0	1 x 90 cm	Violet	CV 0374	
	3.5	0	1 x 90 cm	Violet	CV 0375	
	4	1	1 x 90 cm	Violet	CV 0376	
	<hr/>					
	 <p>HR 40</p>	3	2/0	1 x 75 cm	Violet	CV 0275
3.5		0	1 x 75 cm	Violet	CV 0276	
4		1	1 x 75 cm	Violet	CV 0277	
 <p>HR 48</p>	3	2/0	1 x 90 cm	Violet	CV 1475	
	3.5	0	1 x 90 cm	Violet	CV 1476	
	4	1	1 x 90 cm	Violet	CV 1477	
	5	2	1 x 90 cm	Violet	CV 1478	
	3.5	0	1 x 75 cm	Violet	CV 1402	
	4	1	1 x 75 cm	Violet	CV 1396	
	5	2	1 x 75 cm	Violet	CV 1400	
	3.5	0	1 x 90 cm	Violet	CV 1403	
	4	1	1 x 90 cm	Violet	CV 1401	
	5	2	1 x 90 cm	Violet	CV 1399	




CHIRASORB Plus

UK03
89



Absorbable - Chirisorb Plus

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
HR 50 	3.5	0	1 x 90 cm	Violet	CV 1392	
HR 60 	3 3.5	2/0 0	1 x 75 cm 1 x 75 cm	Violet Violet	CV 0278 CV 0279	
HR 65 	4 5	1 2	1 x 75 cm 1 x 75 cm	Violet Violet	CV 1285 CV 1286	
HR 75 	5	2	1 x 75 cm	Violet	CV 0496	
HR s: ½ circle, round body needle - slim						
HR s 18 	1.5 2	4/0 3/0	1 x 75 cm 1 x 75 cm	Violet Violet	CV 1256 CV 1257	Slim needle
HR s 27 Black 	3	2/0	1 x 75 cm	Violet	CV 2267	Slim, black needle

CHIRASORB Plus

UK03
90



Absorbable - Chirasorb Plus

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRX: ½ circle, round body needle - extra strong						
HRX 22 	2	3/0	1 x 75 cm	Violet	CV 0420	
HRX 27 	2	3/0	1 x 75 cm	Violet	CV 0429	
	3	2/0	1 x 75 cm	Violet	CV 0430	
	3.5	0	1 x 75 cm	Violet	CV 0431	
	4	1	1 x 75 cm	Violet	CV 0432	
HRX 30 	2	3/0	1 x 75 cm	Violet	CV 0426	
	4	1	1 x 90 cm	Violet	CV 1438	
HRX 35 	3	2/0	1 x 75 cm	Violet	CV 0280	
	3.5	0	1 x 75 cm	Violet	CV 0281	
	4	1	1 x 75 cm	Violet	CV 0282	
	3.5	0	1 x 90 cm	Violet	CV 0381	
HRX 40 	3.5	0	1 x 75 cm	Violet	CV 0285	
	4	1	1 x 75 cm	Violet	CV 0286	
	3.5	0	1 x 90 cm	Violet	CV 0385	
	4	1	1 x 90 cm	Violet	CV 0386	
	5	2	1 x 90 cm	Violet	CV 0387	
	3.5	0	1 x 150 cm	Violet	CV 4035	Loop
HRX 45 	4	1	1 x 75 cm	Violet	CV 0290	
	5	2	1 x 75 cm	Violet	CV 0291	
	4	1	1 x 90 cm	Violet	CV 0390	
	5	2	1 x 90 cm	Violet	CV 0391	
HRX 50 	3.5	0	1 x 75 cm	Violet	CV 0294	
	4	1	1 x 75 cm	Violet	CV 0295	
	5	2	1 x 75 cm	Violet	CV 0296	
	4	1	1 x 90 cm	Violet	CV 0395	
HRX 50 	5	2	1 x 90 cm	Violet	CV 0396	
	5	2	1 x 110 cm	Violet	CV 0393	
	4	1	1 x 150 cm	Violet	CV 4540	Loop

CHIRASORB Plus

UK03
91

Absorbable - Chirasorb Plus

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRT: ½ circle, round body needle - tapercut						
HRT 17 	1.5	4/0	1 x 45 cm	Violet	CV 0566	
	1.5 2	4/0 3/0	1 x 75 cm 1 x 75 cm	Violet Violet	CV 0561 CV 0562	
HRT 25 	2	3/0	1 x 75 cm	Violet	CV 0581	
	3	2/0	1 x 75 cm	Violet	CV 0582	
	3.5	0	1 x 75 cm	Violet	CV 0583	
HRT 30 	3.5	0	1 x 75 cm	Violet	CV 0586	
	4	1	1 x 75 cm	Violet	CV 0587	
HRT 31 	4	1	1 x 90 cm	Violet	CV 0589	
	3	2/0	1 x 75 cm	Violet	CV 0606	
HRT 37 	3	2/0	1 x 75 cm	Violet	CV 0592	
	3.5	0	1 x 75 cm	Violet	CV 0593	
	3 3.5 4	2/0 0 1	1 x 90 cm 1 x 90 cm 1 x 90 cm	Violet Violet Violet	CV 0597 CV 0598 CV 0599	
HRT 40 	3.5	0	1 x 75 cm	Violet	CV 0653	
	4	1	1 x 90 cm	Violet	CV 0651	
HRT 45 	3	2/0	1 x 75 cm	Violet	CV 0663	
	5	2	1 x 75 cm	Violet	CV 0665	
HRT 48 	3	2/0	1 x 75 cm	Violet	CV 1671	
	5	2	1 x 90 cm	Violet	CV 1668	


CHIRASORB Plus

UK03
92

Needle	EP	USP	Length	Colour	Cat. No	Other information
--------	----	-----	--------	--------	---------	-------------------

HRN: ½ circle, round body needle - blunt point

HRN 50



4	1	1 x 90 cm	Violet	CV 0686
---	---	-----------	--------	---------

HRF: ½ circle, flattened body needle, taper point

HRF 17



1	5/0	1 x 75 cm	Violet	CV 9242
1.5	4/0	1 x 75 cm	Violet	CV 9256

HRF 22



1.5	4/0	1 x 75 cm	Violet	CV 9344
2	3/0	1 x 75 cm	Violet	CV 9343

HRF 26



2	3/0	1 x 75 cm	Violet	CV 9266
3	2/0	1 x 75 cm	Violet	CV 9267

2	3/0	1 x 90 cm	Violet	CV 9366
---	-----	-----------	--------	---------

HRF 31



2	3/0	1 x 75 cm	Violet	CV 9269
3	2/0	1 x 75 cm	Violet	CV 9273

HRF 36



2	3/0	1 x 75 cm	Violet	CV 9577
3	2/0	1 x 75 cm	Violet	CV 9236
3.5	0	1 x 75 cm	Violet	CV 9235

DR: 3/8 circle, round body needle

DR 15



1.5	4/0	1 x 75 cm	Violet	CV 0724
-----	-----	-----------	--------	---------

DR 21



3	2/0	1 x 75 cm	Violet	CV 0744
---	-----	-----------	--------	---------

DR 30



2	3/0	1 x 75 cm	Violet	CV 0746
---	-----	-----------	--------	---------


CHIRASORB Plus

UK03
93

Absorbable - Chirasorb Plus

Needle	EP	USP	Length	Colour	Cat. No	Other information
--------	----	-----	--------	--------	---------	-------------------

DRN: 3/8 circle, round body needle - blunt point

<p>DRN 60</p> 	4	1	1 x 100 cm	Violet	CV 0699	
--	---	---	------------	--------	---------	--

FR: 5/8 circle, round body needle

<p>FR 27</p> 	3 3.5	2/0 0	1 x 75 cm 1 x 75 cm	Violet Violet	CV 0610 CV 0611	
--	----------	----------	------------------------	------------------	--------------------	--

<p>FR 38</p> 	3 3.5	2/0 0	1 x 75 cm 1 x 75 cm	Violet Violet	CV 0620 CV 0621	
---	----------	----------	------------------------	------------------	--------------------	--

GR: Straight, round body needle

	1.5 2	4/0 3/0	1 x 45 cm 1 x 45 cm	Violet Violet	CV 0406 CV 0407	
<p>GR 20</p> 	1.5	4/0	1 x 75 cm	Violet	CV 0402	
	3	2/0	1 x 90 cm	Violet	CV 0399	

ART: Asymmetric, round body needle - tapercut

<p>ART 30 Black</p> 	1.5	4/0	1 x 45 cm	Violet	CV 0630	Black needle
---	-----	-----	-----------	--------	---------	--------------

JRT: Hook shaped, round body needle - tapercut

<p>JRT 30</p> 	4	1	1 x 75 cm	Violet	CV 0694	
---	---	---	-----------	--------	---------	--

CHIRASORB Plus

UK03
94



Absorbable - Chirascorb Plus

Needle	EP	USP	Length	Colour	Cat. No	Other information
--------	----	-----	--------	--------	---------	-------------------

KR: Ski, round body needle

KR 27 	2	3/0	1 x 40 cm	Violet	CV 4630	Loop
--	---	-----	-----------	--------	---------	------

HS: 1/2 circle, reverse cutting needle

HS 18 	2 3	3/0 2/0	1 x 75 cm 1 x 75 cm	Violet Violet	CV 0424 CV 0425	
HS 21 	1.5	4/0	1 x 75 cm	Violet	CV 0300	
HS 25 	2 3 3.5	3/0 2/0 0	1 x 75 cm 1 x 75 cm 1 x 75 cm	Violet Violet Violet	CV 0441 CV 0442 CV 0443	
HS 30 	3 3.5 4	2/0 0 1	1 x 75 cm 1 x 75 cm 1 x 75 cm	Violet Violet Violet	CV 0445 CV 0446 CV 0447	
HS 37 	3.5 4	0 1	1 x 75 cm 1 x 75 cm	Violet Violet	CV 1451 CV 1452	
HS 37 	3.5 5	0 2	1 x 90 cm 1 x 90 cm	Violet Violet	CV 0451 CV 0453	
HS 40 	3 3.5 4 5	2/0 0 1 2	1 x 75 cm 1 x 75 cm 1 x 75 cm 1 x 75 cm	Violet Violet Violet Violet	CV 0455 CV 0456 CV 0457 CV 0458	
HS 40 	3.5 4	0 1	1 x 90 cm 1 x 90 cm	Violet Violet	CV 0460 CV 0461	
HS 48 	3.5	0	1 x 75 cm	Violet	CV 0466	
HS 48 	4 5	1 2	1 x 90 cm 1 x 90 cm	Violet Violet	CV 0469 CV 0470	

CHIRASORB Plus

Needle	EP	USP	Length	Colour	Cat. No	Other information
DS: 3/8 circle, reverse cutting needle						
DS 12 	0.7	6/0	1 x 45 cm	Violet	CV 0201	
	1	5/0	1 x 45 cm	Violet	CV 0202	
	1.5	4/0	1 x 45 cm	Violet	CV 0203	
DS 15 	1	5/0	1 x 45 cm	Violet	CV 0205	
	1.5	4/0	1 x 45 cm	Violet	CV 0206	
	2	3/0	1 x 45 cm	Violet	CV 0207	
DS 19 	1	5/0	1 x 45 cm	Violet	CV 0214	
	1.5	4/0	1 x 45 cm	Violet	CV 0215	
	2	3/0	1 x 45 cm	Violet	CV 0216	
	1.5	4/0	1 x 75 cm	Violet	CV 0211	
	2	3/0	1 x 75 cm	Violet	CV 0212	
	3	2/0	1 x 75 cm	Violet	CV 0213	
DS 25 	1.5	4/0	1 x 75 cm	Violet	CV 0222	
	2	3/0	1 x 75 cm	Violet	CV 0223	
	3	2/0	1 x 75 cm	Violet	CV 0224	
DS 30 	3	2/0	1 x 75 cm	Violet	CV 0234	
	3.5	0	1 x 75 cm	Violet	CV 0235	
DSM: 3/8 circle, reverse cutting needle - micropoint						
DSM 12 	1	5/0	1 x 45 cm	Violet	CV 0302	Advanced cutting needle
GS: Straight, reverse cutting needle						
GS 13 	1.5	4/0	1 x 45 cm	Violet	CV 1198	



CHIRASORB Plus natural ligatures



UK03
96

Pre-cut threads

EP	USP	Length	Cat. No
1	5/0	3 x 45 cm	CN 0012
1.5	4/0	3 x 45 cm	CN 0013
2	3/0	3 x 45 cm	CN 0014
3	2/0	3 x 45 cm	CN 0015
3.5	0	3 x 45 cm	CN 0016
2	3/0	10 x 45 cm	CN 0064
3	2/0	10 x 45 cm	CN 0065
1	5/0	2 x 75 cm	CN 0022
1.5	4/0	2 x 75 cm	CN 0023
2	3/0	2 x 75 cm	CN 0024
3	2/0	2 x 75 cm	CN 0025
1.5	4/0	5 x 75 cm	CN 0053
2	3/0	5 x 75 cm	CN 0054
3	2/0	5 x 75 cm	CN 0055
3.5	0	5 x 75 cm	CN 0056
4	1	5 x 75 cm	CN 0057
5	2	5 x 75 cm	CN 0058

Chirapack - 20 pcs./box

EP	USP	Length	Cat. No
1.5	4/0	1 x 300 cm	CN 0003
2	3/0	1 x 300 cm	CN 0004
3	2/0	1 x 300 cm	CN 0005
3.5	0	1 x 300 cm	CN 0006



Absorbable - Chirasorb Plus

CHIRASORB Plus

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
HR 13 	1 1.5	5/0 4/0	1 x 45 cm 1 x 45 cm	Natural Natural	CN 0252 CN 0253	
HR 18 	1.5 2	4/0 3/0	1 x 75 cm 1 x 75 cm	Natural Natural	CN 0256 CN 0257	
HR 22 	2 3	3/0 2/0	1 x 75 cm 1 x 75 cm	Natural Natural	CN 0262 CN 0263	
HR 27 	2	3/0	1 x 75 cm	Natural	CN 0266	
HR 30 	3.5	0	1 x 75 cm	Natural	CN 0274	
HRT: ½ circle, round body needle - tapercut						
HRT 17 	1.5 2	4/0 3/0	1 x 75 cm 1 x 75 cm	Natural Natural	CN 0561 CN 0562	
HRT 22 	2	3/0	1 x 75 cm	Natural	CN 0572	
HRT 40 	3 3.5	2/0 0	1 x 75 cm 1 x 75 cm	Natural Natural	CN 0694 CN 0653	



CHIRASORB Plus

UK03
98



Absorbable - Chirisorb Plus

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRF: ½ circle, flattened body needle, taper point						
HRF 13 	1	5/0	1 x 45 cm	Natural	CN 9252	
HRF 22 	2	3/0	1 x 75 cm	Natural	CN 9262	
HRF 26 	3	2/0	1 x 75 cm	Natural	CN 9267	
DS: 3/8 circle, reverse cutting needle						
DS 12 	0.7 1	6/0 5/0	1 x 45 cm 1 x 45 cm	Natural Natural	CN 1201 CN 1202	
DS 15 	1.5 2	4/0 3/0	1 x 45 cm 1 x 45 cm	Natural Natural	CN 1206 CN 1207	
DS 19 	1 1.5 2	5/0 4/0 3/0	1 x 45 cm 1 x 45 cm 1 x 45 cm	Natural Natural Natural	CN 1215 CN 1216 CN 1217	
	1.5 2 3	4/0 3/0 2/0	1 x 75 cm 1 x 75 cm 1 x 75 cm	Natural Natural Natural	CN 1211 CN 1212 CN 1213	

Chirlac

chirlac

Chirlac

Chirlac



CHIRLAC

CHIRLAC is a mid-term absorbable, synthetic, multifilament suture material of extremely high tensile strength made of the homopolymer of hydroacetic acid (glycolic acid).

The thread has a homogeneous structure and surface, which gives it low frictional properties and holds firmly without slipping when knotted.

Chirlac is coated with an absorbable substance, which reduces the capillary effect, provides excellent holding properties and diminishes the tendency to cut into the flesh, while ensuring that the knots hold firmly.

Chirlac is degraded in the body by a simple hydrolytic process. The degraded products are metabolised completely in the body. Absorption begins as a loss of tensile strength without appreciable loss of mass.



Indication: General surgery, Gynaecology, Ophthalmic, Skin Closure (e.g. intracutaneous, subcutaneous), Muscular layer and Ligation.

Contraindication: Chirlac should not be used in tissues subject to expansion, stretching or distension, or which might require long-term mechanical support.

Sizes: 9/0 to 6 USP

Structure: Multifilament

Absorption method: Hydrolysis

Strength Retention Profile:

MID-TERM HEALING (28 - 35 days)	
1 Day	100%
14 Days	65% - 75%
21 Days	45% - 50%

Mass Absorption profile: 80 - 90 days

CHIRLAC violet ligatures

UK03
102



Absorbable - Chirlac

Pre-cut threads


EP	USP	Length	Cat. No
1	5/0	3 x 45 cm	PG 0012
1.5	4/0	3 x 45 cm	PG 0013
2	3/0	3 x 45 cm	PG 0014
3	2/0	3 x 45 cm	PG 0015
3.5	0	3 x 45 cm	PG 0016
4	1	3 x 45 cm	PG 0017
5	2	3 x 45 cm	PG 0018
1	5/0	6 x 45 cm	PG 0042
1.5	4/0	6 x 45 cm	PG 0043
2	3/0	6 x 45 cm	PG 0044
3	2/0	6 x 45 cm	PG 0045
3.5	0	6 x 45 cm	PG 0046
4	1	6 x 45 cm	PG 0047
5	2	6 x 45 cm	PG 0048
2	3/0	10 x 45 cm	PG 0064
3	2/0	10 x 45 cm	PG 0065
3.5	0	10 x 45 cm	PG 0066
4	1	10 x 45 cm	PG 0067
5	2	10 x 45 cm	PG 0068
2	3/0	12 x 45 cm	PG 4564
1	5/0	2 x 75 cm	PG 0022
1.5	4/0	2 x 75 cm	PG 0023
2	3/0	2 x 75 cm	PG 0024
3	2/0	2 x 75 cm	PG 0025
3.5	0	2 x 75 cm	PG 0026
4	1	2 x 75 cm	PG 0027
5	2	2 x 75 cm	PG 0028
1.5	4/0	5 x 75 cm	PG 0053
2	3/0	5 x 75 cm	PG 0054
3	2/0	5 x 75 cm	PG 0055
3.5	0	5 x 75 cm	PG 0056
4	1	5 x 75 cm	PG 0057
5	2	5 x 75 cm	PG 0058
2	3/0	1 x 140 cm	PG 4084

Pre-cut threads

EP	USP	Length	Cat. No
1.5	4/0	1 x 150 cm	PG 0083
2	3/0	1 x 150 cm	PG 0084
3	2/0	1 x 150 cm	PG 0085
3.5	0	1 x 150 cm	PG 0089
4	1	1 x 150 cm	PG 0087
5	2	1 x 150 cm	PG 0090

Chirapack - 20 pcs./box

EP	USP	Length	Cat. No
2	3/0	1 x 250 cm	PG 1004
3	2/0	1 x 250 cm	PG 1005
3.5	0	1 x 250 cm	PG 1006
4	1	1 x 250 cm	PG 1007
5	2	1 x 250 cm	PG 1008
1.5	4/0	1 x 300 cm	PG 0003
2	3/0	1 x 300 cm	PG 0004
3	2/0	1 x 300 cm	PG 0005
3.5	0	1 x 300 cm	PG 0006
4	1	1 x 300 cm	PG 0007
5	2	1 x 300 cm	PG 0008
6	3 + 4	1 x 300 cm	PG 0009



CHIRLAC

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
HR 10 	0.7	6/0	1 x 45 cm	Violet	PG 0246	
	1	5/0	1 x 75 cm	Violet	PG 0248	
HR 13 	0.7	6/0	1 x 45 cm	Violet	PG 0251	
	1	5/0	1 x 45 cm	Violet	PG 0252	
	1.5	4/0	1 x 45 cm	Violet	PG 0253	
	2	3/0	1 x 45 cm	Violet	PG 0250	
	0.7	6/0	1 x 75 cm	Violet	PG 1251	
HR 17 	1	5/0	1 x 75 cm	Violet	PG 0249	
	1.5	4/0	1 x 45 cm	Violet	PG 0258	
	2	3/0	1 x 45 cm	Violet	PG 0259	
HR 18 	0.7	6/0	1 x 75 cm	Violet	PG 0241	
	1	5/0	1 x 75 cm	Violet	PG 0242	
	1.5	4/0	1 x 75 cm	Violet	PG 0256	
	2	3/0	1 x 75 cm	Violet	PG 0257	
	3	2/0	1 x 75 cm	Violet	PG 0260	
	2 x HR 18 	1.5	4/0	1 x 90 cm	Violet	PG 2265
HR 20 	1.5	4/0	1 x 75 cm	Violet	PG 1261	
HR 22 	1.5	4/0	6 x 45 cm	Violet	PG 0526	
	2	3/0	6 x 45 cm	Violet	PG 0527	
	1.5	4/0	1 x 75 cm	Violet	PG 0261	
	2	3/0	1 x 75 cm	Violet	PG 0262	
	3	2/0	1 x 75 cm	Violet	PG 0263	
	2	3/0	4 x 75 cm	Violet	PG 0322	
HR 22 A 	3	2/0	4 x 75 cm	Violet	PG 0323	
	2	3/0	6 x 45 cm	Violet	PG 0327	Take-off needle
2 x HR 22 	3	2/0	1 x 75 cm	Violet	PG 0532	
	2	3/0	1 x 90 cm	Violet	PG 0533	
	3	2/0	1 x 90 cm	Violet	PG 0534	
HR 26 	2	3/0	1 x 75 cm	Violet	PG 0535	
	3	2/0	1 x 75 cm	Violet	PG 0536	

CHIRLAC


UK03
104



Absorbable - Chirlac

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
 <p>HR 27</p>	1.5	4/0	1 x 45 cm	Violet	PG 0245	
	3	2/0	6 x 45 cm	Violet	PG 0538	
	1.5	4/0	1 x 75 cm	Violet	PG 0265	
	2	3/0	1 x 75 cm	Violet	PG 0266	
	3	2/0	1 x 75 cm	Violet	PG 0267	
 <p>HR 30</p>	3.5	0	1 x 75 cm	Violet	PG 0268	
	2	3/0	1 x 90 cm	Violet	PG 0366	
	3	2/0	1 x 90 cm	Violet	PG 0367	
	2	3/0	1 x 45 cm	Violet	PG 0169	
 <p>2 x HR 30</p>	2	3/0	1 x 75 cm	Violet	PG 0269	
	3	2/0	1 x 75 cm	Violet	PG 0273	
	3.5	0	1 x 75 cm	Violet	PG 0274	
	4	1	1 x 75 cm	Violet	PG 1275	
	3	2/0	1 x 90 cm	Violet	PG 0364	
 <p>HR 30 A</p>	3.5	0	4 x 75 cm	Violet	PG 3274	Take-off needle
	3.5	0	1 x 90 cm	Violet	PG 0371	
 <p>HR 35</p>	2	3/0	1 x 75 cm	Violet	PG 0270	
	3	2/0	1 x 75 cm	Violet	PG 0271	
	3.5	0	1 x 75 cm	Violet	PG 0272	
	4	1	1 x 75 cm	Violet	PG 1273	
	5	2	1 x 75 cm	Violet	PG 1272	
	2	3/0	1 x 90 cm	Violet	PG 0373	
 <p>HR 35 A</p>	3	2/0	4 x 75 cm	Violet	PG 3374	Take-off needle
	3	2/0	1 x 90 cm	Violet	PG 0374	
	3.5	0	1 x 90 cm	Violet	PG 0375	
	4	1	1 x 90 cm	Violet	PG 0376	

CHIRLAC

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
 <p>HR 40</p>	3	2/0	1 x 75 cm	Violet	PG 0275	
	3.5	0	1 x 75 cm	Violet	PG 0276	
	4	1	1 x 75 cm	Violet	PG 0277	
 <p>HR 48</p>	3	2/0	1 x 90 cm	Violet	PG 1475	
	3.5	0	1 x 90 cm	Violet	PG 1476	
	4	1	1 x 90 cm	Violet	PG 1477	
	5	2	1 x 90 cm	Violet	PG 1478	
 <p>HR 50</p>	3.5	0	1 x 75 cm	Violet	PG 1402	
	5	2	1 x 75 cm	Violet	PG 1400	
	3.5	0	1 x 90 cm	Violet	PG 1403	
 <p>HR 50</p>	4	1	1 x 90 cm	Violet	PG 1401	
	5	2	1 x 90 cm	Violet	PG 1399	
	3	2/0	1 x 75 cm	Violet	PG 1291	
 <p>HR 60</p>	3	2/0	1 x 90 cm	Violet	PG 1391	
	3.5	0	1 x 90 cm	Violet	PG 1392	
	4	1	1 x 90 cm	Violet	PG 1393	
 <p>HR 65</p>	3	2/0	1 x 75 cm	Violet	PG 0278	
	3.5	0	1 x 75 cm	Violet	PG 0279	
 <p>HR 75</p>	5	2	1 x 75 cm	Violet	PG 1286	
	4	1	1 x 75 cm	Violet	PG 1289	
 <p>HR 75</p>	5	2	1 x 75 cm	Violet	PG 0496	

CHIRLAC

UK03
106



Absorbable - Chirlac

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR s: ½ circle, round body needle - slim						
HR s 18 	1.5 2	4/0 3/0	1 x 75 cm 1 x 75 cm	Violet Violet	PG 1256 PG 1257	Slim needle
HR s 18 black 	1.5 2	4/0 3/0	1 x 75 cm 1 x 75 cm	Violet Violet	PG 2256 PG 2257	Slim, black needle
HR s 22 	2	3/0	1 x 45 cm	Violet	PG 1362	Slim needle
	2 3	3/0 2/0	1 x 75 cm 1 x 75 cm	Violet Violet	PG 1262 PG 1263	
HR s 22 black 	2 3	3/0 2/0	1 x 75 cm 1 x 75 cm	Violet Violet	PG 2262 PG 2263	Slim, black needle
HR s 27 	2	3/0	1 x 45 cm	Violet	PG 1366	Slim needle
	3	2/0	1 x 75 cm	Violet	PG 1267	



CHIRLAC

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRX: ½ circle, round body needle - extra strong						
 HRX 22	2	3/0	1 x 75 cm	Violet	PG 0420	
	3	2/0	1 x 75 cm	Violet	PG 0421	
	3.5	0	1 x 75 cm	Violet	PG 0422	
 HRX 27	2	3/0	1 x 75 cm	Violet	PG 0429	
	3	2/0	1 x 75 cm	Violet	PG 0430	
	3.5	0	1 x 75 cm	Violet	PG 0431	
	4	1	1 x 75 cm	Violet	PG 0432	
	5	2	1 x 75 cm	Violet	PG 0433	
	3	2/0	1 x 90 cm	Violet	PG 0435	
3.5	0	1 x 90 cm	Violet	PG 0436		
4	1	1 x 90 cm	Violet	PG 0437		
 HRX 30	4	1	1 x 75 cm	Violet	PG 1436	
 HRX 35	3	2/0	1 x 75 cm	Violet	PG 0280	
	3.5	0	1 x 75 cm	Violet	PG 0281	
	4	1	1 x 75 cm	Violet	PG 0282	
 HRX 35	3	2/0	1 x 90 cm	Violet	PG 0380	
	3.5	0	1 x 90 cm	Violet	PG 0381	
	4	1	1 x 90 cm	Violet	PG 0382	
 HRX 36	3.5	0	1 x 90 cm	Violet	PG 1381	
	4	1	1 x 90 cm	Violet	PG 1382	
 HRX 40	3.5	0	1 x 75 cm	Violet	PG 0285	
	4	1	1 x 75 cm	Violet	PG 0286	
	5	2	1 x 75 cm	Violet	PG 0287	
 HRX 40	3.5	0	1 x 90 cm	Violet	PG 0385	
	4	1	1 x 90 cm	Violet	PG 0386	
	5	2	1 x 90 cm	Violet	PG 0387	
 HRX 40	3.5	0	1 x 150 cm	Violet	PG 4035	Loop
	4	1	1 x 150 cm	Violet	PG 4040	Loop
 HRX 45	4	1	1 x 75 cm	Violet	PG 0290	
	5	2	1 x 75 cm	Violet	PG 0291	
 HRX 45	4	1	1 x 90 cm	Violet	PG 0390	
	5	2	1 x 90 cm	Violet	PG 0391	



CHIRLAC

UK03
108

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRX: ½ circle, round body needle - extra strong						
 <p>HRX 50</p>	3.5	0	1 x 75 cm	Violet	PG 0294	
	4	1	1 x 75 cm	Violet	PG 0295	
	5	2	1 x 75 cm	Violet	PG 0296	
 <p>HRX 50</p>	4	1	1 x 90 cm	Violet	PG 0395	
	5	2	1 x 90 cm	Violet	PG 0396	
	5	2	1 x 110 cm	Violet	PG 0393	
	4	1	1 x 150 cm	Violet	PG 4540	Loop
 <p>HRX 70</p>	5	2	1 x 240 cm	Violet	PG 6550	Loop
HRT: ½ circle, round body needle - tapercut						
 <p>HRT 13</p>	1.5	4/0	1 x 75 cm	Violet	PG 0551	
	1.5	4/0	1 x 45 cm	Violet	PG 0566	
	2	3/0	1 x 45 cm	Violet	PG 0567	
 <p>HRT 17</p>	1	5/0	1 x 75 cm	Violet	PG 0560	
	1.5	4/0	1 x 75 cm	Violet	PG 0561	
	2	3/0	1 x 75 cm	Violet	PG 0562	
	3	2/0	1 x 90 cm	Violet	PG 0569	
 <p>HRT 17 A</p>	1.5	4/0	6 x 45 cm	Violet	PG 0633	Take-off needle
 <p>HRT 22</p>	3	2/0	1 x 75 cm	Violet	PG 0573	
 <p>HRT 25</p>	1.5	4/0	1 x 75 cm	Violet	PG 0580	
	2	3/0	1 x 75 cm	Violet	PG 0581	
	3	2/0	1 x 75 cm	Violet	PG 0582	

CHIRLAC

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRT: ½ circle, round body needle - tapercut						
HRT 30	3.5	0	1 x 75 cm	Violet	PG 0586	
	4	1	1 x 90 cm	Violet	PG 0589	
HRT 31	3	2/0	1 x 75 cm	Violet	PG 0606	
						
HRT 37	3	2/0	1 x 75 cm	Violet	PG 0592	
	3.5	0	1 x 75 cm	Violet	PG 0593	
	3	2/0	1 x 90 cm	Violet	PG 0597	
	3.5	0	1 x 90 cm	Violet	PG 0598	
	4	1	1 x 90 cm	Violet	PG 0599	
HRT 45	5	2	1 x 75 cm	Violet	PG 0665	
						

UK03
109



Absorbable - Chirlac



CHIRLAC

UK03
110



Absorbable - Chirlac

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRTX: ½ circle, round body needle - tapercut, extra strong						
HRTX 45 	5	2	1 x 90 cm	Violet	PG 0643	
HRTX 50 	4	1	1 x 90 cm	Violet	PG 0656	
	4	1	1 x 150 cm	Violet	PG 5040	Loop
HRN: ½ circle, round body needle - blunt point						
HRN 50 	4	1	1 x 90 cm	Violet	PG 0686	
HRF: ½ circle, flattened body needle, taper point						
HRF 26 	2 3	3/0 2/0	1 x 75 cm 1 x 75 cm	Violet Violet	PG 9266 PG 9267	



CHIRLAC

UK03
111



Absorbable - Chirlac

Needle	EP	USP	Length	Colour	Cat. No	Other information
DR: 3/8 circle, round body needle						
DR 21 	3	2/0	1 x 75 cm	Violet	PG 0744	
DR 30 	2	3/0	1 x 75 cm	Violet	PG 0746	
DRN: 3/8 circle, round body needle - blunt point						
DRN 60 	4	1	1 x 100 cm	Violet	PG 0699	
FR: 5/8 circle, round body needle						
FR 27 	3 4	2/0 1	1 x 75 cm 1 x 75 cm	Violet Violet	PG 0610 PG 0612	
FR 38 	3 3.5	2/0 0	1 x 75 cm 1 x 75 cm	Violet Violet	PG 0620 PG 0621	
GR: Straight, round body needle						
GR 20 	1.5	4/0	1 x 45 cm	Violet	PG 0406	
	2	3/0	1 x 45 cm	Violet	PG 0407	
	1.5	4/0	1 x 75 cm	Violet	PG 0402	
	3	2/0	1 x 90 cm	Violet	PG 0399	

CHIRLAC

UK03
112




Absorbable - Chirlac

Needle	EP	USP	Length	Colour	Cat. No	Other information
--------	----	-----	--------	--------	---------	-------------------

ART: Asymmetric, round body needle - tapercut

ART 30 black



1.5	4/0	1 x 45 cm	Violet	PG 0630	Black needle
-----	-----	-----------	--------	---------	--------------

KR: Ski, round body needle

KR 27



2	3/0	1 x 40 cm	Violet	PG 4630	Loop
---	-----	-----------	--------	---------	------

JRT: Hook shaped, round body needle - tapercut

JRT 30



4	1	1 x 75 cm	Violet	PG 0694	
---	---	-----------	--------	---------	--

J-RX: Hook shaped, round body needle

J-RX 22



3 3.5	2/0 0	1 x 75 cm 1 x 75 cm	Violet Violet	PG 0397 PG 0398	
----------	----------	------------------------	------------------	--------------------	--



CHIRLAC



Needle	EP	USP	Length	Colour	Cat. No	Other information
HS: ½ circle, reverse cutting needle						
 HS 18	2	3/0	1 x 75 cm	Violet	PG 0424	
	3	2/0	1 x 75 cm	Violet	PG 0425	
 HS 21	1.5	4/0	1 x 45 cm	Violet	PG 0298	
	2	3/0	1 x 45 cm	Violet	PG 0299	
 2 x HS 21	1.5	4/0	1 x 75 cm	Violet	PG 0300	
	2	3/0	1 x 75 cm	Violet	PG 0301	
 2 x HS 21	1.5	4/0	1 x 75 cm	Violet	PG 1300	
 HS 25	2	3/0	1 x 75 cm	Violet	PG 0441	
	3	2/0	1 x 75 cm	Violet	PG 0442	
	3.5	0	1 x 75 cm	Violet	PG 0443	
 HS 30	3	2/0	1 x 75 cm	Violet	PG 0445	
	3.5	0	1 x 75 cm	Violet	PG 0446	
	4	1	1 x 75 cm	Violet	PG 0447	
 HS 36	3.5	0	1 x 75 cm	Violet	PG 0351	
 HS 37	2	3/0	1 x 90 cm	Violet	PG 0449	
 HS 40	3	2/0	1 x 75 cm	Violet	PG 0455	
	3.5	0	1 x 75 cm	Violet	PG 0456	
	4	1	1 x 75 cm	Violet	PG 0457	
	5	2	1 x 75 cm	Violet	PG 0458	
	3.5	0	1 x 90 cm	Violet	PG 0460	
	4	1	1 x 90 cm	Violet	PG 0461	

CHIRLAC

UK03
114

Needle	EP	USP	Length	Colour	Cat. No	Other information
HS: ½ circle, reverse cutting needle						



Absorbable - Chirlac

HS 45 	3.5	0	1 x 75 cm	Violet	PG 0472	
HS 48 	3.5	0	1 x 75 cm	Violet	PG 0466	
	3.5	0	1 x 90 cm	Violet	PG 0468	
	4	1	1 x 90 cm	Violet	PG 0469	
	5	2	1 x 90 cm	Violet	PG 0470	
HS 55 	4	1	1 x 75 cm	Violet	PG 0476	
	5	2	1 x 75 cm	Violet	PG 0477	

DS: 3/8 circle, reverse cutting needle						
---	--	--	--	--	--	--

DS 12 	0.7	6/0	1 x 45 cm	Violet	PG 0201	
	1	5/0	1 x 45 cm	Violet	PG 0202	
	1.5	4/0	1 x 45 cm	Violet	PG 0203	
	0.7	6/0	1 x 75 cm	Violet	PG 0196	
DS 15 	1	5/0	1 x 45 cm	Violet	PG 0205	
	1.5	4/0	1 x 45 cm	Violet	PG 0206	
	2	3/0	1 x 45 cm	Violet	PG 0207	
	1.5	4/0	1 x 75 cm	Violet	PG 0209	
DS 19 	1	5/0	1 x 45 cm	Violet	PG 0214	
	1.5	4/0	1 x 45 cm	Violet	PG 0215	
	2	3/0	1 x 45 cm	Violet	PG 0216	
	3	2/0	1 x 45 cm	Violet	PG 0217	
	1.5	4/0	1 x 75 cm	Violet	PG 0211	
	2	3/0	1 x 75 cm	Violet	PG 0212	
	3	2/0	1 x 75 cm	Violet	PG 0213	
DS 24 	2	3/0	1 x 75 cm	Violet	PG 1225	
	3	2/0	1 x 75 cm	Violet	PG 0228	
DS 25 	1.5	4/0	1 x 75 cm	Violet	PG 0222	
	2	3/0	1 x 75 cm	Violet	PG 0223	
	3	2/0	1 x 75 cm	Violet	PG 0224	

CHIRLAC

UK03
115

Absorbable - Chirlac

Needle	EP	USP	Length	Colour	Cat. No	Other information
--------	----	-----	--------	--------	---------	-------------------

DS: 3/8 circle, reverse cutting needle

 <p>DS 26</p>	2	3/0	1 x 45 cm	Violet	PG 1221	
 <p>DS 30</p>	3 3.5 4	2/0 0 1	1 x 75 cm 1 x 75 cm 1 x 75 cm	Violet Violet Violet	PG 0234 PG 0235 PG 0229	
 <p>DS 37</p>	5	2	1 x 75 cm	Violet	PG 1237	
 <p>DS 38</p>	3 4	2/0 1	1 x 75 cm 1 x 75 cm	Violet Violet	PG 0244 PG 1246	
 <p>DS 50</p>	3.5	0	1 x 75 cm	Violet	PG 1247	

DSM: 3/8 circle, reverse cutting needle - micropoint

 <p>DSM 12</p>	1	5/0	1 x 45 cm	Violet	PG 0302	Advanced cutting needle
 <p>DSM 19</p>	2	3/0	1 x 45 cm	Violet	PG 1213	Advanced cutting needle
 <p>DSM 25</p>	2	3/0	1 x 75 cm	Violet	PG 1223	Advanced cutting needle

GS: Straight, reverse cutting needle

 <p>GS 13</p>	1.5	4/0	1 x 45 cm	Violet	PG 1198	
 <p>GS 60</p>	2 3	3/0 2/0	1 x 75 cm 1 x 75 cm	Violet Violet	PG 0236 PG 0237	

CHIRLAC natural ligatures

UK03
116

Pre-cut threads

EP	USP	Length	Cat. No
1	5/0	3 x 45 cm	PN 0012
1.5	4/0	3 x 45 cm	PN 0013
2	3/0	3 x 45 cm	PN 0014
3	2/0	3 x 45 cm	PN 0015
3.5	0	3 x 45 cm	PN 0016
2	3/0	10 x 45 cm	PN 0064
3	2/0	10 x 45 cm	PN 0065
1.5	4/0	2 x 75 cm	PN 0023
2	3/0	2 x 75 cm	PN 0024
3	2/0	2 x 75 cm	PN 0025

Chirapack - 20 pcs./box

EP	USP	Length	Cat. No
1.5	4/0	1 x 300 cm	PN 0003
2	3/0	1 x 300 cm	PN 0004
3	2/0	1 x 300 cm	PN 0005
3.5	0	1 x 300 cm	PN 0006



Absorbable - Chirlac



CHIRLAC

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
HR 13 	1 1.5	5/0 4/0	1 x 45 cm 1 x 45 cm	Natural Natural	PN 0252 PN 0253	
HR 13 A 	1	5/0	6 x 45 cm	Natural	PN 0352	Take-off needle
HR 18 	1.5 2	4/0 3/0	1 x 75 cm 1 x 75 cm	Natural Natural	PN 0256 PN 0257	
HR 18 A 	2	3/0	6 x 45 cm	Natural	PN 0357	Take-off needle
HR 22 	2 3	3/0 2/0	1 x 75 cm 1 x 75 cm	Natural Natural	PN 0262 PN 0263	
HR 27 	2 3 3.5	3/0 2/0 0	1 x 75 cm 1 x 75 cm 1 x 75 cm	Natural Natural Natural	PN 0266 PN 0267 PN 0268	
HR 30 	3.5	0	1 x 75 cm	Natural	PN 0274	

CHIRLAC

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRT: ½ circle, round body needle - tapercut						
HRT 17 	1.5 2	4/0 3/0	1 x 75 cm 1 x 75 cm	Natural Natural	PN 0561 PN 0562	
HRT 22 	2	3/0	1 x 75 cm	Natural	PN 0572	
HRF: ½ circle, flattened body needle, taper point						
HRF 13 	0.7	6/0	1 x 45 cm	Natural	PN 9251	
HS: ½ circle, reverse cutting needle						
HS 18 	2	3/0	1 x 75 cm	Natural	PN 0424	
DS: 3/8 circle, reverse cutting needle						
DS 12 	0.7 1	6/0 5/0	1 x 45 cm 1 x 45 cm	Natural Natural	PN 1201 PN 1202	
DS 15 	1.5 2	4/0 3/0	1 x 45 cm 1 x 45 cm	Natural Natural	PN 1206 PN 1207	
DS 19 	1 1.5 2	5/0 4/0 3/0	1 x 45 cm 1 x 45 cm 1 x 45 cm	Natural Natural Natural	PN 1215 PN 1216 PN 1217	
DS 25 	2	3/0	1 x 45 cm	Natural	PN 1332	
DS 26 	3	2/0	1 x 75 cm	Natural	PN 1334	
DS 30 	3 3.5	2/0 0	1 x 75 cm 1 x 75 cm	Natural Natural	PN 1243 PN 1244	
DS 36 	3.5	0	1 x 75 cm	Natural	PN 1245	

CHIRLAC

UK03
118

Absorbable - Chirlac

Needle	EP	USP	Length	Colour	Cat. No	Other information
DSM: 3/8 circle, reverse cutting needle - micropoint						
DSM 12 	0.7	6/0	1 x 45 cm	Natural	PN 0201	Advanced cutting needle
	1	5/0	1 x 45 cm	Natural	PN 0202	
	1.5	4/0	1 x 45 cm	Natural	PN 0203	
DSM 15 	1	5/0	1 x 45 cm	Natural	PN 0206	Advanced cutting needle
	1.5	4/0	1 x 45 cm	Natural	PN 0207	
	2	3/0	1 x 45 cm	Natural	PN 0208	
DSM 19 	1.5	4/0	1 x 45 cm	Natural	PN 0216	Advanced cutting needle
	2	3/0	1 x 45 cm	Natural	PN 0217	
	1.5	4/0	1 x 75 cm	Natural	PN 0212	
		3/0	1 x 75 cm	Natural	PN 0213	
DSM 25 	2	3/0	1 x 75 cm	Natural	PN 0223	Advanced cutting needle
	3	2/0	1 x 75 cm	Natural	PN 0224	
DSMC: 3/8 circle, conventional cutting needle - micropoint						
DSMC 16 	1.5	4/0	1 x 45 cm	Natural	PN 8207	Advanced cutting needle

Polydiox

Polydiox

Polydiox

P

O

I



W

olydiox

X

X

X



POLYDOX

POLYDOX is a long-term absorbable, synthetic, monofilament suture made of polydioxanone. The suture offers good handling combined with great strength due to its slow degradation profile.

Polydox is ideal for all indications, where an extended wound support of more than 4 weeks is desirable. The threads are available in violet and natural.

The suture is very flexible, pliable and easy to knot. In tissue, the thread can cause a slight initial inflammatory reaction.



Indication: General soft tissue approximation, Orthopaedics (e.g. Joint capsule, Tendons), Fascia Closure, Urology, Gynaecology and Plastic surgery.

Contraindication: Polydox should not be used where prolonged approximation of tissues under stress is required.

Sizes: 7/0 to 4 USP

Structure: Monofilament

Absorption method: Hydrolysis

Strength Retention Profile:

LONG-TERM HEALING (90 - 98 days)	
1 Day	100%
45 Days	50%

Mass Absorption profile: 180 - 240 days

POLYDOX

UK03
122



Absorbable - Polydiox

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
HR 10 	0.7	6/0	1 x 75 cm	Violet	DX 0246	
	1	5/0	1 x 75 cm	Violet	DX 0247	
HR 13 	1	5/0	1 x 45 cm	Violet	DX 0248	
	0.7	6/0	1 x 75 cm	Violet	DX 0249	
	1	5/0	1 x 75 cm	Violet	DX 0250	
2 x HR 13 	1.5	4/0	1 x 75 cm	Violet	DX 0251	
	1	5/0	1 x 75 cm	Violet	DX 5245	
HR 17 	1	5/0	1 x 75 cm	Violet	DX 0254	
HR 18 	1	5/0	1 x 75 cm	Violet	DX 0257	
	1.5	4/0	1 x 75 cm	Violet	DX 0258	
	2	3/0	1 x 75 cm	Violet	DX 0259	
2 x HR 18 	1	5/0	1 x 75 cm	Violet	DX 5250	
	1.5	4/0	1 x 90 cm	Violet	DX 5251	
HR 20 	1	5/0	1 x 75 cm	Violet	DX 0360	
	1.5	4/0	1 x 75 cm	Violet	DX 0256	
HR 22 	1.5	4/0	1 x 75 cm	Violet	DX 0261	
	2	3/0	1 x 75 cm	Violet	DX 0262	
	3	2/0	1 x 75 cm	Violet	DX 0263	
2 x HR 22 	1.5	4/0	1 x 75 cm	Violet	DX 5261	
	2	3/0	1 x 75 cm	Violet	DX 5262	
2 x HR 22 	2	3/0	1 x 90 cm	Violet	DX 5264	
	HR 26 	2	3/0	1 x 75 cm	Violet	DX 0366
3		2/0	1 x 75 cm	Violet	DX 0365	
2 x HR 26 	3	2/0	1 x 75 cm	Violet	DX 5269	
	1.5	4/0	1 x 90 cm	Violet	DX 5271	
2 x HR 26 	2	3/0	1 x 90 cm	Violet	DX 5270	

POLYDOX

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
 HR 27	1.5	4/0	1 x 75 cm	Violet	DX 0265	
	2	3/0	1 x 75 cm	Violet	DX 0266	
	3	2/0	1 x 75 cm	Violet	DX 0267	
	3.5	0	1 x 75 cm	Violet	DX 0268	
	4	1	1 x 75 cm	Violet	DX 0369	
 2 x HR 27	1.5	4/0	1 x 75 cm	Violet	DX 5265	
	2	3/0	1 x 75 cm	Violet	DX 5266	
	3	2/0	1 x 75 cm	Violet	DX 5267	
 HR 27	1.5	4/0	1 x 90 cm	Violet	DX 5295	
	2	3/0	1 x 120 cm	Violet	DX 5297	
 HR 30	1.5	4/0	1 x 75 cm	Violet	DX 0269	
	2	3/0	1 x 75 cm	Violet	DX 0270	
	3	2/0	1 x 75 cm	Violet	DX 0271	
	3.5	0	1 x 75 cm	Violet	DX 0272	
	4	1	1 x 75 cm	Violet	DX 0273	
 2 x HR 30	2	3/0	1 x 90 cm	Violet	DX 5275	
 HR 37	2	3/0	1 x 75 cm	Violet	DX 0264	
	3.5	0	1 x 75 cm	Violet	DX 0356	
 HR 40	3	2/0	1 x 75 cm	Violet	DX 0275	
	3.5	0	1 x 75 cm	Violet	DX 0276	
	4	1	1 x 75 cm	Violet	DX 0277	
	4	1	1 x 90 cm	Violet	DX 0377	
 HR 40	4	1	1 x 150 cm	Violet	DX 0279	
 HR 45	5	2	1 x 90 cm	Violet	DX 0378	



POLYDOX

UK03
124

Needle	EP	USP	Length	Colour	Cat. No	Other information
--------	----	-----	--------	--------	---------	-------------------

HR: ½ circle, round body needle

	HR 50	4	1	1 x 90 cm	Violet	DX 0380	
	HR 65	4	1	1 x 150 cm	Violet	DX 4547	Loop
		4	1	1 x 250 cm	Violet	DX 4546	Loop

HR s: ½ circle, round body needle - slim

	HR s 18	1.5	4/0	1 x 75 cm	Violet	DX 1256	Slim needle
	HR s 22	2	3/0	1 x 75 cm	Violet	DX 1262	Slim needle
	HR s 27	2	3/0	1 x 75 cm	Violet	DX 1272	Slim needle
		3	2/0	1 x 75 cm	Violet	DX 1273	


HRX: ½ circle, round body needle - extra strong

	HRX 27	3	2/0	1 x 75 cm	Violet	DX 0430	
		3.5	0	1 x 75 cm	Violet	DX 0431	
		4	1	1 x 75 cm	Violet	DX 0432	
	HRX 35	3.5	0	1 x 75 cm	Violet	DX 0281	
		4	1	1 x 75 cm	Violet	DX 0282	
		3.5	0	1 x 90 cm	Violet	DX 0381	
	HRX 37	4	1	1 x 90 cm	Violet	DX 0384	

POLYDOX

UK03
125

Absorbable - Polydox

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRX: ½ circle, round body needle - extra strong						
 <p>HRX 40</p>	3.5	0	1 x 75 cm	Violet	DX 0285	
	3.5	0	1 x 90 cm	Violet	DX 0385	
	4	1	1 x 90 cm	Violet	DX 0386	
	5	2	1 x 90 cm	Violet	DX 0387	
	3.5	0	1 x 150 cm	Violet	DX 4035	Loop
	4	1	1 x 150 cm	Violet	DX 4040	Loop
	4	1	1 x 240 cm	Violet	DX 4024	Loop
 <p>HRX 45</p>	4	1	1 x 90 cm	Violet	DX 0390	
	5	2	1 x 90 cm	Violet	DX 0391	
 <p>HRX 48</p>	3.5	0	1 x 90 cm	Violet	DX 0400	
	3.5	0	1 x 150 cm	Violet	DX 4135	Loop
	4	1	1 x 150 cm	Violet	DX 4140	Loop
 <p>HRX 50</p>	3.5	0	1 x 90 cm	Violet	DX 0394	
	4	1	1 x 90 cm	Violet	DX 0395	
	5	2	1 x 90 cm	Violet	DX 0396	
	3.5	0	1 x 150 cm	Violet	DX 4535	Loop
	4	1	1 x 150 cm	Violet	DX 1392	
	4	1	1 x 150 cm	Violet	DX 4540	Loop
	5	2	1 x 150 cm	Violet	DX 4550	Loop
	4	1	1 x 200 cm	Violet	DX 4542	Loop
 <p>HRX 70</p>	4	1	1 x 150 cm	Violet	DX 4543	Loop

POLYDOX

UK03
126



Absorbable - Polydiox

Needle	EP	USP	Length	Colour	Cat. No	Other information
--------	----	-----	--------	--------	---------	-------------------

HRT: ½ circle, round body needle - tapercut

HRT 25



2

3/0

1 x 75 cm

Violet

DX 0581

HRT 26



3

2/0

1 x 75 cm

Violet

DX 0572

HRT 37



4

1

1 x 90 cm

Violet

DX 0599

HRT 48



4

1

1 x 90 cm

Violet

DX 0608

HRTX: ½ circle, round body needle – tapercut, extra strong

HRTX 40



4

1

1 x 90 cm

Violet

DX 0610

HRF: ½ circle, flattened body needle, taper point

HRF 26



2
3
3.5

3/0
2/0
0

1 x 75 cm
1 x 75 cm
1 x 75 cm

Violet
Violet
Violet

DX 9266
DX 9267
DX 9268

HRF 31



3

2/0

1 x 75 cm

Violet

DX 9271

POLYDOX

UK03
127

Absorbable - Polydiox

Needle	EP	USP	Length	Colour	Cat. No	Other information
DR: 3/8 circle, round body needle						
2 x DR 10 	0.7	6/0	1 x 75 cm	Violet	DX 1120	
DR 12 	1 1.5	5/0 4/0	1 x 75 cm 1 x 75 cm	Violet Violet	DX 0120 DX 0121	
DR 17 	1.5	4/0	1 x 75 cm	Violet	DX 0132	
DR 20 	1.5 3	4/0 2/0	1 x 75 cm 1 x 75 cm	Violet Violet	DX 0141 DX 0143	
DRT: 3/8 circle, round body needle - tapercut						
DRT 13 	1	5/0	1 x 75 cm	Violet	DX 0521	
FR: 5/8 circle, round body needle						
2 x FR 27 	3	2/0	1 x 75 cm	Violet	DX 5367	
FR 65 	3.5	0	1 x 150 cm	Violet	DX 4635	Loop
GR: Straight, round body needle						
2 x GR 60 	2	3/0	1 x 75 cm	Violet	DX 5396	
KR: Ski, round body needle						
KR 27 	2	3/0	1 x 40 cm	Violet	DX 4630	Loop
JRT: Hook shaped, round body needle - tapercut						
JRT 30 	4	1	1 x 75 cm	Violet	DX 0694	

POLYDOX

UK03
128



Absorbable - Polydiox

Needle	EP	USP	Length	Colour	Cat. No	Other information
HS: ½ circle, reverse cutting needle						
<p>HS 18</p>	1.5	4/0	1 x 75 cm	Violet	DX 0155	
<p>HS 21</p>	2	3/0	1 x 75 cm	Violet	DX 0166	
<p>HS 22</p>	2 3	3/0 2/0	1 x 45 cm 1 x 45 cm	Violet Violet	DX 0167 DX 0170	
<p>HS 27</p>	3	2/0	1 x 75 cm	Violet	DX 0168	
<p>HS 30</p>	2 3	3/0 2/0	1 x 75 cm 1 x 75 cm	Violet Violet	DX 0171 DX 0172	
<p>HS 37</p>	3 3.5 4	2/0 0 1	1 x 75 cm 1 x 75 cm 1 x 75 cm	Violet Violet Violet	DX 0177 DX 0178 DX 0179	
<p>HS 37</p>	3.5	0	1 x 90 cm	Violet	DX 1178	
<p>HS 40</p>	4	1	1 x 75 cm	Violet	DX 0180	
<p>HS 48</p>	4	1	1 x 90 cm	Violet	DX 0184	

POLYDOX

Needle	EP	USP	Length	Colour	Cat. No	Other information	
DS: 3/8 circle, reverse cutting needle							
DS 15 	1	5/0	1 x 45 cm	Violet	DX 0190		
	1.5 2	4/0 3/0	1 x 75 cm 1 x 75 cm	Violet Violet	DX 0196 DX 0197		
DS 16 	1	5/0	1 x 45 cm	Violet	DX 0198		
	1 1.5 2	5/0 4/0 3/0	1 x 45 cm 1 x 45 cm 1 x 45 cm	Violet Violet Violet	DX 0210 DX 0211 DX 0212		
DS 19 	1.5 2	4/0 3/0	1 x 75 cm 1 x 75 cm	Violet Violet	DX 0215 DX 0216		
	2 3	3/0 2/0	1 x 45 cm 1 x 45 cm	Violet Violet	DX 0222 DX 0223		
DS 25 	2 3 3.5 4	3/0 2/0 0 1	1 x 75 cm 1 x 75 cm 1 x 75 cm 1 x 75 cm	Violet Violet Violet Violet	DX 0227 DX 0228 DX 0229 DX 0230		
	DS 30 	3 3.5 4	2/0 0 1	1 x 75 cm 1 x 75 cm 1 x 75 cm	Violet Violet Violet	DX 0232 DX 0234 DX 0235	
		DSC: 3/8 circle, conventional cutting needle					
		DSC 25 	3	2/0	1 x 75 cm	Violet	DX 8228

UK03
129

Absorbable - Polydiox

POLYDOX

UK03
130



Absorbable - Polydiox

Needle	EP	USP	Length	Colour	Cat. No	Other information
DSM: 3/8 circle, reverse cutting needle - micropoint						
DSM 12 	0.7 1	6/0 5/0	1 x 45 cm 1 x 45 cm	Violet Violet	DX 0201 DX 0202	Advanced cutting needle
DSM 15 	1.5	4/0	1 x 45 cm	Violet	DX 0206	Advanced cutting needle
DSM 19 	1.5 2	4/0 3/0	1 x 75 cm 1 x 75 cm	Violet Violet	DX 0207 DX 0208	Advanced cutting needle
GS: Straight, reverse cutting needle						
2 x GS 65 	2	3/0	1 x 75 cm	Violet	DX 3242	

Non-absorbable

Chiraflon

Chiralen

Silon monofilament

Tervalon

Silon braided

Orsilon braided

Silk

Steel Wire

Non-absorbable

Non-absorbable



Chiraflon

Chiraflon

chiraflon

Chiraflon



CHIRAFLON

CHIRAFLON is a non-absorbable, monofilament suture made of PVDF (similar to teflon). Physiologically this thread is completely inert.

Chiraflon offers optimal smoothness and superior knotting properties. The thread is characterised by its high tensile strength. No splitting occurs even under considerable strain.



Indication: All surgical specialities, especially Cardiovascular, Plastic, Reconstructive, Ophthalmic, Orthopaedic and Traumatology (e.g. tendons).

Contraindication: Chiraflon must not penetrate into the inner parts of the digestive and urinary tracts.

Chemical Composition: Polyvinylidene Fluoride (PVDF)

Sizes: 8/0 to 3 USP

Structure: Monofilament

Absorption method: Non-absorbable

Colour: Blue



CHIRAFLON

UK03
136



Non-absorbable - Chiraflon

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
HR 13 	1	5/0	1 x 45 cm	Blue	TF 6202	
HR 18 	1	5/0	1 x 45 cm	Blue	TF 6211	
	1.5	4/0	1 x 45 cm	Blue	TF 6214	
2 x HR 18 	1	5/0	1 x 75 cm	Blue	TF 6216	
	1.5	4/0	1 x 75 cm	Blue	TF 6217	
HR 27 	3	2/0	1 x 75 cm	Blue	TF 6221	
HRX: ½ circle, round body needle - extra strong						
HRX 40 	3.5	0	1 x 90 cm	Blue	TF 6251	
HRX 50 	4	1	1 x 150 cm	Blue	TF 4540	Loop
	3.5	0	1 x 200 cm	Blue	TF 5035	Loop
	4	1	1 x 200 cm	Blue	TF 5040	Loop
HRT: ½ circle, round body needle - tapercut						
HRT 13 	0.7	6/0	1 x 45 cm	Blue	TF 6421	
2 x HRT 17 	1	5/0	1 x 75 cm	Blue	TF 6332	
DRT: 3/8 circle, round body needle - tapercut						
2 x DRT 10 	0.7	6/0	1 x 75 cm	Blue	TF 6301	
2 x DRT 13 	0.7	6/0	1 x 75 cm	Blue	TF 6311	
	1	5/0	1 x 75 cm	Blue	TF 6312	

CHIRAFLON

UK03
137



Non-absorbable - Chiraflon

Needle	EP	USP	Length	Colour	Cat. No	Other information
GR: Straight, round body needle						
GR 20 	1.5	4/0	1 x 45 cm	Blue	TF 6411	
	2	3/0	1 x 45 cm	Blue	TF 6412	
GR 60 	2	3/0	1 x 75 cm	Blue	TF 6462	
DS: 3/8 circle, reverse cutting needle						
DS 12 	0.7	6/0	1 x 45 cm	Blue	TF 6101	
	1	5/0	1 x 45 cm	Blue	TF 6102	
DS 15 	1	5/0	1 x 45 cm	Blue	TF 6105	
	1.5	4/0	1 x 45 cm	Blue	TF 6106	
DS 19 	1	5/0	1 x 45 cm	Blue	TF 6110	
	1.5	4/0	1 x 45 cm	Blue	TF 6115	
	2	3/0	1 x 45 cm	Blue	TF 6116	
	3	2/0	1 x 45 cm	Blue	TF 6117	
DS 25 	1.5	4/0	1 x 75 cm	Blue	TF 6111	
	2	3/0	1 x 75 cm	Blue	TF 6112	
	3	2/0	1 x 75 cm	Blue	TF 6113	
DS 25 	2	3/0	1 x 45 cm	Blue	TF 6126	
	3.5	0	1 x 45 cm	Blue	TF 6128	
DS 30 	2	3/0	1 x 75 cm	Blue	TF 6122	
	3	2/0	1 x 75 cm	Blue	TF 6123	
	3	2/0	1 x 75 cm	Blue	TF 6131	
DS 30 	3.5	0	1 x 75 cm	Blue	TF 6132	
	4	1	1 x 75 cm	Blue	TF 6133	
	4	1	1 x 75 cm	Blue	TF 6133	
DS 38 	2	3/0	1 x 75 cm	Blue	TF 6140	

CHIRAFLON

UK03
138

Needle	EP	USP	Length	Colour	Cat. No	Other information
DSM: 3/8 circle, reverse cutting needle - micropoint						
 DSM 12	0.7	6/0	1 x 45 cm	Blue	TF 6001	Advanced cutting needle
	1	5/0	1 x 45 cm	Blue	TF 6002	
	1.5	4/0	1 x 45 cm	Blue	TF 6003	
 DSM 15	1	5/0	1 x 45 cm	Blue	TF 6005	Advanced cutting needle
	1.5	4/0	1 x 45 cm	Blue	TF 6006	
 DSM 19	1	5/0	1 x 75 cm	Blue	TF 6008	Advanced cutting needle
	1.5	4/0	1 x 75 cm	Blue	TF 6011	
	2	3/0	1 x 75 cm	Blue	TF 6012	Advanced cutting needle
GS: Straight, reverse cutting needle						
 GS 60	3	2/0	1 x 75 cm	Blue	TF 6167	



Non-absorbable - Chiraflon









Chiralen

Chiralen

Chiralen

Chiralen

Chiralen



CHIRALEN

UK03
141

CHIRALEN is a non-absorbable, monofilament suture made of polypropylene. The thread is hydrophobic and is not subject to degradation by tissue enzymes.

Chiralen has a smooth, non-porous surface, that guarantees a clean and smooth passage through the tissue, optimal elasticity, good elongation properties and high visibility. In tissues, it can cause a minimal acute inflammatory reaction.



Indication: All surgical specialities, especially Cardiovascular, Plastic, Reconstructive, Ophthalmic, Orthopaedic and Traumatology (e.g. tendons).

Contraindication: Chiralen must not penetrate into the inner parts of the digestive and urinary tracts.

Chemical Composition: Polypropylene

Sizes: 10/0 to 2 USP

Structure: Monofilament

Absorption method: Non-absorbable

Colour: Blue and Black

CHIRALEN

UK03
142



Non-absorbable - Chiralen

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
2 x HR 8 	0.5	7/0	1 x 45 cm	Blue	PP 5205	
HR 10 	0.7	6/0	1 x 45 cm	Blue	PP 5142	
	0.7 1	6/0 5/0	1 x 45 cm 1 x 45 cm	Blue Blue	PP 5210 PP 5211	
2 x HR 10 	1	5/0	1 x 60 cm	Blue	PP 6211	
	0.5	7/0	1 x 75 cm	Blue	PP 5209	
HR 13 	0.7 1	6/0 5/0	1 x 75 cm 1 x 75 cm	Blue Blue	PP 5145 PP 5146	
	1	5/0	1 x 45 cm	Blue	PP 5220	
2 x HR 13 	0.7 1 1.5	6/0 5/0 4/0	1 x 75 cm 1 x 75 cm 1 x 75 cm	Blue Blue Blue	PP 5219 PP 5222 PP 5223	
2 x HR 15 	1.5	4/0	1 x 90 cm	Blue	PP 5228	
HR 17 	1	5/0	1 x 90 cm	Blue	PP 5558	
2 x HR 17 	1	5/0	1 x 75 cm	Blue	PP 5236	
	1 1.5	5/0 4/0	1 x 90 cm 1 x 90 cm	Blue Blue	PP 5248 PP 5244	
HR 18 	1 1.5 2	5/0 4/0 3/0	1 x 75 cm 1 x 75 cm 1 x 75 cm	Blue Blue Blue	PP 5152 PP 5153 PP 5154	
2 x HR 18 	1 1.5 2	5/0 4/0 3/0	1 x 75 cm 1 x 75 cm 1 x 75 cm	Blue Blue Blue	PP 5233 PP 5234 PP 5235	
	1 1.5 2	5/0 4/0 3/0	1 x 90 cm 1 x 90 cm 1 x 90 cm	Blue Blue Blue	PP 5237 PP 5238 PP 5239	

CHIRALEN

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
2 x HR 20 	1.5	4/0	1 x 90 cm	Blue	PP 5251	
HR 22 	1	5/0	1 x 75 cm	Blue	PP 5158	
	1.5	4/0	1 x 75 cm	Blue	PP 5157	
	2	3/0	1 x 75 cm	Blue	PP 5156	
	3	2/0	1 x 75 cm	Blue	PP 5155	
	1.5	4/0	1 x 90 cm	Blue	PP 5557	
	2	3/0	1 x 75 cm	Blue	PP 5249	
2 x HR 22 	3	2/0	1 x 75 cm	Blue	PP 5250	
	1.5	4/0	1 x 90 cm	Blue	PP 5246	
	2	3/0	1 x 90 cm	Blue	PP 5247	
	HR 25 	3	2/0	1 x 90 cm	Blue	PP 5559
HR 26 	3	2/0	1 x 75 cm	Blue	PP 5563	
	3.5	0	1 x 75 cm	Blue	PP 5562	
	3	2/0	1 x 90 cm	Blue	PP 5564	
2 x HR 26 	1.5	4/0	1 x 90 cm	Blue	PP 5351	
	2	3/0	1 x 90 cm	Blue	PP 5257	
	2	3/0	1 x 120 cm	Blue	PP 5362	
	2 x HR 26 black 	1.5	4/0	1 x 90 cm	Blue	PP 5258



CHIRALEN

UK03
144



Non-absorbable - Chiralen

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
HR 27 	2	3/0	1 x 75 cm	Blue	PP 5163	
	3	2/0	1 x 75 cm	Blue	PP 5164	
	3.5	0	1 x 90 cm	Blue	PP 5561	
	1.5	4/0	1 x 90 cm	Blue	PP 5266	
2 x HR 27 	2	3/0	1 x 90 cm	Blue	PP 5267	
	3	2/0	1 x 90 cm	Blue	PP 5255	
	2	3/0	1 x 120 cm	Blue	PP 5354	
HR 30 	3	2/0	1 x 75 cm	Blue	PP 5165	
	3.5	0	1 x 75 cm	Blue	PP 5166	
	4	1	1 x 75 cm	Blue	PP 5167	
	3	2/0	1 x 90 cm	Blue	PP 5261	
	3.5	0	1 x 100 cm	Blue	PP 5169	
2 x HR 30 	3	2/0	1 x 90 cm	Blue	PP 5271	
HR 37 	4	1	1 x 75 cm	Blue	PP 5181	
HR 40 	3	2/0	1 x 75 cm	Blue	PP 5185	
	4	1	1 x 100 cm	Blue	PP 5968	
HR s: ½ circle, round body needle - slim						
2 x HR s 18 Black 	1	5/0	1 x 75 cm	Blue	PP 5965	Slim black needles
2 x HR s 22 	1.5	4/0	1 x 75 cm	Blue	PP 5975	Slim needles

CHIRALEN

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRX: ½ circle, round body needle - extra strong						
HRX 35 	3.5 4	0 1	1 x 75 cm 1 x 75 cm	Blue Blue	PP 5373 PP 5372	
HRX 40 	4	1	1 x 90 cm	Blue	PP 5386	
2 x HRX 50 	3	2/0	1 x 120 cm	Blue	PP 5491	

UK03
145




Non-absorbable - Chiralen



CHIRALEN

UK03
146

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRT: ½ circle, round body needle - tapercut						
2 x HRT 10 	0.5	7/0	1 x 60 cm	Blue	PP 5336	
HRT 11 	1	5/0	1 x 60 cm	Blue	PP 5421	
HRT 13 	1	5/0	1 x 75 cm	Blue	PP 5422	
2 x HRT 13 	0.7 1	6/0 5/0	1 x 75 cm 1 x 75 cm	Blue Blue	PP 5321 PP 5322	
HRT 17 	1.5	4/0	1 x 75 cm	Blue	PP 5433	
2 x HRT 17 	1 1.5 2	5/0 4/0 3/0	1 x 75 cm 1 x 75 cm 1 x 75 cm	Blue Blue Blue	PP 5332 PP 5333 PP 5334	
	1 1.5 2	5/0 4/0 3/0	1 x 90 cm 1 x 90 cm 1 x 90 cm	Blue Blue Blue	PP 5932 PP 5933 PP 5934	



CHIRALEN

UK03
147



Non-absorbable - Chiralen

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRT: ½ circle, round body needle - tapercut						
HRT 25 	3	2/0	1 x 75 cm	Blue	PP 5445	
2 x HRT 25 	3	2/0	1 x 75 cm	Blue	PP 5345	
	2 3	3/0 2/0	1 x 90 cm 1 x 90 cm	Blue Blue	PP 5944 PP 5945	
2 x HRT 27 	2	3/0	1 x 90 cm	Blue	PP 5952	
HRT 37 	3.5	0	1 x 75 cm	Blue	PP 5466	
DR: 3/8 circle, round body needle						
2 x DR 8 	0.4	8/0	1 x 75 cm	Blue	PP 5282	
2 x DR 9.3 	0.5	7/0	1 x 60 cm	Blue	PP 5921	
2 x DR 9.3 Black 	0.7	6/0	1 x 60 cm	Blue	PP 5920	Black needles
DR 10 	0.5	7/0	1 x 75 cm	Blue	PP 5212	
	0.7	6/0	1 x 75 cm	Blue	PP 5213	
2 x DR 10 	0.7	6/0	1 x 45 cm	Blue	PP 5215	
	0.5	7/0	1 x 60 cm	Blue	PP 5913	
	0.7	6/0	1 x 60 cm	Blue	PP 5915	
DR 12 	0.7	6/0	1 x 75 cm	Blue	PP 5413	
	1.5	4/0	1 x 75 cm	Blue	PP 5415	

CHIRALEN

UK03
148



Non-absorbable - Chiralen

Needle	EP	USP	Length	Colour	Cat. No	Other information
DR: 3/8 circle, round body needle						
2 x DR 12 	0.5	7/0	1 x 60 cm	Blue	PP 5227	
	0.7	6/0	1 x 60 cm	Blue	PP 5229	
DR 13 	0.7	6/0	1 x 75 cm	Blue	PP 5224	
	1	5/0	1 x 75 cm	Blue	PP 5225	
DR 13 	0.7	6/0	1 x 75 cm	Blue	PP 5416	
2 x DR 13 	0.5	7/0	1 x 45 cm	Blue	PP 5404	
	0.7	6/0	1 x 60 cm	Blue	PP 5328	
	1	5/0	1 x 60 cm	Blue	PP 5403	
2 x DR 13 	0.7	6/0	1 x 75 cm	Blue	PP 5324	
	1	5/0	1 x 90 cm	Blue	PP 5331	
2 x DR 15 	1	5/0	1 x 75 cm	Blue	PP 5226	
DR 17 	1	5/0	1 x 75 cm	Blue	PP 5431	
DR 18 	1.5	4/0	1 x 90 cm	Blue	PP 5535	
2 x DR 18 	1	5/0	1 x 75 cm	Blue	PP 5923	
	1.5	4/0	1 x 75 cm	Blue	PP 5924	
2 x DR 18 	1.5	4/0	1 x 90 cm	Blue	PP 5928	
	2	3/0	1 x 75 cm	Blue	PP 5426	
2 x DR 30 	2	3/0	1 x 75 cm	Blue	PP 5426	



CHIRALEN

UK03
149



Non-absorbable - Chiralen

Needle	EP	USP	Length	Colour	Cat. No	Other information
DRT: 3/8 circle, round body needle - tapercut						
2 x DRT 10 	0.5	7/0	1 x 60 cm	Blue	PP 5298	
	0.7	6/0	1 x 75 cm	Blue	PP 5301	
2 x DRT 10 Black 	0.5	7/0	1 x 75 cm	Blue	PP 5902	Black needles
2 x DRT 11 	0.7	6/0	1 x 60 cm	Blue	PP 5296	
DRT 13 	0.7 1	6/0 5/0	1 x 75 cm 1 x 75 cm	Blue Blue	PP 5309 PP 5310	
2 x DRT 13 	0.7	6/0	1 x 60 cm	Blue	PP 5911	
	0.7 1	6/0 5/0	1 x 75 cm 1 x 75 cm	Blue Blue	PP 5311 PP 5312	
2 x DRT 17 	1	5/0	1 x 90 cm	Blue	PP 5314	
GR: Straight, round body needle						
GR 45  	3	2/0	1 x 45 cm	Blue	PP 5545	
2 x GR 60  	3	2/0	1 x 75 cm	Blue	PP 6245	
	3	2/0	1 x 100 cm	Blue	PP 5575	

CHIRALEN

UK03
150

Needle	EP	USP	Length	Colour	Cat. No	Other information
--------	----	-----	--------	--------	---------	-------------------

HS: ½ circle, reverse cutting needle

2 x HS 25



3

2/0

1 x 90 cm

Blue

PP 6445

HS 40



3.5

0

1 x 75 cm

Blue

PP 6845

HS 50



3.5

0

1 x 75 cm

Blue

PP 6855

4

1

1 x 90 cm

Blue

PP 6955

DS: 3/8 circle, reverse cutting needle

DS 12



0.7
1

6/0
5/0

1 x 45 cm
1 x 45 cm

Blue
Blue

PP 5101
PP 5102

DS 15



0.7
1
1.5

6/0
5/0
4/0

1 x 45 cm
1 x 45 cm
1 x 45 cm

Blue
Blue
Blue

PP 5105
PP 5106
PP 5107

DS 19



1
1.5
2

5/0
4/0
3/0

1 x 75 cm
1 x 75 cm
1 x 75 cm

Blue
Blue
Blue

PP 4108
PP 4109
PP 4110

DS 24



1
1.5
2

5/0
4/0
3/0

1 x 45 cm
1 x 45 cm
1 x 45 cm

Blue
Blue
Blue

PP 5509
PP 5510
PP 5511

1
1.5
2
3
3.5

5/0
4/0
3/0
2/0
0

1 x 75 cm
1 x 75 cm
1 x 75 cm
1 x 75 cm
1 x 75 cm

Blue
Blue
Blue
Blue
Blue

PP 5109
PP 5110
PP 5111
PP 5112
PP 5113

2

3/0

1 x 90 cm

Blue

PP 5519

DS 25



1.5
2
3

4/0
3/0
2/0

1 x 45 cm
1 x 45 cm
1 x 45 cm

Blue
Blue
Blue

PP 5515
PP 5516
PP 5517

1.5
2
3
3.5


4/0
3/0
2/0
0

1 x 75 cm
1 x 75 cm
1 x 75 cm
1 x 75 cm

Blue
Blue
Blue
Blue

PP 5115
PP 5116
PP 5117
PP 5118

CHIRALEN

Needle	EP	USP	Length	Colour	Cat. No	Other information
DS: 3/8 circle, reverse cutting needle						
 <p>DS 30</p>	2	3/0	1 x 75 cm	Blue	PP 5120	
	3	2/0	1 x 75 cm	Blue	PP 5121	
	3.5	0	1 x 75 cm	Blue	PP 5122	
 <p>DS 38</p>	2	3/0	1 x 75 cm	Blue	PP 5130	
	3	2/0	1 x 75 cm	Blue	PP 5131	
	3.5	0	1 x 75 cm	Blue	PP 5132	
	4	1	1 x 75 cm	Blue	PP 5133	
 <p>DS 40</p>	3	2/0	1 x 100 cm	Blue	PP 5139	
DSC: 3/8 circle, conventional cutting needle						
 <p>DSC 15</p>	0.7	6/0	1 x 45 cm	Blue	PP 8105	
 <p>DSC 26</p>	1.5	4/0	1 x 45 cm	Blue	PP 8515	
	2	3/0	1 x 45 cm	Blue	PP 8516	
 <p>DSC 90</p>	4	1	1 x 100 cm	Blue	PP 8033	

UK03
151



Non-absorbable - Chiralen

CHIRALEN

UK03
152



Non-absorbable - Chiralen

Needle	EP	USP	Length	Colour	Cat. No	Other information
DSM: 3/8 circle, reverse cutting needle - micropoint						
 DSM 12	0.5	7/0	1 x 45 cm	Blue	PP 5000	Advanced cutting needle
	0.7	6/0	1 x 45 cm	Blue	PP 5001	
	1	5/0	1 x 45 cm	Blue	PP 5002	
	1.5	4/0	1 x 45 cm	Blue	PP 5003	
 DSM 15	0.7	6/0	1 x 45 cm	Blue	PP 5005	Advanced cutting needle
	1	5/0	1 x 45 cm	Blue	PP 5006	
	1.5	4/0	1 x 45 cm	Blue	PP 5007	
 DSM 19	1	5/0	1 x 45 cm	Blue	PP 5009	Advanced cutting needle
	1.5	4/0	1 x 45 cm	Blue	PP 5010	
	2	3/0	1 x 45 cm	Blue	PP 5011	
	3	2/0	1 x 45 cm	Blue	PP 5012	
 DSM 25	1.5	4/0	1 x 75 cm	Blue	PP 5014	Advanced cutting needle
	2	3/0	1 x 75 cm	Blue	PP 5015	
	2	3/0	1 x 45 cm	Blue	PP 5023	
	3	2/0	1 x 45 cm	Blue	PP 5024	
 DSM 25	3.5	0	1 x 45 cm	Blue	PP 5025	Advanced cutting needle
	2	3/0	1 x 75 cm	Blue	PP 5027	
	3	2/0	1 x 75 cm	Blue	PP 5028	
DSMC: 3/8 circle, conventional cutting needle - micropoint						
 DSMC 16	0.7	6/0	1 x 45 cm	Blue	PP 8005	
GS: Straight, reverse cutting needle						
 GS 60	2	3/0	1 x 75 cm	Blue	PP 5140	
	3	2/0	1 x 75 cm	Blue	PP 5141	

Silon

**Silon
monofilament**






Silon monofilament

Silon

monofilament

Silon

monofilament



SILON monofilament

SILON monofilament is a non-absorbable, synthetic suture made of polyamide, which has a smooth non-porous surface.

The thread is extremely useful for skin sutures and has a high tensile strength in knots, and knots tie down easily. The thread is available in blue and black.

Silon monofilament can cause a slight initial inflammatory tissue reaction.



UK03
155

Indication: Skin Closure (e.g. intracutaneous, subcutaneous), Plastic, General, Reconstructive and Ophthalmic surgery.

Contraindication: Silon monofilament should not be used when long-term mechanical support is required.

Chemical Composition: Polyamide 6

Sizes: 10/0 to 2 USP

Structure: Monofilament

Absorption method: Non-absorbable

Colour: Blue and Black



Non-absorbable - Silon monofilament

SILON monofilament blue ligatures



UK03
156

Pre-cut threads


EP	USP	Length	Cat. No
1	5/0	3 x 45 cm	405
1.5	4/0	3 x 45 cm	406
2	3/0	3 x 45 cm	407
3	2/0	3 x 45 cm	408
1	5/0	10 x 45 cm	409
1.5	4/0	10 x 45 cm	410
2	3/0	10 x 45 cm	411
3	2/0	10 x 45 cm	412
1	5/0	5 x 75 cm	413
1.5	4/0	5 x 75 cm	414

Pre-cut threads

EP	USP	Length	Cat. No
2	3/0	5 x 75 cm	415
3	2/0	5 x 75 cm	416
2	3/0	2 x 100 cm	1485
3.5	0	2 x 100 cm	1487
1	5/0	1 x 150 cm	401
1.5	4/0	1 x 150 cm	402
2	3/0	1 x 150 cm	403
3	2/0	1 x 150 cm	404
3.5	0	1 x 150 cm	1405
4	1	1 x 150 cm	1406



Non-absorbable - Silon monofilament



SILON monofilament

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
HR 10 	1	5/0	1 x 45 cm	Blue	SM 2453	
HR 13 	1.5	4/0	1 x 45 cm	Blue	4125	
HR 18 	1.5	4/0	1 x 75 cm	Blue	4104	
HR 22 	2	3/0	1 x 75 cm	Blue	4106	
HR 27 	3	2/0	1 x 75 cm	Blue	SM 4137	
HR 30 	3	2/0	1 x 75 cm	Blue	SM 4142	
HR 35 	3.5	0	1 x 75 cm	Blue	SM 4149	
HR 40 	4	1	1 x 75 cm	Blue	SM 4152	
	3.5	0	1 x 100 cm	Blue	SM 4154	
HRX: ½ circle, round body needle - extra strong						
HRX 35 	3.5	0	1 x 75 cm	Blue	SM 4241	
HRX 40 	4	1	1 x 150 cm	Blue	SM 5033	Loop

UK03
157



Non-absorbable - Silon monofilament

SILON monofilament

UK03
158



Non-absorbable - Silon monofilament

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRX: ½ circle, round body needle - extra strong						
HRX 48 	4	1	1 x 150 cm	Blue	SM 5034	Loop
	4	1	1 x 200 cm	Blue	SM 5037	Loop
HRX 50 	3.5	0	1 x 150 cm	Blue	SM 5035	Loop
	4	1	1 x 200 cm	Blue	SM 5040	Loop
HRT: ½ circle, round body needle - tapercut						
HRT 13 	0.7 1 1.5	6/0 5/0 4/0	1 x 45 cm 1 x 45 cm 1 x 45 cm	Blue Blue Blue	SM 2552 SM 2553 SM 2554	
DR: 3/8 circle, round body needle						
DR 8 	0.5 1	7/0 5/0	1 x 45 cm 1 x 45 cm	Blue Blue	SM 2448 450	
DR 10 	1	5/0	1 x 45 cm	Blue	451	
DR 15 	1.5	4/0	1 x 45 cm	Blue	4024	
DR 18 	1.5 2	4/0 3/0	1 x 75 cm 1 x 75 cm	Blue Blue	4003 4021	
DR 21 	2 3	3/0 2/0	1 x 75 cm 1 x 75 cm	Blue Blue	SM 2325 SM 2326	
DR 31 	2 3	3/0 2/0	1 x 75 cm 1 x 75 cm	Blue Blue	SM 2330 SM 2331	

SILON monofilament

UK03
159



Non-absorbable - Silon monofilament

Needle	EP	USP	Length	Colour	Cat. No	Other information
--------	----	-----	--------	--------	---------	-------------------

GR: Straight, round body needle

2 x GR 60



3.5

0

1 x 90 cm

Blue

SM 2358



HS: 1/2 circle, reverse cutting needle

HS 12



1

5/0

1 x 45 cm

Blue

SM 2305

1

5/0

1 x 75 cm

Blue

SM 2306

HS 15



2

3/0

1 x 45 cm

Blue

SM 2308

HS 18



1.5

4/0

1 x 75 cm

Blue

4420

HS 25



2

3/0

1 x 75 cm

Blue

SM 2311

HS 40



4

1

1 x 75 cm

Blue

SM 2413

3.5

0

1 x 100 cm

Blue

SM 3382

DS: 3/8 circle, reverse cutting needle

DS 10



0.5

7/0

1 x 45 cm

Blue

SM 2251

0.7

6/0

1 x 45 cm

Blue

SM 2252

DS 11



0.5

7/0

1 x 45 cm

Blue

SM 2254

DS 12



0.5

7/0

1 x 45 cm

Blue

4328

0.7

6/0

1 x 45 cm

Blue

SM 2255

1

5/0

1 x 45 cm

Blue

SM 2256

1.5

4/0

1 x 45 cm

Blue

SM 2257

DS 13



1

5/0

1 x 45 cm

Blue

SM 2248

SILON monofilament

UK03
160



Non-absorbable - Silon monofilament

Needle	EP	USP	Length	Colour	Cat. No	Other information
DS: 3/8 circle, reverse cutting needle						
 DS 15	0.7	6/0	1 x 45 cm	Blue	SM 2260	
	1	5/0	1 x 45 cm	Blue	SM 2261	
	1.5	4/0	1 x 45 cm	Blue	SM 2262	
	2	3/0	1 x 45 cm	Blue	SM 2263	
	1.5	4/0	1 x 75 cm	Blue	SM 2362	
 DS 16	0.7	6/0	1 x 45 cm	Blue	SM 2460	
	0.7	6/0	1 x 60 cm	Blue	SM 2459	
 DS 19	0.7	6/0	1 x 45 cm	Blue	4332	
	1	5/0	1 x 45 cm	Blue	SM 2266	
	1.5	4/0	1 x 45 cm	Blue	SM 2367	
	2	3/0	1 x 45 cm	Blue	SM 2368	
	3	2/0	1 x 45 cm	Blue	SM 2369	
	1	5/0	1 x 75 cm	Blue	SM 2366	
	1.5	4/0	1 x 75 cm	Blue	SM 2267	
	2	3/0	1 x 75 cm	Blue	SM 2268	
	3	2/0	1 x 75 cm	Blue	SM 2269	
	1	5/0	1 x 90 cm	Blue	SM 2466	
 DS 20	1.5	4/0	1 x 75 cm	Blue	SM 2471	
 DS 24	1.5	4/0	1 x 75 cm	Blue	SM 3277	
	2	3/0	1 x 75 cm	Blue	SM 3276	
	3	2/0	1 x 75 cm	Blue	SM 3275	
	1.5	4/0	1 x 45 cm	Blue	SM 2372	
	2	3/0	1 x 45 cm	Blue	SM 2373	
	3	2/0	1 x 45 cm	Blue	SM 2374	
 DS 25	1	5/0	1 x 75 cm	Blue	SM 2375	
	1.5	4/0	1 x 75 cm	Blue	SM 2273	
	2	3/0	1 x 75 cm	Blue	SM 2270	
	3	2/0	1 x 75 cm	Blue	SM 2271	
	3.5	0	1 x 75 cm	Blue	SM 2272	
	4	1	1 x 75 cm	Blue	SM 2274	
	2	3/0	1 x 90 cm	Blue	SM 3272	
	3	2/0	1 x 90 cm	Blue	SM 3271	

SILON monofilament

UK03
161



Non-absorbable - Silon monofilament

Needle	EP	USP	Length	Colour	Cat. No	Other information
DS: 3/8 circle, reverse cutting needle						
	3	2/0	1 x 45 cm	Blue	SM 2384	
	3.5	0	1 x 45 cm	Blue	SM 2385	
 <p>DS 30</p>	1.5	4/0	1 x 75 cm	Blue	SM 4324	
	2	3/0	1 x 75 cm	Blue	4325	
	3	2/0	1 x 75 cm	Blue	SM 2284	
	3.5	0	1 x 75 cm	Blue	SM 2285	
	2	3/0	1 x 90 cm	Blue	SM 4325	
	3	2/0	1 x 90 cm	Blue	SM 4326	
	4	1	1 x 90 cm	Blue	SM 4327	
 <p>DS 37</p>	3.5	0	1 x 75 cm	Blue	SM 2386	
 <p>DS 38</p>	2	3/0	1 x 75 cm	Blue	SM 2290	
	3	2/0	1 x 75 cm	Blue	SM 2291	
	3.5	0	1 x 75 cm	Blue	SM 2292	
	4	1	1 x 75 cm	Blue	SM 2293	
	3	2/0	1 x 90 cm	Blue	SM 3291	
 <p>DS 40</p>	3	2/0	1 x 90 cm	Blue	SM 2497	
	3.5	0	1 x 90 cm	Blue	SM 2495	
 <p>DS 50</p>	4	1	1 x 100 cm	Blue	SM 2693	
 <p>DS 60</p>	3	2/0	1 x 100 cm	Blue	SM 3388	
	4	1	1 x 100 cm	Blue	SM 3390	
 <p>DS 90</p>	4	1	1 x 100 cm	Blue	SM 3395	

SILON monofilament

UK03
162

Needle	EP	USP	Length	Colour	Cat. No	Other information
--------	----	-----	--------	--------	---------	-------------------

DSC: 3/8 circle, conventional cutting needle

DSC 90



3
3.5

2/0
0

1 x 100 cm
1 x 100 cm

Blue
Blue

SM 8397
SM 8394

DSM: 3/8 circle, reverse cutting needle - micropoint

DSM 12



0.7
1
1.5

6/0
5/0
4/0

1 x 45 cm
1 x 45 cm
1 x 45 cm

Blue
Blue
Blue

SM 2055
SM 2056
SM 2057

Advanced
cutting
needle

1

5/0

1 x 75 cm

Blue

SM 2053

DSM 15



1
1.5
2

5/0
4/0
3/0

1 x 45 cm
1 x 45 cm
1 x 45 cm

Blue
Blue
Blue

SM 2060
SM 2061
SM 2062

Advanced
cutting
needle

DSM 19



1
1.5
2

5/0
4/0
3/0

1 x 45 cm
1 x 45 cm
1 x 45 cm

Blue
Blue
Blue

SM 2064
SM 2065
SM 2066

Advanced
cutting
needle

1

5/0

1 x 75 cm

Blue

SM 2166

2

3/0

1 x 75 cm

Blue

SM 2168

DSM 25



1.5
2

4/0
3/0

1 x 45 cm
1 x 45 cm

Blue
Blue

SM 2073
SM 2074

Advanced
cutting
needle

2

3/0

1 x 75 cm

Blue

SM 2078

3

2/0

1 x 75 cm

Blue

SM 2079

DSM 40



3

2/0

1 x 75 cm

Blue

SM 2090

Advanced
cutting
needle

DSMC: 3/8 circle, conventional cutting needle - micropoint

DSMC 19



1

5/0

1 x 100 cm

Blue

SM 8081

Advanced
cutting
needle

DSMC 20



1.5

4/0

1 x 45 cm

Blue

SM 8065

Advanced
cutting
needle

SILON monofilament

Needle	EP	USP	Length	Colour	Cat. No	Other information
GS: Straight, reverse cutting needle						
GS 50 	3	2/0	1 x 90 cm	Blue	SM 3295	
GS 60 	2 3	3/0 2/0	1 x 75 cm 1 x 75 cm	Blue Blue	SM 3283 SM 3284	
GS 75 	3.5	0	1 x 90 cm	Blue	SM 3294	
			1 x 75 cm	Blue	SM 3299	

UK03
163



Non-absorbable - Silon monofilament






Tervalon

Tervalon

Tervalon

Tervalon



TERVALON

TERVALON is a braided, polyester suture material with a special coating, which guarantees effortless tissue passage and excellent knot-tying characteristics.

The extreme tensile strength and tissue compatibility are further positive characteristics of Tervalon, which is used in all areas of surgery due to its favourable features.

Tervalon is applicable wherever non-absorbability and low-friction passage through tissues are needed.

The suture elicits a minimal inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissue.



UK03
167



Non-absorbable - Tervalon

Indication: Cardiovascular, Plastic, Reconstructive and General surgery.

Contraindication: Tervalon must not penetrate into the inner parts of the digestive and urinary tracts.

Chemical Composition: Polyester

Sizes: 6/0 to 6 USP

Structure: Multifilament

Absorption method: Non-absorbable

Colour: Green and White

TERVALON green ligatures

UK03
168



Non-absorbable - Tervalon

Pre-cut threads

EP	USP	Length	Cat. No
2	3/0	10 x 30 cm	TG 4064
1	5/0	3 x 45 cm	TG 4032
1.5	4/0	3 x 45 cm	TG 4033
2	3/0	3 x 45 cm	TG 4034
3	2/0	3 x 45 cm	TG 4035
3.5	0	3 x 45 cm	TG 4036
4	1	3 x 45 cm	TG 4037
5	2	3 x 45 cm	TG 4038
1.5	4/0	6 x 45 cm	TG 4043
2	3/0	6 x 45 cm	TG 4044
3	2/0	6 x 45 cm	TG 4045
3.5	0	6 x 45 cm	TG 4046
4	1	6 x 45 cm	TG 4047
5	2	6 x 45 cm	TG 4048
1.5	4/0	10 x 45 cm	TG 4053
2	3/0	10 x 45 cm	TG 4054
3	2/0	10 x 45 cm	TG 4055
3.5	0	10 x 45 cm	TG 4056
4	1	10 x 45 cm	TG 4057
5	2	10 x 45 cm	TG 4058
3.5	0	15 x 45 cm	TG 4142
6	3 + 4	4 x 60 cm	TG 4099
2	3/0	13 x 60 cm	TG 4086
2	3/0	2 x 70 cm	TG 4014
3	2/0	2 x 70 cm	TG 4015
3.5	0	2 x 70 cm	TG 4016
4	1	2 x 70 cm	TG 4017
5	2	2 x 70 cm	TG 4018

Pre-cut threads

EP	USP	Length	Cat. No
1.5	4/0	5 x 70 cm	TG 4080
3	2/0	5 x 70 cm	TG 4082
4	1	5 x 70 cm	TG 4084
5	2	5 x 70 cm	TG 4085
3	2/0	10 x 70 cm	TG 3045
3.5	0	10 x 70 cm	TG 3046
2	3/0	6 x 75 cm	TG 3035
3.5	0	6 x 75 cm	TG 3037
4	1	6 x 75 cm	TG 3038
3.5	0	1 x 150 cm	TG 4076
4	1	1 x 150 cm	TG 4077
5	2	1 x 150 cm	TG 4078
6	3 + 4	1 x 150 cm	TG 4079
7	5	1 x 150 cm	TG 4180

Chirapack - 20 pcs./box

EP	USP	Length	Cat. No
4	1	1 x 250 cm	TG 3097
5	2	1 x 250 cm	TG 3098
1.5	4/0	1 x 300 cm	TG 4093
2	3/0	1 x 300 cm	TG 4094
3	2/0	1 x 300 cm	TG 4095
3.5	0	1 x 300 cm	TG 4096
4	1	1 x 300 cm	TG 4097
5	2	1 x 300 cm	TG 4098
6	3 + 4	1 x 300 cm	TG 4100

TERVALON

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
HR 13 	1	5/0	1 x 45 cm	Green	TG 4452	
	1.5	4/0	1 x 45 cm	Green	TG 4453	
2 x HR 13 	1.5	4/0	1 x 75 cm	Green	TG 4495	
HR 18 	1.5	4/0	1 x 75 cm	Green	TG 4456	
	2	3/0	1 x 75 cm	Green	TG 4457	
	3	2/0	1 x 75 cm	Green	TG 4458	
HR 18 A 	2	3/0	1 x 45 cm	Green	TG 5447	Take-off needle
HR 22 	1.5	4/0	1 x 75 cm	Green	TG 4465	
	2	3/0	1 x 75 cm	Green	TG 4466	
	3	2/0	1 x 75 cm	Green	TG 4467	
2 x HR 22 	3	2/0	1 x 75 cm	Green	TG 4653	
	3	2/0	1 x 90 cm	Green	TG 4652	
HR 25 	2	3/0	1 x 75 cm	Green	TG 4471	
2 x HR 26 	3	2/0	1 x 90 cm	Green	TG 4650	
HR 26 A 	3	2/0	5 x 45 cm	Green	TG 5452	Take-off needle
HR 27 	2	3/0	1 x 75 cm	Green	TG 4461	
	3	2/0	1 x 75 cm	Green	TG 4462	
	3.5	0	1 x 75 cm	Green	TG 4463	
	4	1	1 x 75 cm	Green	TG 4464	
HR 30 	2	3/0	1 x 75 cm	Green	TG 4315	
	3	2/0	1 x 75 cm	Green	TG 4316	
	3.5	0	1 x 75 cm	Green	TG 4317	
HR 35 	4	1	1 x 75 cm	Green	TG 4373	



TERVALON

UK03
170



Non-absorbable - Tervalon

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
HR 37 	3.5 4	0 1	1 x 75 cm 1 x 75 cm	Green Green	TG 4377 TG 4374	
HR 40 	5	2	1 x 90 cm	Green	TG 4376	
HR s: ½ circle, round body needle - slim						
HR s 18 	1.5	4/0	1 x 75 cm	Green	TG 2474	Slim needle
HR s 22 	2	3/0	1 x 75 cm	Green	TG 2485	Slim needle
HRX: ½ circle, round body needle - extra strong						
HRX 26 	3.5	0	1 x 75 cm	Green	TG 4362	
HRX 27 	5	2	1 x 75 cm	Green	TG 4365	
HRX 35 	2 3.5	3/0 0	1 x 75 cm 1 x 75 cm	Green Green	TG 4318 TG 4320	
HRX 50 	6	3	1 x 90 cm	Green	TG 4673	
HRX 70 	6	3 + 4	1 x 180 cm	Green	TG 7060	Loop

TERVALON

UK03
171



Non-absorbable - Tervalon

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRT: ½ circle, round body needle - tapercut						
HRT 17 	1	5/0	1 x 45 cm	Green	TG 4622	
	1.5	4/0	1 x 45 cm	Green	TG 4623	
	2	3/0	1 x 45 cm	Green	TG 4624	
2 x HRT 17 	1.5	4/0	1 x 75 cm	Green	TG 4627	
	2	3/0	1 x 75 cm	Green	TG 4628	
	1	5/0	1 x 90 cm	Green	TG 4631	
2 x HRT 17 	1.5	4/0	1 x 90 cm	Green	TG 4632	
	2	3/0	1 x 90 cm	Green	TG 4633	
	HRT 17 A	2	3/0	6 x 45 cm	Green	TG 4734
HRT 20 	3	2/0	1 x 75 cm	Green	TG 4635	
2 x HRT 25 	3	2/0	1 x 90 cm	Green	TG 4645	
HRT 37 	5	2	1 x 75 cm	Green	TG 4758	
HRT 45 	5	2	4 x 75 cm	Green	TG 5426	
HRT 50 	7	5	4 x 75 cm	Green	TG 5438	

TERVALON

UK03
172



Non-absorbable - Tervalon

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRTX: 1/2 circle, round body needle – tapercut, extra strong						
 <p>HRTX 50</p>	7	5	1 x 75 cm	Green	TG 4438	
DR: 3/8 circle, round body needle						
 <p>DR 12</p>	1	5/0	1 x 45 cm	Green	TG 4326	
 <p>DR 18</p>	2	3/0	1 x 75 cm	Green	TG 4331	
 <p>DR 21</p>	3	2/0	1 x 75 cm	Green	TG 4336	
DRT: 3/8 circle, round body needle - tapercut						
 <p>2 x DRT 17</p>	2	3/0	1 x 90 cm	Green	TG 4649	
FR: 5/8 circle, round body needle						
 <p>FR 27</p>	2	3/0	1 x 45 cm	Green	TG 4499	
GR: Straight, round body needle						
 <p>GR 20</p>	1 1.5 2 3	5/0 4/0 3/0 2/0	1 x 45 cm 1 x 45 cm 1 x 45 cm 1 x 45 cm	Green Green Green Green	TG 4484 TG 4485 TG 4486 TG 4487	
JRT: Hook shaped, round body needle - tapercut						
 <p>JRT 30</p>	4	1	1 x 75 cm	Green	TG 5768	

TERVALON

UK03
173

Non-absorbable - Tervalon

Needle	EP	USP	Length	Colour	Cat. No	Other information
HS: ½ circle, reverse cutting needle						
 HS 12	1	5/0	1 x 45 cm	Green	TG 4502	
 HS 18	1	5/0	1 x 45 cm	Green	TG 4505	
	1.5	4/0	1 x 45 cm	Green	TG 4506	
	2	3/0	1 x 45 cm	Green	TG 4507	
 2 x HS 18	1.5	4/0	1 x 75 cm	Green	TG 4508	
	2	3/0	1 x 75 cm	Green	TG 4509	
 2 x HS 18	1.5	4/0	1 x 90 cm	Green	TG 4618	
	2	3/0	1 x 90 cm	Green	TG 4618	
 HS 21	1.5	4/0	1 x 75 cm	Green	TG 4608	
	2	3/0	1 x 75 cm	Green	TG 4609	
 HS 25	2	3/0	1 x 75 cm	Green	TG 4511	
 2 x HS 25	3	2/0	1 x 90 cm	Green	TG 4612	
 HS 36	4	1	1 x 75 cm	Green	TG 4542	
 HS 37	3.5	0	1 x 75 cm	Green	TG 4534	
 HS 50	5	2	1 x 75 cm	Green	TG 4556	
 HS 55	6	3 + 4	1 x 75 cm	Green	TG 4567	

TERVALON

UK03
174



Non-absorbable - Tervalon

Needle	EP	USP	Length	Colour	Cat. No	Other information
DS: 3/8 circle, reverse cutting needle						
 DS 15	1	5/0	1 x 45 cm	Green	TG 4547	
	1.5	4/0	1 x 45 cm	Green	TG 4548	
	1.5	4/0	1 x 70 cm	Green	TG 6544	
	1.5	4/0	1 x 45 cm	Green	TG 6550	
 DS 19	2	3/0	1 x 45 cm	Green	TG 4551	
	3	2/0	1 x 70 cm	Green	TG 9555	
	1.5	4/0	1 x 75 cm	Green	TG 4553	
	2	3/0	1 x 75 cm	Green	TG 4554	
 DS 25	2	3/0	1 x 75 cm	Green	TG 4564	
	3	2/0	1 x 75 cm	Green	TG 4565	
	3.5	0	1 x 75 cm	Green	TG 5564	
 DS 30	2	3/0	1 x 75 cm	Green	TG 4574	
	3	2/0	1 x 75 cm	Green	TG 4575	
	3.5	0	1 x 75 cm	Green	TG 4576	
	4	1	1 x 75 cm	Green	TG 4577	
	4	1	1 x 120 cm	Green	TG 4578	
 DS 38	3.5	0	1 x 75 cm	Green	TG 4586	
	4	1	1 x 75 cm	Green	TG 4587	
	5	2	1 x 75 cm	Green	TG 4588	
 2 x DS 38	4	1	2 x 75 cm	Green	TG 5587	
DSC: 3/8 circle, conventional cutting needle						
 DSC 90	4	1	1 x 100 cm	Green	TG 8589	

TERVALON

Needle	EP	USP	Length	Colour	Cat. No	Other information
DSM: 3/8 circle, reverse cutting needle - micropoint						
DSM 12 	1 1.5	5/0 4/0	1 x 45 cm 1 x 45 cm	Green Green	TG 4513 TG 4514	Advanced cutting needle
DSM 15 	1.5	4/0	1 x 45 cm	Green	TG 4518	Advanced cutting needle
DSM 19 	2	3/0	1 x 75 cm	Green	TG 4524	Advanced cutting needle
DSM 25 	3	2/0	1 x 75 cm	Green	TG 4531	Advanced cutting needle
GS: Straight, reverse cutting needle						
GS 60 	2 3	3/0 2/0	1 x 75 cm 1 x 75 cm	Green Green	TG 4594 TG 4595	

UK03
175



Non-absorbable - Tervalon



TERVALON white ligatures

UK03
176

Pre-cut threads

EP	USP	Length	Cat. No
3	2/0	10 x 45 cm	TG 4055W
4	1	10 x 45 cm	TG 4057W
5	2	10 x 45 cm	TG 4058W

Chirapack - 20 pcs./box

EP	USP	Length	Cat. No
2	3/0	1 x 250 cm	TG 3094W
3	2/0	1 x 250 cm	TG 3095W
3.5	0	1 x 250 cm	TG 3096W
4	1	1 x 250 cm	TG 3097W



Non-absorbable - Tervalon

TERVALON

Needle	EP	USP	Length	Colour	Cat. No	Other information
DSM: 3/8 circle, reverse cutting needle - micropoint						
DSM 25 	3	2/0	1 x 75 cm	White	TG 4531W	Advanced cutting needle

Silon braided
Orsilon braided

Silon braided
Orsilon braided

Silon braided

Orsilon



SILON braided/ORSILON braided

SILON braided/ORSILON braided are non-absorbable, multifilament sutures made of polyamide. Braided sutures made of polyamide have excellent tissue biocompatibility, minimal resistance to tissue passage and high knot security. The sutures are characterized by a strong knot-pull strength.

Characteristics of Orsilon braided are higher tensile strength, lower extension, reduced capillary effect and a smoother surface.

The sutures elicit a minimal initial inflammatory reaction in tissues, which is followed by gradual encapsulation of the sutures by fibrous connective tissue.



UK03
179

Indication: Skin closure, Subcutaneous tissue, Transversely cut muscles, Fascia, Peritoneum, Gastrointestinal tract and Blood vessels ligation.

Contraindication: Silon braided/Orsilon braided should not be used for suturing infected wounds. It should not be used when long-term mechanical support is required.

Chemical Composition: Polyamide

Sizes: 5/0 to 5 USP

Structure: Multifilament

Absorption method: Non-absorbable

Colour: White and Black



SILON braided white ligatures

UK03
180



Non-absorbable - Silon braided / Orsilon braided

Pre-cut threads

EP	USP	Length	Cat. No
1	5/0	3 x 45 cm	526
1.5	4/0	3 x 45 cm	527
2	3/0	3 x 45 cm	528
2.5	-	3 x 45 cm	529
3	2/0	3 x 45 cm	530
3.5	0	3 x 45 cm	531
4	1	3 x 45 cm	532
5	2	3 x 45 cm	557
6	3 + 4	3 x 45 cm	558
1.5	4/0	6 x 45 cm	SB 2032
3.5	0	6 x 45 cm	SB 2036
4	1	6 x 45 cm	SB 2037
1	5/0	10 x 45 cm	533
1.5	4/0	10 x 45 cm	534
2	3/0	10 x 45 cm	535
2.5	-	10 x 45 cm	536
3	2/0	10 x 45 cm	537
3.5	0	10 x 45 cm	538
4	1	10 x 45 cm	539
5	2	10 x 45 cm	547
6	3 + 4	10 x 45 cm	548
1	5/0	10 x 70 cm	566
1.5	4/0	10 x 70 cm	567
2	3/0	10 x 70 cm	568
2.5	-	10 x 70 cm	569
3	2/0	10 x 70 cm	570
3.5	0	10 x 70 cm	571
4	1	10 x 70 cm	572
5	2	10 x 70 cm	573
6	3 + 4	10 x 70 cm	574
1.5	4/0	2 x 75 cm	581
2	3/0	2 x 75 cm	582
2.5	-	2 x 75 cm	583
3	2/0	2 x 75 cm	584

Pre-cut threads

EP	USP	Length	Cat. No
3.5	0	2 x 75 cm	585
4	1	2 x 75 cm	586
5	2	2 x 75 cm	587
6	3 + 4	2 x 75 cm	588
2	3/0	5 x 75 cm	592
2.5	-	5 x 75 cm	593
3	2/0	5 x 75 cm	594
3.5	0	5 x 75 cm	595
4	1	5 x 75 cm	596
5	2	5 x 75 cm	597
6	3 + 4	5 x 75 cm	598
1	5/0	1 x 150 cm	516
1.5	4/0	1 x 150 cm	517
2	3/0	1 x 150 cm	518
2.5	-	1 x 150 cm	519
3	2/0	1 x 150 cm	520
3.5	0	1 x 150 cm	521
4	1	1 x 150 cm	522
5	2	1 x 150 cm	523
6	3 + 4	1 x 150 cm	524
7	5	1 x 150 cm	525

Chirapack - 20 pcs./box

EP	USP	Length	Cat. No
1	5/0	1 x 300 cm	SB 2053
1.5	4/0	1 x 300 cm	SB 2054
2	3/0	1 x 300 cm	SB 2055
2.5	-	1 x 300 cm	SB 2056
3	2/0	1 x 300 cm	SB 2057
3.5	0	1 x 300 cm	SB 2058
4	1	1 x 300 cm	SB 2059
5	2	1 x 300 cm	SB 2060
6	3 + 4	1 x 300 cm	SB 2061
7	5	1 x 300 cm	SB 2062

SILON braided

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
HR 37 	4	1	1 x 75 cm	White	SB 2618	
HRT: ½ circle, round body needle - tapercut						
HRT 37 	3.5	0	1 x 90 cm	White	SB 2570	
J-RX: Hook shaped, round body needle						
J-RX 22 	3.5	0	1 x 75 cm	White	SB 2580	
HS: ½ circle, reverse cutting needle						
HS 37 	4	1	1 x 75 cm	White	SB 2555	
HS 55 	5	2	1 x 75 cm	White	SB 2520	
DS: 3/8 circle, reverse cutting needle						
DS 30 	3.5	0	1 x 75 cm	White	SB 4576	
DS 38 	4	1	2 x 75 cm	White	SB 5587	



SILON braided black ligatures

UK03
182



Non-absorbable - Silon braided / Orsilon braided

Pre-cut threads

EP	USP	Length	Cat. No
2	3/0	3 x 45 cm	SC 0528
3	2/0	3 x 45 cm	SC 0530
3.5	0	3 x 45 cm	SC 0531
4	1	3 x 45 cm	SC 0532
5	2	3 x 45 cm	SC 0557
6	3 + 4	3 x 45 cm	SC 0558
2	3/0	6 x 45 cm	SC 2033
3.5	0	6 x 45 cm	SC 2036
4	1	6 x 45 cm	SC 2037
1.5	4/0	10 x 45 cm	SC 0534
2	3/0	10 x 45 cm	SC 0535
3	2/0	10 x 45 cm	SC 0537
3.5	0	10 x 45 cm	SC 0538
4	1	10 x 45 cm	SC 0539
5	2	10 x 45 cm	SC 0547
2	3/0	17 x 45 cm	575
3	2/0	17 x 45 cm	576
4	1	17 x 45 cm	577
2	3/0	10 x 70 cm	SC 0568
3	2/0	10 x 70 cm	SC 0570
3.5	0	10 x 70 cm	SC 0571
4	1	10 x 70 cm	SC 0572
5	2	10 x 70 cm	SC 0573
1	5/0	2 x 75 cm	SC 0580
2	3/0	2 x 75 cm	SC 0582
3	2/0	2 x 75 cm	SC 0584
2	3/0	5 x 75 cm	SC 0592
3	2/0	5 x 75 cm	SC 0594
3.5	0	5 x 75 cm	SC 0595
4	1	5 x 75 cm	SC 0596
5	2	5 x 75 cm	SC 0597
6	3 + 4	1 x 150 cm	SC 0524

Chirapack - 20 pcs./box

EP	USP	Length	Cat. No
2	3/0	1 x 300 cm	SC 2055
3	2/0	1 x 300 cm	SC 2057
3.5	0	1 x 300 cm	SC 2058
4	1	1 x 300 cm	SC 2059
5	2	1 x 300 cm	SC 2060
6	3 + 4	1 x 300 cm	SC 2061

SILON braided

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
HR 25 	2	3/0	1 x 75 cm	Black	SC 2420	
HR 27 	2	3/0	1 x 75 cm	Black	SC 2427	
	3	2/0	1 x 75 cm	Black	SC 2426	
HR 27 A 	2	3/0	5 x 75 cm	Black	SC 5451	Take-off needle
	2	3/0	5 x 75 cm	Black	SC 2428	
HR 37 	4	1	1 x 75 cm	Black	SC 2618	
HRX: ½ circle, round body needle - extra strong						
HRX 35 A 	4	1	1 x 45 cm	Black	SC 2756	Take-off needle
HS: ½ circle, reverse cutting needle						
HS 37 	3	2/0	1 x 75 cm	Black	SC 2534	
	4	1	1 x 75 cm	Black	SC 2555	



ORSILON braided

UK03
184



Non-absorbable - Silon braided / Orsilon braided

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
HR 10 	1	5/0	1 x 45 cm	White	653	
HR 13 	1.5	4/0	1 x 45 cm	White	6114	
HR 15 	1	5/0	1 x 45 cm	White	OS 2407	
HR 18 	1.5 2	4/0 3/0	1 x 75 cm 1 x 75 cm	White White	6102 OS 2410	
HR 22 	2 2.5 3	3/0 - 2/0	1 x 75 cm 1 x 75 cm 1 x 75 cm	White White White	OS 2415 6103 OS 2416	
	2	3/0	1 x 120 cm	White	OS 2417	
HR 27 	3	2/0	1 x 75 cm	White	OS 2426	
HR 30 	2.5 3	- 2/0	1 x 75 cm 1 x 75 cm	White White	6106 6107	
DR: 3/8 circle, round body needle						
DR 15 	1.5	4/0	1 x 45 cm	White	6013	
DR 18 	1.5	4/0	1 x 75 cm	White	6001	
DR 21 	2.5	-	1 x 75 cm	White	6011	
DR 31 	2.5	-	1 x 75 cm	White	6012	

ORSILON braided

Needle	EP	USP	Length	Colour	Cat. No	Other information
HS: ½ circle, reverse cutting needle						
HS 18 	1.5	4/0	1 x 45 cm	White	6340	
HS 21 	1 1.5 2	5/0 4/0 3/0	1 x 75 cm 1 x 75 cm 1 x 75 cm	White White White	6344 6345 6346	
HS 25 	2	3/0	1 x 75 cm	White	6356	
HS 30 	3	2/0	1 x 75 cm	White	6367	
DS: 3/8 circle, reverse cutting needle						
DS 15 	1.5	4/0	1 x 45 cm	White	6317	
DS 19 	1.5	4/0	1 x 75 cm	White	6316	
DS 25 	2.5 3	- 2/0	1 x 75 cm 1 x 75 cm	White White	6318 6320	
DS 30 	2.5	-	1 x 75 cm	White	6319	






Non-absorbable - Silon braided / Orsilon braided

Silk

Silk Silk



SILK

SILK is a natural and non-absorbable, multifilament suture material. To give Silk a certain stiffness a special technology is applied.

The braided suture is stiffened by inserting a twisted stem from the same material. The Silk is coated to achieve water resistance and avoid its high capillary effect.

Silk is characterised by its constant thread diameter and high tensile strength. The advantages of Silk sutures are their extreme pliability, ease of handling and excellent knot strength.



Indication: Skin Closure, Subcutaneous tissue, Fascia, Peritoneum and Gastrointestinal Tract.

Contraindication: Silk should not be used for suturing infected wounds. Silk cannot be used when permanent suture stability is required.

Chemical Composition: Natural silk

Sizes: 9/0 to 5 USP

Structure: Multifilament

Absorption method: Non-absorbable

Colour: Black and Blue

SILK black ligatures

UK03
190

Pre-cut threads

EP	USP	Length	Cat. No
3	2/0	5 x 30 cm	HB 7055
1.5	4/0	3 x 45 cm	HB 4033
2	3/0	3 x 45 cm	HB 4034
3	2/0	3 x 45 cm	HB 4035
5	2	3 x 45 cm	HB 4038
1.5	4/0	12 x 45 cm	HB 3043
2	3/0	12 x 45 cm	HB 3044
3	2/0	12 x 45 cm	HB 3045
3.5	0	12 x 45 cm	HB 3046
1	5/0	1 x 150 cm	317
1.5	4/0	1 x 150 cm	318
2	3/0	1 x 150 cm	319
3	2/0	1 x 150 cm	321
3.5	0	1 x 150 cm	322
4	1	1 x 150 cm	323
5	2	1 x 150 cm	324
6	3 + 4	1 x 150 cm	325
7	5	1 x 150 cm	326

Chirapack - 20 pcs./box

EP	USP	Length	Cat. No
2	3/0	1 x 300 cm	HB 3053
3	2/0	1 x 300 cm	HB 3054
3.5	0	1 x 300 cm	HB 3055
4	1	1 x 300 cm	HB 3056



Non-absorbable - Silk

SILK

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
 HR 10	1	5/0	1 x 45 cm	Black	355	
 HR 13	1	5/0	1 x 45 cm	Black	3102	
 HR 15	0.7 2	6/0 3/0	1 x 75 cm 1 x 75 cm	Black Black	3110 3113	
 HR 16	1.5	4/0	1 x 75 cm	Black	3124	
 HR 18	1.5	4/0	1 x 75 cm	Black	3104	
 HR 22	2	3/0	1 x 75 cm	Black	3108	
 HR 27	3 3.5	2/0 0	1 x 75 cm 1 x 75 cm	Black Black	3115 3116	
 HR 30	2 3	3/0 2/0	1 x 75 cm 1 x 75 cm	Black Black	3118 3119	
 HR 35	4	1	1 x 75 cm	Black	3131	

UK03
191



Non-absorbable - Silk

SILK

UK03
192

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRT: ½ circle, round body needle - tapercut						
 HRT 17	1.5	4/0	1 x 45 cm	Black	3153	
	2	3/0	1 x 45 cm	Black	3154	
DR: 3/8 circle, round body needle						
 DR 8	0.5	7/0	1 x 45 cm	Black	350	
	1	5/0	1 x 45 cm	Black	352	
 DR 10	0.5	7/0	1 x 45 cm	Black	351	
 DR 12	1	5/0	1 x 45 cm	Black	3001	
 DR 15	1.5	4/0	1 x 45 cm	Black	3012	
	1	5/0	1 x 75 cm	Black	3010	
 DR 18	1.5	4/0	1 x 75 cm	Black	3003	
 DR 21	1.5	4/0	1 x 75 cm	Black	3011	
GR: Straight, round body needle						
 GR 20	1.5	4/0	1 x 75 cm	Black	3205	
	2	3/0	1 x 75 cm	Black	3206	
 GR 45	3	2/0	1 x 75 cm	Black	3217	



Non-absorbable - Silk

SILK

UK03
193

Non-absorbable - Silk

Needle	EP	USP	Length	Colour	Cat. No	Other information
HS: ½ circle, reverse cutting needle						
 HS 18	1.5	4/0	1 x 45 cm	Black	3402	
	1.5	4/0	1 x 75 cm	Black	3405	
	2	3/0	1 x 75 cm	Black	3406	
	3	2/0	1 x 75 cm	Black	3407	
 HS 21	3	2/0	1 x 75 cm	Black	3427	
 HS 24	2	3/0	1 x 75 cm	Black	3434	
 HS 40	4	1	1 x 75 cm	Black	3465	
DS: 3/8 circle, reverse cutting needle						
 DS 12	0.7	6/0	1 x 45 cm	Black	3312	
	1	5/0	1 x 45 cm	Black	3313	
 DS 15	1.5	4/0	1 x 45 cm	Black	3320	
 DS 16	1.5	4/0	1 x 75 cm	Black	3331	
 DS 19	1.5	4/0	1 x 45 cm	Black	3335	
	2	3/0	1 x 75 cm	Black	3336	
	3	2/0	1 x 75 cm	Black	3337	

SILK

UK03
194



Non-absorbable - Silk

Needle	EP	USP	Length	Colour	Cat. No	Other information
--------	----	-----	--------	--------	---------	-------------------

DS: 3/8 circle, reverse cutting needle

	1.5	4/0	1 x 45 cm	Black	3343	
	2	3/0	1 x 45 cm	Black	3344	
	3	2/0	1 x 45 cm	Black	3345	

	3	2/0	1 x 75 cm	Black	3348	
	4	1	1 x 75 cm	Black	3347	

	3	2/0	1 x 75 cm	Black	3351	
	3.5	0	1 x 75 cm	Black	3352	
	4	1	1 x 75 cm	Black	3353	

	3.5	0	1 x 75 cm	Black	3366	
---	-----	---	-----------	-------	------	--

DSC: 3/8 circle, conventional cutting needle

	0.7	6/0	1 x 45 cm	Black	8326	
	1	5/0	1 x 75 cm	Black	8327	
	3	2/0	1 x 75 cm	Black	8325	

	3	2/0	1 x 100 cm	Black	8365	
---	---	-----	------------	-------	------	--

	4	1	1 x 75 cm	Black	8377	
--	---	---	-----------	-------	------	--

DSM: 3/8 circle, reverse cutting needle - micropoint

	2	3/0	1 x 45 cm	Black	3534	Advanced cutting needle
---	---	-----	-----------	-------	------	-------------------------

GS: Straight, reverse cutting needle

	3	2/0	1 x 75 cm	Black	3225	
--	---	-----	-----------	-------	------	--

	2	3/0	1 x 75 cm	Black	3223	
	3.5	0	1 x 75 cm	Black	3224	



Steel Wire



STEEL Wire

STEEL Wire is a non-absorbable, monofilament and sterile surgical suture composed of 316LVM stainless steel for surgical implants. Steel Wire material is supplied as combinations of stainless steel needles.

Steel Wire sutures are the most inert of all suture materials and provide maximum tensile strength. Steel Wire is available with fixed or rotating (axial) needles. All boxes are supplied in long packs.



UK03
197



Non-absorbable - Steel Wire

Indication: Abdominal wound closure, Sternal closure, Orthopaedic procedures and Tendons.

Contraindication: Steel Wire is contraindicated in patients with known sensitivities or allergies to steel and/or its principal metallic components, such as chromium and nickel. The presence of steel may interfere with certain radio-diagnostics and its use is contraindicated where radiotransparency of suture material is required.

Chemical Composition: 316LVM Stainless steel

Sizes: 1 to 7 USP

Structure: Monofilament

Absorption method: Non-absorbable

Colour: Steel

STEEL Wire

UK03
198



Non-absorbable - Steel Wire

Needle	EP	USP	Length	Colour	Cat. No	Other information
--------	----	-----	--------	--------	---------	-------------------

HRX: 1/2 circle, round body needle - extra strong

HRX 48 rotating



8

6

4 x 45 cm

Steel

SW 2122

Long pack

HRTX: 1/2 circle, round body needle - tapercut, extra strong

HRTX 50



7

5

4 x 45 cm

Steel

SW 2200

Long pack

8

6

4 x 45 cm

Steel

SW 2210

HS: 1/2 circle, reverse cutting needle

HS 48



4

1

4 x 45 cm

Steel

SW 2000

Long pack

6

4

4 x 45 cm

Steel

SW 2010

7

5

4 x 45 cm

Steel

SW 2100

8

6

4 x 45 cm

Steel

SW 2110

9

7

4 x 45 cm

Steel

SW 2120

HS 55 rotating



7

5

4 x 45 cm

Steel

SW 2310

Long pack

8

6

4 x 45 cm

Steel

SW 2330

7

5

2 x 75 cm

Steel

SW 2300

8

6

2 x 75 cm

Steel

SW 2320

HSC: 1/2 circle, conventional cutting needle

HSC 48



7

5

4 x 45 cm

Steel

SW 2400

Long pack

8

6

4 x 45 cm

Steel

SW 2410

9

7


4 x 45 cm

Steel

SW 2420

The background is a solid dark blue color. It features several large, overlapping, semi-transparent shapes in a lighter shade of blue. These shapes are organic and flowing, resembling stylized leaves or petals. Some are circular or semi-circular, while others are more elongated and pointed. The overall effect is a layered, textured look.

Special products



Fix-clip

Fix-clip Cat. No. FS 0030

Intracutaneous suture with Fix-clip for an optimal cosmetic result.

Pierce the skin at two points only


No suture knots

Fixation with two clips

Easy to remove after wound healing: simply cut off one of the clips

Easy suture removal by pulling out the second clip

For threads EP 1-2 (USP 5/0-3/0)



UK03
201

Special Products

Minimal invasive surgery

Material	EP	USP	Length	Cat. No
CHIRASORB violet, Loop + Applicator	3	2/0	1 x 50 cm	MS 0143
CHIRASORB violet, Loop + Applicator	3.5	0	1 x 50 cm	MS 0144
CHIRASORB violet, Loop + Applicator, KR 27	4	1	1 x 120 cm	MS 1125
CHIRLAC violet, Loop + Applicator	4	1	1 x 70 cm	MS 0125
CHIRLAC violet, Loop + Applicator	5	2	1 x 70 cm	MS 0126
POLYDOX violet, Loop + Applicator	3.5	0	1 x 70 cm	MS 0134
POLYDOX violet, Loop + Applicator	4	1	1 x 70 cm	MS 0135
SILON monofilament blue, Loop + Applicator	3.5	0	1 x 70 cm	MS 1124

Ophthalmic - Microsurgery

UK03
202



Special Products


Type of needles

- DR - 3/8 circle, round body needle
- DSM - 3/8 circle, reverse cutting needle, micropoint
- VSM - 1/4 circle, reverse cutting needle, micropoint
- HLZ - 1/2 circle, spatulated
- DLZ - 3/8 circle, spatulated
- VLZ - 1/4 circle, spatulated
- GLZ - Straight, spatulated



Material	Needle	EP	USP	Length	Cat. No
CHIRASORB violet	2 x DLZ 6.4 / 200 µm	0.4	8/0	1 x 45 cm	MI 2422
		0.5	7/0	1 x 30 cm	MI 2421
CHIRASORB violet	2 x VLZ 8 / 350 µm	0.5	7/0	1 x 45 cm	MI 2609
		0.7	6/0	1 x 45 cm	MI 2610
CHIRLAC violet	1 x DLZ 6.4 / 200 µm	0.4	8/0	1 x 30 cm	MI 1419
		0.5	7/0	1 x 30 cm	MI 1420
CHIRLAC violet	2 x DLZ 6.4 / 200 µm	0.4	8/0	1 x 30 cm	MI 1426
		0.5	7/0	1 x 30 cm	MI 1427
CHIRLAC violet	2 x VLZ 8 / 350 µm	0.7	6/0	1 x 45 cm	MI 1610
		1	5/0	1 x 45 cm	MI 1710

Ophthalmic - Microsurgery

Material	Needle	EP	USP	Length	Cat. No
CHIRALEN black	1 x DR 6 	0.4	8/0	1 x 45 cm	MI 3241
CHIRALEN black	2 x DR 6.5 	0.4	8/0	1 x 60 cm	MI 3242
CHIRALEN black	1 x DLZ 6.4 / 200 µm 	0.2	10/0	1 x 15 cm	MI 1219
CHIRALEN black	2 x DLZ 6.4 / 200 µm 	0.2	10/0	1 x 23 cm	MI 1220
CHIRALEN black	2 x GLZ 16 / 150 µm 	0.2	10/0	1 x 23 cm	MI 1280
SILON black monofilament	1 x DR 3.75 / 75 µm 	0.2	10/0	1 x 13 cm	MI 2248
SILON black monofilament	1 x DR 5 Black / 130 µm 	0.4	8/0	1 x 13 cm	MI 2241
SILON black monofilament	1 x DR 5 / 100 µm 	0.3	9/0	1 x 15 cm	MI 2344
SILON black monofilament	1 x DR 5.1 / 100 µm 	0.3	9/0	1 x 13 cm	MI 2345
SILON black monofilament	1 x DR 5.1 / 75 µm 	0.2	10/0	1 x 13 cm	MI 2246

UK03
203



Special Products

Ophthalmic - Microsurgery

UK03
204



Special Products

Type of needles

- DR - 3/8 circle, round body needle
- DSM - 3/8 circle, reverse cutting needle, micropoint
- VSM - 1/4 circle, reverse cutting needle, micropoint
- HLZ - 1/2 circle, spatulated
- DLZ - 3/8 circle, spatulated
- VLZ - 1/4 circle, spatulated
- GLZ - Straight, spatulated

Material	Needle	EP	USP	Length	Cat. No
SILON black monofilament	1 x DR 6.5 / 150 μ m 	0.4	8/0	1 x 15 cm	MI 1240
SILON black monofilament	1 x DSM 6.5 / 200 μ m 	0.4	8/0	1 x 13 cm	MI 2448
SILON black monofilament	1 x DLZ 6.2 / 150 μ m 	0.2	10/0	1 x 30 cm	MI 2223
SILON black monofilament	2 x DLZ 6.2 / 150 μ m 	0.2 0.3	10/0 9/0	1 x 30 cm 1 x 30 cm	MI 1225 MI 1226
TERVALON green	2 x VLZ 8 / 350 μ m 	1	5/0	1 x 45 cm	MI 8710
SILK black twisted	2 x DSM 8 / 230 μ m 	0.5	7/0	1 x 45 cm	MI 1550
SILK black twisted	2 x HLZ 8 / 200 μ m 	0.4	8/0	1 x 30 cm	MI 1430

Ophthalmic - Microsurgery

Material	Needle	EP	USP	Length	Cat. No
SILK black	2 x VSM 8 / 200 μ m 	0.7	6/0	1 x 45 cm	MI 1690
SILK black	2 x DLZ 6.4 / 200 μ m 	0.4	8/0	1 x 30 cm	MI 1431
		0.5	7/0	1 x 45 cm	MI 1432

UK03
205



Special Products



Eye-needles

Surgical Eye-needles are made of high-quality stainless steel. Different types of needles are available to meet a full range of requirements.

The needles are always packed in foil pouches with 12 pieces.

The surgical Eye-needle range consists of Eye-needles with a double spring eye and with a regular eye.

Eye-needles
Double spring



Eye-needles
Regular



Eye-needles

Double spring



Needle	Needle no.	Wire diameter (mm)	Cat. No
PB HR: ½ circle, round body needle			
<p>HR 61</p>	0	0.65	EY 8002
<p>HR 52</p>	1	0.65	EY 8003
<p>HR 42</p>	2	0.65	EY 8004
<p>HR 36</p>	3	0.55	EY 8005
<p>HR 30</p>	4	0.55	EY 8006
<p>HR 24</p>	5	0.55	EY 8007
<p>HR 20</p>	7	0.55	EY 8009
<p>HR 13</p>	9	0.45	EY 8011

UK03
207



Special Products - Eye-needles

Eye-needles

Double spring



UK03
208

Needle	Needle no.	Wire diameter (mm)	Cat. No
E HRX: 1/2 circle, round body needle - extra strong			



Special Products - Eye-needles

HRX 100



0000

1.20

EY 8700

HRX 72




000

1.10

EY 8701

HRX 63



00

1.10

EY 8703

HRX 54



0

1.00

EY 8705



Eye-needles

Double spring



Needle	Needle no.	Wire diameter (mm)	Cat. No
E HRX: ½ circle, round body needle - extra strong			
HRX 46 	1	1.00	EY 8707
HRX 38 	2	0.85	EY 8709
HRX 30 	3	0.75	EY 8711
HRX 24 	4	0.65	EY 8713
HRX 22 	5	0.55	EY 8715
HRX 17.5 	6	0.55	EY 8717
PC GR: Straight, round body needle			
GR 70 	00	0.65	EY 8801
GR 52 	1	0.65	EY 8803
GR 36 	3	0.55	EY 8804
GR 24 	5	0.55	EY 8805

UK03
209



Special Products - Eye-needles

Eye-needles

Double spring



UK03
210

Needle	Needle no.	Wire diameter (mm)	Cat. No
GA HS: 1/2 circle, reverse cutting needle			









Special Products - Eye-needles

<p>HS 72</p>	0	1.80	EY 8623
<p>HS 60</p>	2	1.60	EY 8625
<p>HS 54</p>	3	1.40	EY 8626
<p>HS 46</p>	4	1.40	EY 8627
<p>HS 38</p>	5	1.30	EY 8628
<p>HS 32</p>	6	1.10	EY 8629
<p>HS 23</p>	8	1.00	EY 8631

Eye-needles

Double spring



Needle	Needle no.	Wire diameter (mm)	Cat. No
GE HS: 1/2 circle, reverse cutting needle			
 <p>HS 80</p>	1	1.30	EY 8601
 <p>HS 68</p>	3	1.10	EY 8603
 <p>HS 58</p>	5	1.10	EY 8605
 <p>HS 50</p>	8	1.00	EY 8608
 <p>HS 42</p>	11	0.85	EY 8611
 <p>HS 34</p>	14	0.75	EY 8614
 <p>HS 28</p>	17	0.65	EY 8617
 <p>HS 24.5</p>	19	0.65	EY 8619

UK03
211



Special Products - Eye-needles

Eye-needles

Double spring



UK03
212



Special Products - Eye-needles

Needle	Needle no.	Wire diameter (mm)	Cat. No
GE HS: 1/2 circle, reverse cutting needle			
HS 20 	21	0.55	EY 8621
HS 17.5 	22	0.55	EY 8622
BE DS: 3/8 circle, reverse cutting needle			
DS 58 	5	1.10	EY 8405
DS 50 	8	1.00	EY 8408
DS 44 	10	0.85	EY 8410
DS 36 	13	0.75	EY 8413
DS 32 	15	0.75	EY 8415
DS 28 	17	0.65	EY 8417
DS 24.5 	19	0.65	EY 8419
DS 20 	21	0.55	EY 8421

Eye-needles

Double spring



Needle	Needle no.	Wire diameter (mm)	Cat. No
BE DS: 3/8 circle, reverse cutting needle			
DS 17.5 	22	0.55	EY 8422
DS 12 	24	0.55	EY 8424

UK03
213




Special Products - Eye-needles

Eye-needles

Regular



Needle	Needle no.	Wire diameter (mm)	Cat. No
PB HR: 1/2 circle, round body needle			
HR 24 	5	0.55	EY 9007
E HRX: 1/2 circle, round body needle - extra strong			
HRX 72 	000	1.10	EY 9701
HRX 38 	2	0.85	EY 9709
GA HS: 1/2 circle, reverse cutting needle			
HS 32 	6	1.10	EY 9629

Eye-needles

Regular



UK03
214



Special Products - Eye-needles

Needle	Needle no.	Wire diameter (mm)	Cat. No
GE HS: 1/2 circle, reverse cutting needle			
<p>HS 42</p>	11	0.85	EY 9611
<p>HS 34</p>	14	0.75	EY 9614
<p>HS 24.5</p>	19	0.65	EY 9619
<p>HS 20</p>	21	0.55	EY 9621
BE DS: 3/8 circle, reverse cutting needle			
<p>DS 58</p>	5	1.10	EY 9405
<p>DS 44</p>	10	0.85	EY 9410
<p>DS 32</p>	15	0.75	EY 9415
<p>DS 28</p>	17	0.65	EY 9417
<p>DS 24.5</p>	19	0.65	EY 9419
<p>DS 20</p>	21	0.55	EY 9421

Chiralen

Meshes

Capromes

apromesh

Chiralen

Caprom

SO

Surgical Meshes

Surgical Meshes have become standard for the repair of abdominal wall hernias. Surgical mesh is a sterile woven material designed for permanent implantation within the body during open or laparoscopic procedures.

A wide range of mesh implants are available, with two main functions: to stabilise and strengthen soft tissue defects, and to act as a sling to support prolapsed organs and viscera. In tissue repair applications, the mesh serves to mechanically strengthen the weakened area whilst simultaneously promoting long-term stability by acting as a scaffold for new tissue growth. It is most commonly used in pelvic organ prolapse, urinary stress incontinence and hernia repairs.

CAPROMESH meshes


Partially absorbable

The sterile, partially absorbable, composite mesh is made from two components: one absorbable and one non-absorbable.

Capromesh is composed of an absorbable segmented copolymer of glycolide and ϵ -caprolactone and non-absorbable polypropylene monofilament fibres. The monofilament structure of the mesh decreases the risk of infection and inflammation.

After absorption of the absorbable part, only the polypropylene mesh remains in the body. The structure and size of the residual mesh are optimally designed for the physiological stress to which the abdominal wall is subject.

Indication: Capromesh is indicated in the operative treatment of weakened abdominal wall and fascial defects - in particular hernias, the healing of which requires long-term bridging or strengthening of the damaged structures..



Technical information

Fibre thickness: Approx. 500 μ m


Pore size: 2 - 4 mm

Weight: Initial: approx. 85 g/m². PP - contents approx. 28 g/m²

Stiffness: Initial: approx. 20 mg. After absorption: approx. 1 mg

Burst strength of the mesh: Initial: approx. 300 N. After absorption: approx. 175 N

Absorption time: PGA/PCL part: approx. 90 - 120 days



Advantages


Advantages for the surgeon:


- Low inflammation risk due to monofilament structure
- Mesh can be cut into any shape and remains stable
- Optimal and safe handling during surgery
- Suitable for all operation techniques
- Made from proven materials: Monolac (PGA/PCL) and Chiralen (PP)
- Good balance between price and quality

Advantages for the patient:

- Lower risk of infection
- Minimal foreign body component after partial absorption
- Optimal elongation in both directions
- Maximal comfort for the patient due to optimal softness
- Residual weight of the PP-net is only 25 g/m²

Size	Cat. No
6 cm x 11 cm	PC 0611
10 cm x 15 cm	PC 1015
15 cm x 15 cm	PC 1515
15 cm x 30 cm	PC 1530
30 cm x 30 cm	PC 3030





CHIRALEN meshes

CHIRALEN meshes

Non-absorbable


The sterile and surgical meshes are made from undyed polypropylene filaments. The knitted filaments have a smooth and hydrophobic surface, which biologically are absolutely inert. High elasticity in both directions guarantees good elasticity and flexibility. This construction enables the mesh to be cut to any shape or size without unravelling.

The meshes' tensile strength and tear resistance are in line with the natural properties of the body, so the meshes can be used for various indications. The special structure, increased pore size and particularly fine monofilaments provide a macroporous, very light and flexible mesh implant. Its greater porosity allows more rapid ingrowth of the body's cells, which stabilises the tissues.

The mesh can be fixed in the location of the wound. For this purpose non-absorbable suture material should be used at a distance of at least 6.5 mm from the edge of the mesh.

Chiralen meshes are available in different sizes, all offering high strength and high flexibility. Chiralen meshes are available in a weight of approx. 125 g/m² or in the light version with a weight of approx. 55 g/m².

Indication: Chiralen meshes are indicated in the operative treatment of fascial defects, in particular hernias, which require long-term bridging of damaged structures.



Technical information - Chiralen mesh

Weight: 125 g/m²

**Thickness of
filament:** 0.420 mm

Pore density: 48 pores/1 cm²


Technical information - Chiralen mesh light

Weight: 55 g/m²

**Thickness of
filament:** 0.320 mm

Pore density: 48 pores/1 cm²

Pore size gross: 218 µm
(net: 180 µm)



CHIRALEN meshes light

Advantages:


- Minimises foreign-body reactions
- Security combined with great comfort for the patient
- Optimised scar tissue formation

CHIRALEN meshes

Size	Cat. No
6 cm x 11 cm	PM 0611
7.6 cm x 15 cm	PM 0715
10 cm x 15 cm	PM 1015
15 cm x 15 cm	PM 0151
30 cm x 30 cm	PM 0303

CHIRALEN meshes light

Size	Cat. No
6 cm x 11 cm	PL 0611
10 cm x 15 cm	PL 1015
15 cm x 15 cm	PL 1515
15 cm x 30 cm	PL 1530
30 cm x 30 cm	PL 3030



Veterinary



Veterinary

UK03
222



Veterinary

CHIRASORB violet

Bottles

EP	USP	Length	Cat. No
2	3/0	1 x 50 m	LV 1013
3	2/0	1 x 50 m	LV 1014
3.5	0	1 x 50 m	LV 1015
4	1	1 x 50 m	LV 1016
5	2	1 x 50 m	LV 1017
6	3 + 4	1 x 50 m	LV 1018
2	3/0	1 x 25 m	LV 1043
3	2/0	1 x 25 m	LV 1044
3.5	0	1 x 25 m	LV 1045
2	3/0	1 x 15 m	LV 1023
3	2/0	1 x 15 m	LV 1024
3.5	0	1 x 15 m	LV 1025
4	1	1 x 15 m	LV 1026
5	2	1 x 15 m	LV 1027
6	3 + 4	1 x 15 m	LV 1028

CHIRLAC violet


Bottles

EP	USP	Length	Cat. No
2	3/0	1 x 50 m	PG 1013
3	2/0	1 x 50 m	PG 1014
3.5	0	1 x 50 m	PG 1015
4	1	1 x 50 m	PG 1016
5	2	1 x 50 m	PG 1017
6	3 + 4	1 x 50 m	PG 1018
2	3/0	1 x 25 m	PG 1043
3	2/0	1 x 25 m	PG 1044
3.5	0	1 x 25 m	PG 1045
4	1	1 x 25 m	PG 1046
5	2	1 x 25 m	PG 1047
2	3/0	1 x 15 m	PG 1023
3	2/0	1 x 15 m	PG 1024
3.5	0	1 x 15 m	PG 1025
4	1	1 x 15 m	PG 1026
5	2	1 x 15 m	PG 1027
6	3 + 4	1 x 15 m	PG 1028

POLYDOX violet

Bottles

EP	USP	Length	Cat. No
2	3/0	1 x 50 m	DX 1013
4	1	1 x 50 m	DX 1016
2	3/0	1 x 25 m	DX 1043
3	2/0	1 x 25 m	DX 1044
3.5	0	1 x 25 m	DX 1045



Veterinary

SILON monofilament blue

Bottles

EP	USP	Length	Cat. No
1	5/0	1 x 50 m	SM 0702
1.5	4/0	1 x 50 m	SM 0703
2	3/0	1 x 50 m	SM 0704
3	2/0	1 x 50 m	SM 0705
3.5	0	1 x 50 m	SM 0706
4	1	1 x 50 m	SM 0707
5	2	1 x 50 m	SM 0708
3	2/0	1 x 100 m	SM 0725

TERVALON green

Bottles

EP	USP	Length	Cat. No
1.5	4/0	1 x 100 m	TG 0503
2	3/0	1 x 100 m	TG 0504
3	2/0	1 x 100 m	TG 0505
3.5	0	1 x 100 m	TG 0506
4	1	1 x 50 m	TG 0507
5	2	1 x 50 m	TG 0508
6	3 + 4	1 x 50 m	TG 0509
7	5	1 x 25 m	TG 0510
8	6	1 x 25 m	TG 0511

SILON braided white

Bottles

EP	USP	Length	Cat. No
1	5/0	2 x 200 m	501
1.5	4/0	2 x 200 m	502
2	3/0	2 x 140 m	503
2.5	-	2 x 140 m	504
3	2/0	2 x 100 m	505
3.5	0	2 x 70 m	506
4	1	2 x 50 m	507
5	2	2 x 40 m	508
6	3 + 4	2 x 30 m	509
7	5	2 x 25 m	510


SILK black

Bottles

EP	USP	Length	Cat. No
1	5/0	1 x 100 m	302 N
1.5	4/0	1 x 100 m	303 N
2	3/0	1 x 100 m	304 N
3	2/0	1 x 100 m	306 N
4	1	1 x 50 m	308 N
5	2	1 x 50 m	309 N
7	5	1 x 25 m	311 N

UK03
223

Veterinary





The CTI group



The CTI group

A constellation of 3 European manufacturers of high quality products for the healthcare sector

The CTI group consists of: CHIRANA T. Injecta in Slovakia, CHIRANA T. Injecta in Czech Republic and Vitrex Medical in Denmark.

The intention of CTI group is to serve as an “umbrella-organization”, where the individual companies collaborate to be “Committed to Innovation” and reach the joint goal:
- to develop, produce and distribute high quality, healthcare products.





CHIRANA T. Injecta,
Slovakia

CHIRANA T. Injecta is a global, state of the art manufacturer of healthcare products, which was founded in 1935.

In 1947 CHIRANA T. Injecta started manufacturing glass syringes and they have ever since focused on development and production of high quality, disposable products for the healthcare sector. Through continuously updating production facilities and expanding their product range, they are today supplying:

- Hypodermic syringes
- Perfusion syringes
- Pen and Blood collection needles
- Insulin syringes
- Cannula and lancet needles

CHIRANA T. Injecta,
Czech Republic

CHIRANA T. Injecta was founded in 1930. Since the early 1990s they have made substantial investments in upgrading the production resources, enabling the company to offer products in an unique quality.

Their more than 85 years of experience with manufacturing sutures has resulted in top priority within quality control, product development and sutures manufactured in line with state of the art technology. Their present range consist of:

- Surgical Sutures, Ligatures and Meshes
- Veterinary bottles and Eye-needles
- Microsurgery and Cardiovascular sutures
- Antibacterial sutures

VITREX Medical,
Denmark

Vitrex Medical was founded in 1965 and has since the 1970s mainly focused on supplying to the laboratory and analytical sector.

With more than 50 years of experience Vitrex Medical has become the largest manufacturer of capillary tubes in the world today by developing equipment, establishing technologies, and optimizing production processes. Vitrex Medical currently supplies:

- Capillary tubes and Micro-Hematocrit tubes
- Micro-pipettes
- Safety lancets
- Vacuum tubes and Blood collection needles
- Pen needles and Diabetic lancets



Version UK03

Designed by Vitrex Medical a/s

ref. 20170801JR03

CHIRANA T. Injecta, s.r.o.

**Komoranska 2148
143 00 Prague
Czech Republic**

www.cti-group.com